

**Thermo Scientific
ALPS™ 50V
User's Guide**



Important: This product is designed and constructed to be safe when properly used, in accordance with the supplied documentation, and when the operating precautions outlined in this document are fully observed.

IT IS ESSENTIAL THAT THE USER OF THIS MANUAL IS AWARE OF THE POTENTIAL HAZARDS ASSOCIATED WITH THE UNIT AND ITS ACCESSORIES. ALL OPERATORS SHOULD BE FAMILIAR WITH THE SAFETY PRECAUTIONS AND WARNINGS GIVEN IN THIS SECTION PRIOR TO ATTEMPTING TO OPERATE THE UNIT. IF THE UNIT IS USED IN A MANNER NOT SPECIFIED BY THE MANUFACTURER, THE PROTECTION PROVIDED BY THE EQUIPMENT MAY BE IMPAIRED.

The following symbols are used in this manual:



CAUTION/WARNING



HEATED SURFACE

The information in this user guide including any apparatus, Methods, Techniques and Concepts described herein, are the property of Thermo Fisher Scientific or its Licensees and may not be copied, disclosed or used for any purpose not expressly authorised by the owners thereof.

Since Thermo Fisher Scientific constantly strives to improve all of its products, We reserve the right to modify equipment and any user instruction manual without prior notice. No part of this instruction manual may be reproduced in any form without the prior consent of Thermo Fisher Scientific.

Thermo Fisher Scientific
ABgene House, Blenheim Road
Epsom, Surrey
KT19 9AP, UNITED KINGDOM

© 2007 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc.



WARNING: If the unit is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

The unit is intended for indoor laboratory use only, at an altitude of less than 2200m above sea level, within a temperature range of 18°C to 30°C and a relative humidity range of 20% to 80% non-condensing. If the instrument is stored outside these ranges, it should be left to stand until it equilibrates to within the above limits.

Ensure the voltage selection switch is set to the correct voltage for the required region (See Fig 3).

Ensure the correct fuse for the required voltage setting is fitted (see Specifications section).

Do not operate the unit outside the rated power supply range specified.

There are no user accessible or serviceable parts inside the unit. Do not remove or open the unit's casework.

Before using any cleaning or decontamination method except those recommended by the manufacturer, users should check with the manufacturer that the proposed method will not damage the equipment.

Ensure the unit is only connected to an earthed supply.



WARNING: The hot plate can reach temperatures up to 200°C. Care must be taken not to touch it or serious burns may occur.

The unit will remain hot for a considerable time after it has been turned off and must be allowed to cool down to an acceptable level before cleaning.

Foil sealed plates may remain hot for a number of seconds after sealing and should be handled with care.

Contents

1. General Description of Apparatus.....	Page 5
2. Intended function of Apparatus.....	Page 5
3. Installation.....	Page 5
4. Types of Consumable.....	Page 6
5. Operating Instructions.....	Page 7
6. Cleaning.....	Page 10
7. Specification.....	Page 10
8. Over Temperature Safety Cut Out.....	Page 10
9. Switching Off the Audible Warning (if fitted).....	Page 11
10. Lifting and Moving the Instrument.....	Page 11
11. Warranty and Returns.....	Page 11
12. Declaration of Conformity.....	Page 12

1. General Description of Apparatus

The Thermo Scientific ALPS 50V™ provides a safe, controlled method of sealing sample plates made from various plastics. The sealer is adjustable for both sealing time and temperature and gives a visual countdown display during the sealing process. The variable heat sealer can seal deep well plates with various thermal sealing plastic or foil films and may be adapted to seal Microplates and PCR plates by using a suitable plate carrier.

The variable heat sealer is intended to provide a platform with which to load and thermally seal sample plates to prevent sample loss at high temperature or prolonged periods of storage. Various types of plate may be loaded onto a plate carrier and covered by a thermal sealing medium, which may be foil or clear plastic.

2. Intended Function of Apparatus

The operator may pre-set a precise temperature and seal time for the sealing operation. Sealing is initiated by pulling down a handle which raises the sample plate into contact with the hot plate. The plate sealer detects when the sample plate is in the correct position and starts a countdown timer. The countdown is displayed on a three digit LED display. When the display reads zero, the operator lifts the handle back to the upright position, which moves the sample plate away from the hot plate. The plate may then be manually removed for further processing.

3. Installation

The variable temperature sealer should be removed from its packaging and all internal packaging removed. The instrument should be placed on a level surface, away from direct sunlight and drafts ensuring access to the power switch on the rear of the unit. Ensure that the vents on the back panel are not obstructed.

3.1 - Mains inlet selection:

The mains voltage selector switch should be checked to ensure that the desired voltage has been set and the correct value of fuse has been fitted (see FIG 3, FIG 4 and the specification section). The instrument should now be connected to an earthed mains supply with an IEC 320 connector type C13 mains lead ensuring there is sufficient space around it so that the mains connection is not obscured.

4. Types of Consumable

The ALPS 50V™ is designed to be used with the following consumables:-

Plate Types: (Polypropylene, polyethylene or polystyrene)

- Deep well plates without a carrier.
- PCR (Polymerase chain reaction) plates with carrier.
- Microplates with carrier.

Sealing film types:

- Foil polypropylene laminate
- Clear polyester polypropylene laminate
- Clear polymer
- Thin clear polymer
- Foil laminate
- Foil

Plate Carriers:

- Enclosed Well Plate Carrier S120530
- Open Well Plate Carrier S120529
- Bio Block Plate Carrier S120531

5. Operating Instructions

If the instrument has been removed from storage, it should be left to reach ambient room temperature. For correct operating conditions, please see specification section. Please refer to FIG 2 for location of user interface components and FIG 4 for the location of the power switch.

When the instrument is first turned on, The LED display shows the current setting for sealing temperature which is retained from the last seal. A new setting between 125°C and 200°C may be set by pressing the up and down arrow keys. The “°C” mode LED lights to indicate that the adjustable field is temperature. Pressing the mode switch causes the seal time to be displayed. This is also adjustable via the up and down keys in 0.5 Second increments between 1 – 9 seconds. Pressing the mode button again switches back to the sealing temperature adjustment. If the up or down arrow keys are held for more than 1 second, the key will auto-repeat.

When the desired settings have been set the heat on/off button may be pressed to begin heating the sealing plate to the set temperature. The heat LED flashes at a slow rate to indicate that the heater is on but not yet up to temperature. The heater should heat up to the set point in about ten minutes. When the heater is within two degrees of the set point, the heat LED will be continuously on and the sealer is ready to begin sealing plates.

5. Operating Instructions

A Typical operating sequence is:

1. Switch instrument on and set required sealing time and heater set-point temperature for the desired sealing process, via the keypad and LED display.
2. Enable the heater by pressing the Heat on/off button.
3. Allow sufficient time for the instrument to reach desired temperature. The Heat on/off LED will flash whilst the sealer is coming up to temperature and will remain on permanently when the set point ($\pm 2^{\circ}\text{C}$) is reached.
4. Load sample-plate onto plate carrier and add the sealing film. Care should be taken not to touch the heating surface whilst loading the sample plate.
5. Grasp the handle with one or two hands as appropriate (see Fig. 1a).
6. Lower handle to thermally compress sealing film onto sample plate (see Fig. 1b & 1c).
7. When correct clamping pressure is achieved feedback is provided and the “timer” will count down to zero. Care should be taken not to apply more pressure to the handle than is necessary to operate the micro switch.
8. On zero seconds, the display will flash and an audible warning will sound (if audible sounder fitted) . Release the heater plate from the sealed sample plate by raising the handle to its upper position and review the seal that has been made. Failure to raise the handle could ruin the sample plate.
9. Remove sealed plate for subsequent processing.

The seal temperature and time may vary with each type of plate. It is recommended that the temperature and time settings for the sealer are reviewed empirically to determine optimum seal integrity.

5. Operating Instructions

Fig. 1a



Fig. 1b



Fig. 1c



Fig. 2

1. Heat On/Off LED
2. Heat On/Off Button
3. Down Button
4. Up Button
5. Mode Button
6. Set Seal Time Mode LED
7. Set Temperature Mode LED




Fig. 3



Voltage Selector Switch Set to 115V



Voltage Selector Switch Set to 230V

Fig. 4



Mains Power Switch

Fig. 5



Fuse Position

6. Cleaning

The instrument is to be cleaned only by wiping Virkon solution at the manufacturer's guideline concentration. No parts are to be autoclaved.

7. Specification

Equipment Reference:
ALPS 50V™

Model No:
953-A

Specification:

Temperature Set Point	125°C to 200°C in 1°C increments.
Range	1 to 9 seconds in 0.5 second increments.
Time Set Range	20% to 80% non-condensing.
Operating Humidity Range	18°C to 30°C.
Operating Temperature	220 x 321 x 425 mm
Weight	7.2 kg
Range	100 - 130V AC 50/60Hz or 220 - 240V AC 50/60Hz
Power Supply	350 W Max
Power Rating	100 - 130V - T3.15A (IEC 127) or
Fuse	220 - 240V - T1.6A (IEC 127)
Mains Input Connector	IEC 320
Applicable Standards	FCC Class A, CFR 47 Part 15 BS EN 61326:1998

8. Over Temperature Safety Cut Out

The instrument is protected against a fault condition leading to the sealing plate heating up to a temperature exceeding its maximum set point. Protection is provided by a single use thermal cut out, operating in the range 235°C to 240°C. In the unlikely event that this cut out should operate, it would necessitate the sealer being returned to the manufacturer for repair.

9. Switching Off the Audible Warning

Some units are fitted with a sounder which give an audible warning at the end of the seal time. The sounder is on by default but may be programmed to be on or off. To turn the Audible warning off, place the sealer in Temperature set mode (Set Temperature Mode LED on) and simultaneously press the Mode button and the Heat On/Off buttons together. The display will show “01” to indicate that the audible warning is ON and both Set Temperature Mode and Set Seal Time Mode LEDs will be off. Pressing the Down button turns the audible warning OFF. The display should now show “00” . Pressing the Mode button will save this configuration so that the audible warning is always off and the display will return to the Set Temperature Mode. To turn the audible warning back on repeat this process but instead press the Up button to display “01”.

10. Lifting and Moving the Instrument



WARNING: The instrument should be lifted by holding by the handle. Care should be taken when lifting the instrument due to its weight. Refer to current Health and Safety procedures.

11. Warranty and Returns

Thermo Fisher Scientific
ABgene House, Blenheim Road
Epsom, Surrey
KT19 9AP, UNITED KINGDOM

+44 (0)1372 723456
+44 (0)1372 741414 fax
abgene.info@thermofisher.com
www.thermo.com

12. Declaration of Conformity

The ALPS™ 50V unit has been designed in accordance with, and satisfies the requirements of, article 11 of the Low Voltage Directive 73/23/EEC as realigned by 93/68/EEC on the harmonisation of the laws of the Member States relating to electrical equipment designed for use within certain voltage limits, to the essential requirements of BS EN 61010-1:2001.

The ALPS™ 50V unit has been type tested by EMC Projects Limited (a UKAS and CAA approved test facility and UK appointed Notified Body), and issued a Certificate of Compliance No. 6091/06 to the following EMC standard:

BS EN 61326:1998, Electrical equipment for measurement, control and laboratory use. EMC requirements

Satisfying the EMC Directive(s) 89/336/EEC and 92/31/EEC as realigned by 93/69/EEC.

Declaration:

I declare that as the Authorised representative, that the above information in relation to the supply/manufacture of this product, is in conformity with the stated standards and other related documents, following the provisions of (89/392/EEC)

Signed:

A handwritten signature in black ink, appearing to be 'T. J. ...', written over a faint, illegible printed name.

Thermo Fisher Scientific
United Kingdom

Notes:

Notes:

Notes:

Thermo Fisher Scientific
ABgene House, Blenheim Road
Epsom, Surrey
KT19 9AP, UNITED KINGDOM

www.thermo.com