
Operating Instructions | Betriebsanleitung | Mode d'emploi |  
Instrucciones de funcionamiento | Istruzioni per l'uso

## Sartorius Portable Meter

PT-10, PT-10P, PT-10-P20

Portable pH/mV Meter | Tragbares pH/mV-Meter |  
pH-mètre/millivoltmètre portable | medidor pH/mV portátil | pH/mV-metri portatili

---

**English – page 3**

In cases involving questions of interpretation, the German-language version shall prevail.

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---

# Contents

3	<b>Contents</b>	16	<b>Operating the Meter: pH Mode</b>
4	<b>Intended Use</b>	18	<b>Standardizing the Meter for Millivolt Measurement</b>
5	<b>Warning and Safety Information</b>	19	<b>Operating the Meter: Toggling between Relative mV Mode and Absolute mV Modes</b>
6	<b>General View of the Equipment</b>	20	<b>Temperature Compensation</b>
7	<b>Getting Started</b>	21	<b>Troubleshooting</b>
7	Equipment Supplied	22	<b>Overview</b>
8	Installing the Battery	22	Specifications
9	Connecting the Equipment to AC Power (Optional)	23	Accessories
10	Installing Electrodes	24	Declaration of Conformity
13	<b>Standardizing the Meter for pH Measurement</b>		
16	<b>Operating the Meter: Measurement</b>		

---

## Intended Use

The PT Portable Meters from Sartorius measure pH, mV and temperature within the following ranges:

- pH range: 0.00 to 14.00
- Measurement range in mV: -1800 to +1800
- Measurement range in °C: -5 to +105

PT portable meters ease your workload and speed up your routine with the following features:

- Simple operation
- Rugged, sturdy construction
- Battery operation using a 9-volt battery (optional AC adapter available separately)
- Easy to standardize with 1, 2 or 3 standardization points
- Automatic buffer detection
- Automatic electrode test during standardization
- Automatic temperature compensation
- Clear readout with easy-to-understand symbols and LCD

---

## Warning and Safety Information

The PT portable meter complies with the European Council Directives as well as international regulations and standards for electrical equipment, electromagnetic compatibility, and the stipulated safety requirements. Improper use or handling, however, can result in damage and/or injury.

- To prevent damage to the equipment, read these operating instructions thoroughly before using your portable meter. Keep these instructions in a safe place. Please observe the following to ensure safe and reliable operation of your portable meter:

△ Do not use this equipment in hazardous areas/locations.

△ If you use electrical equipment in installations and under ambient conditions requiring higher safety standards, you must comply with the provisions as specified in the applicable regulations for installation in your country.

△ If you use an AC adapter, make sure the voltage rating printed on the AC adapter is identical to your local line voltage.

- Use Sartorius accessories, as these are optimally designed for use with your portable meter.
- The buffer solution used for standardization must match the stored values exactly.


- Note on installation:  
The operator shall be responsible for any modifications to Sartorius equipment and must check and, if necessary, correct these modifications. On request, Sartorius will provide information on the minimum operating specifications (in accordance with the Standards listed above for defined immunity to interference).

Industrial protection ratings for the housing:

- Portable meter: IP67
  - AC adapter (optional): IP20
- If you have any problems with your portable meter: contact your local Sartorius office, dealer or service center.


# General View of the Equipment

Display and keypad:


- 1 LC display
- 2 [Enter] key: Select options
- 3 [Standardize] key: Activate standardization functions
- 4 On/off
- 5 Jack for optional AC adapter
- 6 "Measuring" symbol
- 7 Temperature
- 8 Mode
- 9 Result
- 10 Texts: settings/error indicators
- 11 Buffer symbols
- 12 Symbols: Standardization in process; slope symbol
- 13 "Stability" symbol
- 14 "Standardize" symbol
- 15 Low-battery indicator
- 16 Connector for ATC probe
- 17 BNC connector for electrode
- 18 [Mode] key: Toggle between pH and mV mode
- 19 [Setup] key: Activate device settings (e.g. to delete standardization data, etc.)

Rear view:

Display:

---

# Getting Started

## Storage and Shipping Conditions

Do not expose the portable meter unnecessarily to extreme temperatures, moisture, shocks, blows or vibration.

## Unpacking the Portable Meter

- After unpacking the portable meter, please check it immediately for any visible damage.
- If any sign of damage is visible, proceed as directed in the chapter entitled “Care and Maintenance,” under the section on “Safety Inspection”.

It is a good idea to save the box and all parts of the packaging until you have successfully installed your portable meter. Only the original packaging provides the best protection for shipment. Before packing your portable meter, unplug all connected cables to prevent damage.

## Equipment Supplied

The equipment supplied includes the components listed below:

### PT-10


- pH meter
- 9V battery



### PT-10P

- pH meter
- 9V battery
- Gel electrode with temperature sensor (PY-P12)
- Carrying case
- Buffer solution; pH = 4.00, 90 ml
- Buffer solution; pH = 7.00, 90 ml
- 2 PE bottles, 60 ml

### PT-10-P20

- pH meter
- 9 V battery
- Gel electrode without temperature sensor
- Carrying case
- Buffer solution; pH = 4.00, 90 ml
- Buffer solution; pH = 7.00, 90 ml
- 2 PE bottles, 60 ml


### Installing the Battery

- Use only commercially available 9 V batteries. If desired, use a rechargeable battery (not included).
- Rechargeable 9 V batteries are available from your local electronic supplies dealer.
- △ If you use a rechargeable battery, recharge it with an external battery charger.
- Lay the portable meter face down.
- △ Open the battery compartment: Remove the compartment cover.
- Install the battery in the compartment.
- Make sure the polarity is correct.
- △ Do not discard used batteries with normal household waste. Rechargeable batteries contain toxic materials and must be disposed of in accordance with local waste disposal regulations.
- Close the battery compartment: Slide the cover into position until it snaps into place.

### Note:

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense. Changes or modifications not expressly approved by Sartorius AG could void the user's authority to operate the equipment.

### Connecting the Equipment to AC Power (Optional)

- Use only original Sartorius equipment; in Europe:

Power adapters for PT Meters

230 V European plug 102078.1

240 V UK plug 102079.1

100 V Japanese plug 102155.1

230 V Australian plug 102156.1

115 V North American plug 102077.1

This AC power adapter has protection rating IP20 in accordance with the European standard EN 60529.

- Plug the AC adapter into the wall outlet (mains supply).
- The pole of the output voltage is connected to the portable meter housing, which can be grounded for operation.



# Installing Electrodes

You can use the following electrodes with the PT portable meter:


- Combination pH electrode with BNC connector
- Separate 2.5 mm phone jack for the ATC probe
- Ion-selective or ORP (redox) electrode with BNC connector

## Installation

1. Remove the protective cover from the electrode.

2. If the pH electrode is dry (e.g., when it is new), soak it overnight in an electrode filling solution, 3-mol/l KCl solution or electrode storage solution before use.


3. Check the device for correct functioning as follows:

Connect the shorting cap to the BNC input on the portable meter.


Select mV mode and record the mV value:

Press the [Mode] key


Make sure the portable meter is in mV mode, not relative mV mode. If the mV value is  $0 \pm 0.3$ , the meter is functioning correctly. Note that the long-term drift specified for this portable meter is 0.1 mV/month since the most recent standardization.

Connect the electrode to the input connector on the meter (push in and twist to lock). For high-precision measurements, you can connect an ATC probe to the ATC port on the portable meter.


4. Rinse and blot-dry electrodes between measurements (do not wipe). Rinse electrodes with distilled water, deionized water, or some of the next solution to be measured.

5. Store pH electrodes in KCl solution or electrode storage solution. Always leave the filling hole open when in use and closed when in storage. Refill with filling solution when the internal solution drops to 1 inch below the filling hole.


⚠ PT-10P/PT-10-P20 models: The electrode supplied is a low-maintenance gel-filled electrode which cannot be re-filled.

# Standardizing the Meter for pH Measurement


Because electrodes vary in their response, you need to standardize your pH meter and electrodes to compensate for variations. The more frequently you standardize, the more accurate your measurements. For best results, standardize at least once a day.

1. Immerse the electrode in standard solution. Stir gently. Allow the electrode to reach a stable value.

2. If the display is not in pH mode: Press the [Mode] mode key to toggle to pH mode.


**[Setup]**

3. Clear the stored standardization data before performing a new 2- or 3-point standardization: Press the [Setup] key and then press [Enter]. Existing standardization data is deleted and you can select a new set of buffers.


**[Enter]**

## [Standardize]


4. Display the current buffer set:  
Press the [Standardize] key.

The current buffer set flashes on the display and the meter automatically detects the current buffer from among this set. Once the signal stabilizes, the buffer data is stored automatically. When performing standardization, older data is deleted automatically.


5. The meter displays the slope percentage of the electrode for the first buffer as 100.0%. When you enter a second or third buffer, the meter performs a diagnostic check on the electrode and displays the slope.


## [Standardize]




6. To enter a second buffer, place the electrode in the second buffer solution, stir, allow to stabilize, and press [Standardize]. The meter detects the buffer.

7. Next, the meter performs a diagnostic test of the electrode. The display shows the slope symbol and the slope % of the electrode.

## [Standardize]

8. If “Error” is displayed next to the slope symbol, your electrode is not functioning properly. The slope must be between 90 and 105%. The portable meter does not accept or store values displayed with the “Error” indicator. Press [Enter] to continue.

If no ATC probe is connected during calibration, the meter uses the nominal buffer values at 25°C. If an ATC probe is connected, the meter adjusts continually for temperature. In this case, the pH value of the buffer at the current temperature is taken into account when calculating the slope of the electrode.


9. Enter a third standard:  
Immerse the electrode in the third buffer solution, stir, allow to stabilize, and press [Standardize]. Refer to steps 7 and 8 before continuing.

10. After each buffer is entered, the “Standardizing” symbol switches off and the “Measuring” or “Stability” symbol is displayed, indicating that the meter has returned to the normal measuring mode.

Standardize your meter and electrode using at least two buffers, with pH values bracketing the expected pH of your samples.

Note: Stirring shortens the electrode response time.

## Operating the Meter: Measurement

Place the electrode in the medium to be measured so that both the glass membrane and the diaphragm are immersed.


## Operating the Meter: pH Mode

Function of the [Setup] key:

Press the [Setup] key to access functions for displaying or deleting standardization data, as well as for selecting the desired set of buffers.


Exiting the Setup mode: Press the [Setup] key again.

### [Setup]

1. To delete all standardization data:  
Press [Setup]; "Clear" is displayed.  
To delete all values, press the [Enter] key.  
The meter deletes all standardization data and returns to the measuring mode.


### [Enter]


2. Display of the electrode slope:  
Following successful standardization with 2 buffers, "OK" is displayed, indicating that the electrode tested is OK. Additionally, the slope between the first and second buffers is displayed (a numeric value and the "Slope" symbol) along with the corresponding buffer symbols.


**2× [Setup]**

**[Setup]**

**[Setup] or [Enter]**

**[Setup]**

1.68 4.01 6.86 9.18 12.46

**[Setup]**

"SET Auto"

**[Mode]**

"SET Man"

**[Enter]**

3. When three buffers are used:  
To view the electrode slope between the second and third buffers, press the [Setup] key again.  
The electrode slope between the second and third buffers is displayed.

4. To view stored buffer sets:  
Press the [Setup] key again.

5. To select the displayed buffer set:  
Press the [Enter] key.  
To select the next set:  
Press the [Setup] key.


6. To select NIST buffer sets: Press the [Enter] key.

7. To select "Always shut off manually":  
Press the [Mode] key and then press [Enter] to confirm or select "Automatic shutoff" (factory setting; turns off the meter if not used for 20 minutes; available only for battery operation).

8. To return to the normal measuring mode:  
Press the [Setup] key.

To confirm all settings: Press the [Enter] key.


# Standardizing the Meter for Millivolt Measurement

Millivolt measurement is used in electrode diagnostics, titration, and in measuring oxidation/reduction (redox) potential (ORP). ORP measurements indicate the oxidizing or reducing capacity of a solution. You can use ORP values to monitor and control solutions that require a defined level of oxidants or reductants.


1. Immerse the electrode in standard solution.

## [Mode]


2. If the display is not in mV mode: Press the [Mode] mode key to toggle to mV mode.

## [Standardize]

Standardizing

3. Set a relative mV value:  
Press the [Standardize] key.


4. When the signal has stabilized:  
The meter automatically stores the measured value.  
The current absolute mV value becomes zero relative millivolts.

---


## Operating the Meter: Toggling between Relative and Absolute mV Modes

**[Setup]**

1. Delete the mV offset and return to the absolute mV mode: Press the [Setup] key. "Clear" is displayed.

**[Enter]**

2. Delete the mV offset:  
Press the [Enter] key.  
The meter then returns to the absolute mode.

---

# Temperature Compensation

△ Note: The temperature is compensated automatically only if an ATC probe is connected.

Temperature has two major effects on pH values:

1. Temperature affects buffer pH values.

The pH value in a buffer is dependent on the temperature of the buffer. The values are usually given on the buffer label. Most technical buffers correspond to the values listed in the table shown here on the left.

When you perform standardization in pH mode, the pH value is adjusted to the current value for the current temperature.

## Standard buffer:

	pH 4.00	pH 7.00	pH 10.00
0 °C	4.005	7.13	10.34
5 °C	4.003	7.10	10.26
10 °C	4.001	7.07	10.19
15 °C	4.002	7.05	10.12
20 °C	4.003	7.02	10.06
25 °C	4.008	7.00	10.00
30 °C	4.010	6.99	9.94
35 °C	4.020	6.98	9.90
40 °C	4.030	6.97	9.85
50 °C	4.061	6.97	9.78


For example, if the buffer has a pH of 7.00 at 25°C and the current buffer temperature is 20°C, the portable meter standardizes the buffer at pH 7.02.

2. Temperature affects electrode efficiency.


The theoretical voltage change per pH unit is approximately 59.17 mV at 25°C. This value changes, however, as the temperature changes.

The meter compensates for this effect by changing the pH reading as the temperature changes


# Troubleshooting

1. If "-----" is displayed, the signal from the electrode is out of range.  
This might indicate that the electrode is not in the solution.

2. If "Error" is displayed with the "Slope" symbol, this indicates an error in the electrode reaction. If this is displayed during standardization, this means the electrode shows under 90% or over 105% of the correct slope. The "Error" message can mean one of the following:
  - bad electrode
  - bad buffer(s)
  - wrong buffer selected

3. If "---" is displayed, the meter has detected an error in the ATC probe. This indicates that no ATC probe is used for temperature compensation (default temperature: 25°C).

## Electrode Test

pH=7  $0 \pm 30$  mV

pH=4 159 to 186 mV  
higher than pH 7

pH=10 159 to 186 mV  
lower than pH 7

4. Testing a pH-electrode:  
Immerse the electrode in a good pH 7 buffer. Press the [Mode] key to activate mV mode and record the displayed value. Make sure the meter is in mV mode and not in relative mV mode. Repeat the same steps with a pH 4 or pH 10 buffer. At a temperature of approximately 25°C, the electrode signal should be within the limits shown here.

---

# Overview

## Specifications

### PT-10 Portable Meter

#### pH Measurement

Measurement range	0.00 to 14.00
Readability	$\pm 0.01$
Accuracy	$\pm 0.005$
Standardization points, max.	3
Automatic buffer detection	16 buffers (NIST, technical buffers)
Automatic temperature compensation	Yes

#### mV Measurement

Measurement range in mV	-1800 to +1800
Readability in mV	$\pm 0.1$
Accuracy in mV	$\pm 0.2$

#### Temperature measurement

Temperature range in °C	-5 to +105
Readability in °C	$\pm 0.1$
Accuracy in °C	$\pm 0.3$

# Accessories

	Order No.
<b>AC adapter</b> for Europe	102078.1
Power adapters for PT Meters	
230 V European plug	102078.1
240 V UK plug	102079.1
100 V Japanese plug	102155.1
230 V Australian plug	102156.1
115 V North American plug	102077.1
<b>pH Combination Electrodes:</b>	
Fiber junction, gel electrolyte, integrated ATC probe	PY-P12
Fiber junction, gel electrolyte	PY-P20
Fiber junction, liquid KCl electrolyte, integrated ATC probe	PY-P10
Platinum junction, liquid KCl electrolyte, integrated ATC probe	PY-P11
Platinum junction, liquid KCl electrolyte	PY-P21
Platinum junction (microelectrode), liquid KCl electrolyte	PY-P22
Annular gap junction (surface electrode), gel electrolyte	PY-P23
Adjustable sleeve junction, liquid KCl electrolyte	PY-P24
<b>Redox Combination Electrode:</b>	
Platinum disk, liquid KCl	PY-R01
<b>Temperature Probe:</b>	
10 k $\Omega$ ATC	PY-T01
<b>pH Buffer Solutions:</b>	
Buffer pH = 4.00 $\pm$ 0.01 at 25°C, 6 x 90 ml	PY-Y21-6
Buffer pH = 7.00 $\pm$ 0.01 at 25°C, 6 x 90 ml	PY-Y22-6
Buffer pH = 4.00 $\pm$ 0.01 at 25°C, pump bottle, 500 ml	PY-Y21
Buffer pH = 7.00 $\pm$ 0.01 at 25°C, pump bottle, 500 ml	PY-Y22
Buffer pH = 10.00 $\pm$ 0.01 at 25°C, pump bottle, 500 ml	PY-Y23
<b>Cleaning solution</b> , Pepsin/hydrochloric acid, 500 ml	PY-Y06
<b>Electrolyte solution</b> , KCl (3mol/l), silver ion-free, 500 ml	PY-Y07

# Declaration of Conformity to Council Directives 89/336/EEC and 73/23/EEC

The electronic pH meter of the series PT-..

meets the requirements of the test standards listed below, in conjunction with the associated power supplies, auxiliary peripheral devices and installation equipment listed in Annex A2 (see Annex A1 for a technical description and variants).

## 1. Electromagnetic Compatibility

1.1 Source for 89/336/EEC: EC Official Journal, No. 2001/C105/03

EN 61326-1 Electrical equipment for measurement, control and laboratory use- EMC requirements  
Part 1: General requirements

Limitation of emissions: Residential areas, Class B

Defined immunity to interference: Minimum test requirements, continuous monitored operation

## 2. Safety of Electrical Equipment

2.1 Source for 73/23/EEC: EC Official Journal, No. 2001/C106/03

EN 61010 Safety requirements for electrical equipment for measurement, control and laboratory use  
Part 1: General requirements

Sartorius AG  
37070 Goettingen, Germany  
2002

J. Rehwald  
(Senior Vice President  
Mechanical Engineering  
Mechatronics Division)

C. Oldendorf  
(Vice President  
R&D, Technological Operations and Innovations  
Mechatronics Division)

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Copyright by Sartorius AG,  
Goettingen, Germany.  
All rights reserved. No part  
of this publication may  
be reprinted or translated in  
any form or by any means  
without the prior written  
permission of Sartorius AG.  
The status of the information,  
specifications and illustrations  
in this manual is indicated  
by the date given below.  
Sartorius AG reserves the  
right to make changes to the  
technology, features,  
specifications and design of the  
equipment without notice.

Status:  
June 2002, Sartorius AG,  
Goettingen, Germany