

# IKAMAG® RET control-visc C

*safety control*

# IKAMAG® RET control-visc

*safety control*

BETRIEBSANLEITUNG	DE	3
OPERATING INSTRUCTIONS	EN	26
MODE D'EMPLOI	FR	49
INDICACIONES DE SEGURIDAD	ES	72
VEILIGHEIDSINSTRUCTIES	NL	74
NORME DI SICUREZZA	IT	76
SÄKERHETSANVISNINGAR	SV	78
SIKKERHETSHENVISNINGER	DA	80
SIKKERHETSHENVISNINGER	NO	82
TURVALLISUUSOHJEET	FI	84
INSTRUÇÕES DE SEGURANÇA	PT	86
WSKAZÓWKI BEZPIECZEŃSTWA	PL	88
BEZPEČNOSTNÍ UPOZORNĚNÍ	CS	90
BIZTONSÁGI UTASÍTÁSOK	HU	92
VARNOSTNA NAVODILA	SL	94
BEZPEČNOSTNÉ POKYNY	SK	96
OHUTUSJUHISED	ET	98
DROŠĪBAS NORĀDES	LV	100
SAUGOS REIKALAVIMAI	LT	102

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DE

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## EG-DECLARATION OF CONFORMITY

EN

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 61010-1 and EN 61326.

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In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

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## 1 Warranty

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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# Contents

	Page
<b>1 Warranty</b>	<b>2</b>
<b>2 Safety instructions</b>	<b>27</b>
<b>3 Proper use</b>	<b>28</b>
<b>4 Unpacking</b>	<b>28</b>
<b>5 Commissioning</b>	<b>29</b>
<b>6 Control elements</b>	<b>29</b>
<b>7 Display</b>	<b>29</b>
<b>8 Switching on</b>	<b>30</b>
<b>9 Default settings</b>	<b>30</b>
9.1 Resetting parameters to factory default settings (RESET)	31
<b>10 Modes of operating</b>	<b>31</b>
10.1 A Mode of operating	31
10.2 B Mode of operating	31
10.3 C Mode of operating	31
10.4 X1 Mode of operating	31
10.5 Reversing the modes of operating	31
<b>11 Function heating</b>	<b>32</b>
11.1 Control of medium temperature with Pt 100 temperature sensor	32
11.2 Medium temperature control with twin Pt 1000 temperature sensors	32
11.3 Medium temperature control without external temperature sensors	33
11.4 Residual heat warning (HOT)	33
11.5 Changing the temperature unit of measure- ment from °C; °F	34
11.6 Setting the safety temperature	34
<b>12 Function stirring</b>	<b>34</b>
12.1 Intermittent mode	34
12.2 Stirring bar breakaway monitoring	35
12.3 Setting the intermittent mode and the stirring bar breakaway monitoring function	35
12.4 Measurement of viscosity trend	35

	Page
<b>13 Tips</b>	<b>36</b>
<b>14 Connections and interfaces</b>	<b>36</b>
14.1 RS 232 C serial interface	37
14.1.1 Configuration	37
14.1.2 Instruction syntax	37
14.1.3 Overview of the NAMUR-instructions	37
14.1.4 "Watchdog" -function, monitoring the serial data flow	40
14.2 Connections between RET control-visc /C and external devices	40
14.2.1 PC 5.1 Adapter IKA - Control	40
14.2.2 PC 2.1 Cable	40
14.2.3 PC 1.5 Cable	40
14.2.4 PC 1.2 Adapter	41
14.2.5 AK 2.1 Cable	41
14.2.6 AK 2.2 Cable	41
14.3 Tips on using labworldsoft	41
<b>15 Faults</b>	<b>41</b>
15.1 Troubleshooting. Error messages, causes and solutions	42
<b>16 Maintenance and cleaning</b>	<b>45</b>
16.1 RET control visc C Technical Enamel	45
16.1.1 Protective conductor RET control visc C	45
<b>17 Associated standards and regulations</b>	<b>45</b>
<b>18 Accessories</b>	<b>46</b>
<b>19 Technical data</b>	<b>47</b>
<b>20 List of spare parts RET control visc /C</b>	<b>106</b>
<b>21 Spare parts diagram RET control visc /C</b>	<b>107</b>

## 2 Safety instructions

EN

### *For your protection*

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Socket must be earthed (protective ground contact).
- **Caution - Magnetism!** Effects of the magnetic field have to be taken into account (e.g. data storage media, cardiac pacemakers ...).
- **Risk of burns!** Exercise caution when touching the housing parts and the heating plate. The heating plate can reach temperatures in excess of 340 °C. Pay attention to the residual heat after switching off.
- Ensure that the mains power supply cable does not touch the heating base plate.
- Wear your personal protective equipment in accordance with the hazard category of the media to be processed. Otherwise there is a risk from:
  - splashing and evaporation of liquids
  - ejection of parts
  - release of toxic or combustible gases.
- Set up the appliance in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The feet of the appliance must be clean and undamaged.
- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- Gradually increase the speed.
- Reduce the speed if
  - the medium splashes out of the vessel because the speed is too high
  - the appliance is not running smoothly
  - the container moves on the base plate.
- **Caution!** Never process and heat up any media that has a flash point higher than the adjusted safe temperature limit that has been set (50 to 350 °C).

The safe temperature limit must always be set to at least 25 °C lower than the fire point of the media used.


- Beware of hazards due to:
  - flammable materials
  - combustible media with a low boiling temperature
  - glass breakage
  - incorrect container size
  - overfilling of media
  - unsafe condition of container.
- The appliance may heat up when in use.
- The base plate can heat up due to the action of the drive magnets at high motor speeds, even if the heater is not operational.
- Process pathogenic materials only in closed vessels under a suitable extractor hood. Please contact **IKA** if you have any questions.
- **Do not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Please observe the operating instructions for any accessories used.
- Ensure that the external temperature sensor (PT 1000, PT 100, ETS-D ...) is inserted in the media to a depth of at least 20 mm.
- The PT 1000 external temperature sensor must always be inserted in the media when connected.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Accessories must be securely attached to the device and cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.
- Always disconnect the plug before fitting accessories.
- The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- The device will automatically restart in mode B following any interruption to the power supply.
- It may be possible for wear debris from rotating accessory parts to reach the material being processed.

- When using PTFE-coated magnetic bars, the following has to be noted: *Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.*

(Source: Römpps Chemie-Lexikon and "Ulmann", Volume 19)

#### For protection of the equipment

- The appliance may only be opened by experts.
- The voltage stated on the type plate must correspond to the mains voltage.
- Do not cover the device, even partially e.g. with metallic plates or film. This results in overheating.
- Ensure that the base plate is kept clean.
- Protect the appliance and accessories from bumps and impacts.
- Only use glass encapsulated temperature sensors for electrolysis procedures.
- Observe the minimum distances between the devices, between the device and the wall and above the assembly (min. 800 mm).

### 3 Proper use

- **Use**
  - For mixing and/or heating liquids
- **Range of use**
  - Laboratories
  - Schools
  - Pharmacies

This device is suitable for use in all areas except:

- Residential areas
- Areas that are connected directly to a low-voltage supply network that also supplies residential areas.

### 4 Unpacking



Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The delivery scope covers:  
an IKAMAG RET control vis /C, an operating instructions, an main cable and a short instruction.

## 5 Commissioning


Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). Insert the temperature measuring sensor (not included in the supply schedule) into the 5-polar DIN-bushing at the rear of the instrument. Put a magnetic bar into a vessel that is permeable for magnetic field lines, and fill the vessel with the medium to be processed. Put the vessel on the heating plate and immerse the temperature measuring sensor into the medium (min. 20mm). With the enterprise without external temperature measuring sensor see chapter 11.3. If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged. Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

## 6 Control elements

## 7 Display

When the functions are switched off, the digital display always shows the set speed value.

In all three modes of operation only function "heating" has been started, the display shifts to rated and real temperature.


In all three modes of operation only function "stirring" has been started, the display shifts to rated and real speed.

If both functions have been started, the display of function „heating“ always has preference. If in this case the speed is changed via rotary knob „Mot“ (not possible in **C** -mode of operation), the display reverses to speed display for a duration of 5 seconds.

The mode of the display can also be changed by pressing foil key „Mode“ once. After 5 seconds the display shifts back to temperature display.

If function „heating“ is switched off, and if the heating plate temperature is higher than 50°C, LED „Temp“ will be blinking green (20% ON ... 80%OFF).

The digital display changes constantly from rated temperature **°C**

to rated temperature **HOT**

## 8 Switching on

The instrument is switched on by the mains switch on the left side of the instrument (see Ill. 1).

A self test cycle for all the safety functions is run and a recognition cycle of the temperature sensors connected up is also run.

- The “Temp” and “Mot” LEDs and the LCD display segments flash for approximately 2 seconds.
- The software version number (e.g. “V000”) is displayed.
- The status of the intermittent mode (**On/OFF Intermit**) and the stirring bar breakaway monitoring (**On/OFF SAFE Rpm**) functions are displayed.
- The message **TEST Probe** is displayed.
- The operating mode (**A, B** or **C**) is displayed.

- The temperature sensors connected up and identified are displayed:

**PT OFF** No external temperature sensor connected up.  
**PT2 OFF** PT100 (single)  
**PT2 On** 2x PT1000

- The safe temperature limit of the heating plate is displayed as **SAFE H xxx °C** (SV flashes).

- **SAFE H OK** is displayed.

- If external temperature sensors (PT100 or PT1000) are connected up, their safe temperature limit is displayed as **SAFE P xxx °C** (SV flashes).

- The message **SAFE P OK** is displayed.

## 9 Default settings

The RET control-visc /C is supplied with the following default settings:

Operating mode	<b>A</b>	
Setpoint temperature	<b>0</b>	°C
SAFE H	<b>350</b>	°C
SAFE P	<b>350</b>	°C
Setpoint speed	<b>50</b>	rpm
Intermittent stirring mode	<b>OFF</b>	
Cycle time (at setpoint speed)	<b>60</b>	sec. (t)
Interruption (at 50 rpm)	<b>10</b>	sec. (t)
Stirring bar breakaway monitoring	<b>OFF</b>	
Error 5 response time	<b>180</b>	sec. (t)
Pt100 temperature sensor offset	<b>0</b>	K (t)
Pt1000 1 temperature sensor offset	<b>0</b>	K (t)
Pt1000 2 temperature sensor offset	<b>0</b>	K (t)

(t = can only be modified via serial interface)

## 9.1 Resetting parameters to factory default settings (RESET)

To reset the parameters to the default settings, proceed as follows:

- Switch off the unit at the main ON/OFF switch.
- Hold the membrane buttons "Visc" and "Temp" pressed down and switch on the main ON/OFF switch. After approx. 10 seconds the display will begin to flash. Release the membrane buttons.
- The unit is now set on the factory default settings.

## 10 Modes of operating

The IKAMAG RET control visc can be operated in the following three modes of operation:

### 10.1 A Mode of operation

The IKAMAG RET control visc /C is switched at our Works to the **A** -mode of operation. In this mode, an operation is only possible when **the Pt temperature measuring sensor** (Pt100 or Pt1000) is inserted. When the instrument is switched on („power“), functions heating and stirring are switched off.

The set rated values are stored and digitally displayed. They are taken over when the respective function is switched on. The set rated values can be changed. If there has been an interruption of the mains, **both functions have to be restarted**.

The safety circuit can be set or modified. When the power supply switch is switched on, **A** is displayed.

### 10.2 B Mode of operation

In this mode of operation can be worked with or without **Pt temperature measuring sensor** (Pt100 or Pt1000). When the instrument is switched on, the state of the functions heating and stirring before the last switching off (mains OFF) as well as

the data set last are taken over. The set rated values can be changed. The safety circuit can be set or modified.

The **B** - mode of operation which is indicated by symbol „O“. When the power supply switch is switched on, **B** is displayed.

### 10.3 C Mode of operation

In this mode of operation can be worked with or without **Pt temperature measuring sensor** (Pt100 or Pt1000). When the instrument is switched on, the state of the functions heating and stirring before the last switching off (mains OFF), as well as the data set last are taken over. The set rated values from **A** -and **B** -modes of operating cannot be changed.

The safety circuit **cannot** be set or modified. When the power supply switch is switched on, **C** is displayed.

### 10.4 X1 Mode of operation

In this mode of operation X1 the equipment behaves as in the mode of operation A with the following exceptions:

- The allocation of the two Pt1000 temperature measuring sensors is firm.
- No Error 5, if the Pt1000 temperature measuring sensor from the reaction container is taken.

### 10.5 Reversing the modes of operation

For reversing the modes of operation, the following steps are required:

- Switch off the instrument by way of mains switch.
- Keep foil key „Mot“ and „Temp“ depressed and switch instrument on via mains switch, release foil keys after about 1 second.

## 11 Function heating

The unit has a built-in 600W heating plate with a stainless steel or coated surface (C). Heating is controlled by a wear-proof, electronic TRIAC element.

The heating plate is kept at a constant temperature by a control circuit. The heating plate temperature can also be monitored from a separate, adjustable safety circuit. The two temperature sensors (Pt1000) for this purpose are built into the heating plate.

The heating plate safety circuit (SAFE H (hotplate)) can be set in the temperature range from 50° to 350°C.

The unit can also be controlled from a single, external Pt100 temperature sensor or from twin Pt1000 temperature sensors. Any external temperature sensors must be plugged in before the unit is switched on so that the automatic self-test/sensor recognition cycle can take place.

The external temperature sensors (medium) are monitored by a separate, adjustable safety circuit. The medium safety circuit (SAFE P (probe)) can be set from 50 to 350°C or to the maximum setting of the heating plate safety circuit.

Use the "Temp" knob to set the setpoint temperature required in the range from 0° to 340°C (max.) The temperature set will be displayed on the upper right-hand side of the digital display. "SV" and "PROBE" will flash.

To switch the heating function on and off, use the "Temp" membrane button.

When the heating function is switched on, the green "Temp" will light up and the digital display will switch to the next display page. The actual temperature is displayed.

While the plate is heating up, the green "Temp" LED will change colour from green to orange (orange = power supply to heating plate).


### 11.1 Control of medium temperature with Pt 100 temperature sensor

When the unit is switched on, the message "PT2 OFF" is displayed.

The "PROBE" symbol is shown on the digital display to indicate that the Pt100 sensor is operating. The setpoint and actual temperatures displayed are for the medium.

The heating plate temperature is limited by the setting of the safety temperature circuit (SAFE H).

The medium temperature is controlled and regulated by a fuzzy logic controller. This ensures that the medium, sensed by the Pt100 sensor, is heated as quickly as possible without exceeding to the set temperature. The fuzzy logic controller automatically adjusts to the various heating capacities of different vessels and their content. This system ensures optimum temperature control without temperature drift or fluctuations.

### 11.2 Medium temperature control with twin Pt 1000 temperature sensors

When the unit is switched on, the message "PT2 OFF" is displayed.

The "PROBE" symbol is shown on the digital display to indicate that two Pt1000 sensors are operating.

Heat transfer medium and reaction medium (2 vessels)

There is a temperature sensor in the heat transfer medium and a temperature sensor in the reaction medium.

Automatic recognition of the heat transfer and reaction medium sensors takes place.

The setpoint temperature displayed is for the reaction medium temperature. The actual temperature of the reaction medium is shown in large figures to the left of the setpoint temperature.


The heating plate temperature and the heat transfer medium temperature are limited by the heating plate safety temperature (SAFE H) setting. The reaction medium temperature is limited by the medium safety temperature (SAFE P) setting.

The reaction medium temperature is controlled and regulated by a fuzzy logic controller. The controller is optimised by the heat transfer medium temperature. This ensures that the reaction medium is heated as quickly as possible without exceeding the set temperature. The fuzzy logic controller automatically adjusts to various heating capacities of different vessels and their content. This system ensures optimum temperature control without temperature drift or fluctuations.

### 11.3 Medium temperature control without external temperature sensors

This is not possible in operating mode A. Medium temperature control using external temperature sensors is more precise than standard heating plate temperature control.

The operation without Pt 100 is digitally shown by the symbol "PROBE" blinking. The display rated and real temperature refer to the temperature of the heating plate.

### 11.4 Residual heat warning (HOT)

In order to prevent the risk of burns from a hot heating plate, the RET control-visc /C has a residual heat (HOT) warning function.

When the heating function is switched off and the heating plate temperature is still above 50°C, the green "Temp" LED will flash. The digital display will switch backwards and forwards between "Setpoint temperature+ °C" and "Setpoint temperature + HOT" as a warning that the plate is still hot.

When the viscosity trend measurement function is enabled, measurement of the heating plate temperature over 50°C is not possible and the HOT warning on the digital display is not possible. However, the green "Temp" LED will continue to flash to indicate this condition.

If the main ON/OFF switch on the unit is switched off, the display will remain active and the message "—" will flash in the area of the actual value display area. In this state the stirring function can be switched on and the setpoint speed can also be modified.

If the heating plate temperature drops to below 50°C, the unit will automatically switch off. The display will be switched off and the unit electronics will be powered down by a relay. If you have to switch the heating plate off suddenly when the heating plate is still hot, you must switch the main ON/OFF switch on and then off again rapidly.

The residual heat warning function will not operate in the event of a mains power failure or if the power plug is pulled out of the socket.

## 11.5 Changing the temperature unit of measurement from °C to °F

To toggle the unit of measurement for temperature readings from °C to °F and back again, proceed as follows:

- Switch off the unit at the main ON/OFF switch. (The residual heat display should no longer be enabled. See, Residual Heat Warning (HOT))
- Hold the "Mode" and "Temp" membrane buttons pressed down and switch on the main ON/OFF switch. After approx. 5 seconds the display will begin to flash. Release the membrane buttons.
- Temperatures will now be displayed in the unit of measurement selected. Temperatures in the memory will be converted accordingly.

## 11.6 Setting the safety temperature

Setting and/or changing the safety temperature is only possible in the „A“-mode of operation.

The safety circuit is set as follows:

Switch instrument on with mains switch on the left side of the instrument.

As soon as the display shifts to SAFE H (safety temperature limiter of the heatin plate) (SV is blinking), press the foil key "Mode" (SV is blinking no longer), keep this foil key depressed and set the requested safety temperature with the TEMP rotary knob. The set safety temperature (e.g. 100°C) is displayed. After releasing the foil key "Mode", SAFE H OK is displayed for about 2 seconds, and the safety temperature is stored this way.

When the display changes to show SAFE P (safety temperature limit for external sensors) (SV flashes), press the "Mode" membrane button. The SV will stop flashing. Hold down the "Mode" membrane button down and turn the TEMP knob to set the safety limit temperature required. The safety temperature (e.g. 60°C) will be displayed. Release the "Mode" membrane button. After approx. 2 seconds, the message SAFE OK

will be displayed to indicate that the safety temperature has now been stored in the memory.

## 12 Function stirring

The motor speed is set on the rotary knob "Mot", (0 and/or. 50 to 1500 1/min in steps of 10 1/min each). "SV" is blinking.

The function "stirring" is started and stopped with foil key "Mot".

If the function stirring is started, the motor runs smoothly to the set speed, and the LED "Mot" beside the imprint "Mot 1/min" shines green.

The display changes to real speed display.

### 12.1 Intermittent mode

If Intermit mode is active, "On IntErmit" will be displayed when the device is switched on.

The intermittent mode interrupts the stirring cycle at regular intervals. For example, if the stirring function is switched on with the setpoint speed set at 1 rpm the stirring cycle will be interrupted for approx. 10 seconds after a cycle time of 60 seconds. The cycle time and the length of the interruption can be modified via the serial interface (see section 14.1).

**The stirring bar breakaway monitoring function cannot be enabled when the intermittent mode is enabled.**

The intermittent mode can also be used with the special small stirring bars. These can be used when the stirring bar breakaway monitoring function is not enabled and will continue to stir the medium even when the stirring bar breaks away. These stirring bars are picked up once again when the speed drops cyclically.

## 12.2 Stirring bar breakaway monitoring

When the unit is switched on, the stirring bar breakaway monitoring function is displayed (On/OFF SAFE Rpm).

Stirring bar breakaway can occur with any magnetic stirrer. Breakaway can occur when stirring at high speeds or when large stirring bars are being driven at speeds close to the magnetically transmissible maximum torque and there are variations in medium viscosity or controller fluctuations.

On the RET control-visc /C you can select the stirring bar breakaway monitoring function which operates in the setpoint speed range between 300 and 1500 rpm.

When stirring bar breakaway is detected, stirring is stopped momentarily to enable the stirring bar to be picked up. After this, stirring starts again and ramps up to the setpoint speed (< flashes). If a further stirring bar breakaway occurs within the next 3 minutes of stirring at a constant setpoint speed, the stirring bar will be picked up again and stirring will restart but this time at the setpoint speed reduced by 100 rpm. The red "Temp" LED flashes and the buzzer will sound for several seconds. After the fourth stirring bar breakaway – with less than 3 minutes normal stirring – the error message ER 41 will be displayed (the motor will continue to run and the heating is switched off) the red "Temp" LED will flash and the buzzer will sound. When further stirring bar breakaways occur, the speed will be reduced each time in 100 rpm steps until the minimum of 200 rpm is reached. The speed is no longer displayed. In this case you should select new test conditions.

### IMPORTANT:

**You must calibrate and test the stirring bar breakaway monitoring function for your particular configuration taking into account the stirring bars, speed, vessel and medium you are using.**

The reliability of stirring bar breakaway monitoring increases as the transmitted torque increases.

In a glass beaker filled with water, the function will detect the breakaway of stirring bars longer than 50 mm (d=8 mm) in the speed range from 300 to 1500 rpm. Increasing the stirring bar

length, the stirring bar surface area (e.g. TRIKA coated bars) or increasing the speed will increase the torque transmitted.

**The stirring bar breakaway monitoring function will not operate if the stirring bar breaks away while the stirring is still ramping up to the setpoint speed.**

## 12.3 Setting the intermittent mode and the stirring bar breakaway monitoring function

- Switch off the unit at the main ON/OFF switch.
- Hold the "Mot" and the "Mode" membrane buttons pressed down and switch on the main ON/OFF switch. After approx. 5 seconds release the membrane buttons.
- After the display self-test (when all the LCD segments flash), the message On or OFF IntErmitt will be displayed.
- Press the "Mot" membrane button within the next 5 seconds to switch the intermittent mode ON or OFF as required.
- If you disable the intermittent mode by selecting OFF, the message On or OFF Safe Rpm will be displayed.
- To switch the stirring bar monitoring function ON or OFF, press the "Mot" membrane button within the next five seconds.

### IMPORTANT:

**These changes will be stored to memory when the self-test (SAFE H or SAFE P OK is displayed) has been completed.**

## 12.4 Measurement of viscosity trend

This measurement is determined with the help of the control software via the controlling variable given to the motor, and is shown in percent in the digital display. The maximum controlling variable (maximum torque of motor) is 100%.

The instruments are **not** designed for measuring the absolute viscosity. Only the relative change in the product in relation to an initial point determined by the user, is measured and displayed. As the controlling variable is strongly dependent on the real speed, this should not be changed during a measurement.

The intermittent mode must be OFF.

After the speed of the motor and of the magnetic bar in the medium have stabilized to the set rated speed (waiting time about 3 minutes), the viscosity measurement is started by pressing foil key „Visc“.

The real controlling variable is stored as reference value 0 %  $\Delta P$ , and is digitally displayed. Further displayed are the set rated speed, as well as the viscosity change in %.

Depending on the increase or decrease of viscosity, the value is displayed as positive or negative value. A simultaneous display of the percentage and of the real speed is not possible.



By pressing foil key „mode“, it is possible to inquire rated and real temperature, but the measurement of the viscosity trend has priority, and the instrument switches the digital display back to trend display after about 5 seconds. When the stirring function is switched off, the trend display is switched off as well. Analog, the trend display can only be switched on with an activated function „stirring“.

**PLEASE NOTE:**

**While the measurement of the viscosity trend is on, the state of the heating plate cannot be seen on the digital display, even if the heating plate is exceeding 50°C.**

- be measured especially in media which have poor conductivity. The tip of the measuring sensor must be at least 5 mm from the vessel bottom; a distance of 10 mm is ideal.
- With highly viscous reaction media, the distance between the heating plate and the stirring bars must be kept as small as possible (to reduce the risk of stirring bar breakaway).

## 14 Connections and interfaces

## 13 Tips

- Ensure that the medium is well mixed. (High speeds and large stirring bars ensure good mixing). A medium which is well mixed will ensure uniform temperature distribution and better temperature control.
- The temperature sensors (PT100, PT1000) should be immersed to the maximum possible depth. The minimum immersion depth is 20 mm; do not go below this figure.
- When using metal vessels, do not place the temperature sensors on the bottom of the vessel. Placing sensors on the vessel bottom can cause excessively high temperatures to

## 14.1 Serial RS 232 C interface

### 14.1.1 Configuration

The unit can be controlled from an external PC (using the dedicated software labworldsoft) via the RS 232 C serial interface fitted to the unit.

**To increase safety when controlling the RET control-visc/C from a PC, enable the watchdog function for monitoring the continuous serial data flow. (see section 14.1.4)**

- The functions of the interface lines between laboratory instrument and automation system are selected from the specified signals of the EIA-standard RS 232 C, corresponding with DIN 66 020 Part 1. The allotment of the bushing can be taken from illustration.
- For the electrical properties of the interface lines and for the allotment of the signal status, standard RS 232 C, corresponding with DIN 66 259 Part 1 applies.
- Transmission method: Asynchronous signal transmission in start-stop-operation.
- Mode of transmission: Fully Duplex.
- Character format: Character imaging acc. to data format DIN 66 022 for start-stop-operation. 1 start bit; 7 character bits; 1 parity bit [straight (even)]; 1 stop bit
- Transmission speed: 9600 bit/s
- Data flow control: Hardware handshake RTS/CTS
  - RTS: (pin 7) LOW (positive tension)  
/ computer may send
  - RTS: (pin 7) HIGH (negative tension)  
/ computer may not send
  - CTS: (pin 8) LOW (positive tension)  
/ computer recipient
  - CTS: (pin 8) HIGH (negative tension)  
/ computer not recipient
- Access method: Data communication from

laboratory instrument to computer is only possible on demand of the computer

### 14.1.2 Instruction syntax

Here applies the following:

- The instructions are generally sent from the processor (master) to the laboratory instrument (slave).
- The laboratory instrument exclusively sends on demand of the processor. Even error codes cannot be spontaneously communicated from the laboratory instrument to the processor (automatic system)
- The instructions are transmitted in capital letters.
- Instructions and parameters as well as subsequent parameters are separated by at least one blank. (Code: hex 0x20)
- Each individual instruction including parameters and data as well as each reply are terminated with CR LF (Code: hex 0x0D and 0x0A) and have a maximum length of 80 characters.
- The decimal separator in a floating point number is the point. (Code: hex 0x2E)

The above statements largely correspond with the recommendations of the NAMUR-Association. (NAMUR-recommendations for the design of electric plug connections for the analog and digital signal transmission to laboratory - MSR individual units. Rev. 1.1)

### 14.1.3 Overview of the NAMUR-instructions

Abbreviations:

- X,y = numbering parameter (integer number)
- m = variable value, integer
- n = value of variable, floating point number
- X = 1 Pt100 or Pt1000 medium temperature
- X = 2 heating plate temperature
- X = 3 heating plate - safety temperature
- X = 4 speed
- X = 5 viscosity trend
- X = 7 Pt1000 heat transfer medium temperature

X = 8 Pt temperature sensor – safety temperature  
 X = 50 Pt100 temperature sensor offset in K  
 (-3.0 <= n <=+3.0 )  
 X = 52 Pt1000 1st temperature sensor offset in K  
 (-3.0 <= n <=+3.0 )  
 X = 53 Pt1000.2nd temperature sensor offset in K  
 (-3.0 <= n <=+3.0 )

X = 54 Error 5 response time in seconds  
 (180 <= n <= 1200)  
 X = 55 Intermittent mode cycle time in seconds  
 (10 <= n <= 600)  
 X = 56 Intermittent mode, duration of interruption in  
 seconds ( 5 <= n <= 60)

NAMUR instructions	Function	Display (additional)
IN_NAME	Input description name	
IN_PV_X X=1;2;3;4;5; 7; 8;	Reading the real value	
IN_SOFTWARE	Input software ID number, date and version	
IN_SP_X X=1;2;3;4;7;8; 12;42;50;52;53; 54;55;56;	Reading the set rated value	
IN_TYPE	Input laboratory unit ID	
OUT_NAME name	Output description name. (Max. 6 characters, default: IKAMAG)	
OUT_SP_12@n	Setting the WD safety temperature with the echo of the set value	
OUT_SP_42@n	Setting the WD safety speed with the echo of the set value	
OUT_SP_X n X=1;2;4;7;50; 52;53;54;55;56	Setting the rated value to n	
OUT_WD1@m	Watchdog mode 1: When a WD1 event occurs, the heating and stirring functions are shutdown and message ER 2 is displayed. Set the watchdog time to m (20...1500) seconds, with echo of the watchdog time. This instruction starts the watchdog function and must be sent within the set watchdog time.	
OUT_WD2@m	Watchdog mode 2: When a WD2 event occurs, the speed setpoint will be set to the WD safety setpoint speed and the temperature setpoint will be set to the WD safety setpoint temperature. The WD warning is displayed. The WD2 event can be reset with OUT_WD2@0 - resetting also blocks the watchdog function. Set the watchdog time to m (20...1500) seconds, with echo of the watchdog time. This command starts the watchdog function and must be sent within the set watchdog time.	
RESET	Switching off the instrument function.	

RMP_CONT_X X=1;4	Continuation of ramp function (After prior RMP_PAUSE_X).	
RMP_IN_X X=1;4	Reading the real segment number of ramp. With ramp not started: 0	
RMP_INX_y X=1;4	Reading the accumulated value and the ramp segment duration of ramp segment y.	
RMP_LOOP_RESET_X X=1;4	Ending of ramp loop.	
RMP_LOOP_SET_X X=1;4	To work off the ramps in one loop.	
RMP_OUT_X_yn X=1;4 hh:mm:ss	Setting the accumulated value (n) and the ramp segment duration (hh:mm:ss) for ramp segment y.	
RMP_PAUSE_X X=1;4	Stopping the ramp function. Freezing of real rated value and real ramp segment time.	
RMP_RESET_X	Switching off ramp functions and deleting of all set ramp segments	
RMP_START_X X=1;4	Starting the ramp function, beginning with ramp function No.1 (Only possible after prior START_X. After RMP_STOP_X START_X is not necessary.	
RMP_STOP_X X=1;4	Switching off ramp function. Rated value = 0. (Ramp is maintained, that means, ramp can be restarted with RMP_START_X).	
START_X X=1;2;4;5;7	Starting the instrument's (remote) function.	REMOTE
START_X X=16 (old function)	Activation of the "watchdog" function. If no rated value is entered given in the next 60 seconds, the heating is switched off.	
STATUS	<p>Display of status</p> <p>1S: mode of operation A  2S: mode of operation B  3S: mode of operation C  S0: manual operation without fault  S1: Automatic operation Start (without fault)  S2: Automatic operation Stop (without fault)</p> <p>&lt;0: error code: (-1)</p> <p>- 1: error 1  - ... (see table)  -31: error 31  -83: wrong parity  -84: unknown instruction  -85: wrong instruction sequence  -86: invalid rated value  -87: not sufficient storage space</p>	
STOP_X X=1;2;4;5;7	Switching off the instrument - (remote) function. Variables set with OUT_SP_X are maintained. Contains the instruction RMP_STOP.	REMOTE

### 14.1.4 "Watchdog" function, monitoring the serial data flow

The following applies to situations where the watchdog function is enabled (see Namur instructions). If no new transmissions of these commands from the PC take place within the preset watchdog time, the heating and stirring functions will be shutdown according to the watchdog mode selected or will be controlled using the preset setpoints. An operating system crash, a PC power failure or a fault in the connecting cable to the RET control-*visc* /C can cause an interruption in data transmission.

#### Watchdog – Mode 1

If an interruption in data transmission occurs which is longer than the preset watchdog time, the heating and stirring functions will be shutdown and the error message ER 2 will be displayed.

#### Watchdog – Mode 2


If an interruption in data transmission occurs which is longer than the preset watchdog time, the speed setpoint value will be set to the WD safety speed setpoint and the temperature setpoint will be set to the WD safety temperature setpoint. The WD warning message will be displayed.

## 14.2 Connections between RET control-*visc* /C and external devices

The RET control-*visc* /C unit can be connected to other external devices using the following IKA adapters and adapter cables.


### 14.2.1 PC 5.1 Adapter IKA - Control

An adaptor cable is available as an accessory part. It splits the analog and serial signals. The analog output signals are allotted to a 7-polar coupling according to the NAMUR-recommendation; the serial signals to a 9-polar Sub-D-bushing plug (RS 232 C).

### 14.2.2 PC 2.1 Cable

Cable PC 2.1 for connection to the 9-pin socket on the computer.

### 14.2.3 PC 1.5 Cable

Cable PC 1.5 for connection to the 15-pin socket on the computer (25-pin plug).

#### 14.2.4 PC 1.2 Adapter

Adapter PC 1.2 for connecting the 9-pin socket with an 8-way serial interface (25-pin plug).

#### 14.2.5 AK 2.1 Kabel

Cable AK 2.1 for connecting the 7-pin socket of the PC 5.1 adapter to a printer (4 mm banana plug).

#### 14.2.6 AK 2.2 Kabel

Cable AK 2.2 for connecting the 15-pin socket to a printer (4 mm banana plug).

### 14.3 Tips on using labworldsoft

- The safety temperatures (SAFE H, SAFE P) must be set to match the test conditions set on the RETcv /C also when you are using labworldsoft.
- If the setpoint settings in labworldsoft exceed the setpoint value range or the safety temperature set on the RETcv /C, the labworldsoft settings will not be accepted. The speed setpoint must not exceed 1500 rpm.
- If the ramp function to gradually increase the medium setpoint temperature is used, the system variation and the control hysteresis can be greater than they would be if a setpoint temperature is entered directly.
- Heating plate temperature control in critical processes should only be used by experienced users. Critical processes in this case refers to the testing of large volumes under poor control conditions with a single medium temperature setpoint setting.
- When a watchdog mode is enabled, the watchdog time must be set at 20 times the sensing time. This is to prevent process delays, caused by parallel programs (e.g. Excel) or devices (e.g. NPort-Server), from triggering an Error 2 (watchdog time timeout).

## 15 Faults

A unit fault is indicated by the red flashing "Temp" LED (flashing at one second intervals) and the display of an error code number "Er n" on the display (see the error code table). The buzzer also sounds.

To switch off the warning buzzer, press the "Mode" button.

The heating will be powered down (exception: WD – see table).

Following a fault in the heating circuit, the stirring function will continue to operate at the speed set before the fault took place (exception: WD - see table).

To reset a fault, switch off the unit at the main ON/OFF switch for a few seconds. The unit should remain off long enough for the temperature to settle. (The residual heat display must be OFF. See the residual heat display).

When the unit is switched on again, the safety circuit test is run within 20 seconds. After this the other test/monitoring functions are enabled.

**If the fault will not clear, you must run a RESET cycle (see section 9.1 Resetting parameters to factory default settings (RESET)).**

**If a RESET cycle does not clear the fault, take the unit to your technical service centre.**

## 15.1 Troubleshooting. Error messages, causes and solutions

Error code	Description	Cause	Effect	Solution
Er 1	In operating mode A; temperature sensor not plugged in	<ul style="list-style-type: none"> <li>When the unit was switched on in the mode A, the temperature sensor was not plugged in</li> </ul>	The heating remains switched off	<ul style="list-style-type: none"> <li>Switch off the unit, plug in the temperature sensor, switch on the unit</li> <li>Switch over the operating mode</li> </ul>
Er 2	In remote operation (PC), with watchdog funktion 1 enabled: No communication between PC and RET control visc /C	<ul style="list-style-type: none"> <li>PC does not send any data during the watchdog time.</li> <li>The connection/cable to the PC is broken.</li> </ul>	Heating OFF. Motor OFF.	<ul style="list-style-type: none"> <li>Change the watchdog time.</li> <li>Send data from the PC within the watchdog time (OUT_WDx@m).</li> <li>Check the connector cable and plug.</li> </ul>
Er 3	Temperature inside unit is higher than 80°C	<ul style="list-style-type: none"> <li>Build up of heat between the heating plate and the casing.</li> <li>Permitted ambient temperature has been exceeded.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>Switch off the unit. Allow it to cool down and then switch on again.</li> <li>Change test process.</li> <li>Restore ambient temperature to within the permitted maximum.</li> </ul>
Er 4	Difference between setpoint and actual speeds is greater than 300 rpm.	<ul style="list-style-type: none"> <li>Motor seized or overloaded.</li> </ul>	Heating OFF. Motor OFF.	<ul style="list-style-type: none"> <li>Reduce the torque load, e.g. use smaller magnetic bars.</li> <li>Reduce the setpoint speed.</li> </ul>
Er 5	No temperature increase measured by temperature sensor where temperature difference remains the same (after approx. 3 minutes).	<ul style="list-style-type: none"> <li>Sensor not in medium.</li> <li>Volume of heat transfer medium too large.</li> <li>Medium heat transfer rating is too low.</li> <li>Heat temperature rating of the vessel. is too low.</li> <li>With indirect heating the total heat transfer resistance is too high.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>Immerse the sensor in the medium.</li> <li>Reduce medium volume.</li> <li>Use a heat transfer oil with a better heat transfer rating.</li> <li>Use a metal vessel in place of a glass vessel.</li> </ul>
Er 8	The calibration value of the temperature sensor is outside the limit value.	<ul style="list-style-type: none"> <li>Fault in calibration procedure.</li> <li>Value was incorrectly stored to memory.</li> <li>EEPROM switching error.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>Repeat the calibration procedure.</li> </ul>
Er 9	Fault storing the safety circuit setpoint temperatures.	<ul style="list-style-type: none"> <li>EEPROM fault</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>Switch off the unit. Allow it to cool down and then switch on again.</li> </ul>
Er 11	Temperature sensor unplugged.	<ul style="list-style-type: none"> <li>Medium temperature sensor 1 or 2 (Pt100 or PT1000) was unplugged during an operating cycle.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>Only unplug the temperature sensor when the unit is switched off.</li> </ul>

Error code	Description	Cause	Effect	Solution
Er 12	Temperature sensor plugged in.	<ul style="list-style-type: none"> <li>• Medium temperature sensor 1 or 2 (Pt100 or PT1000) was unplugged during an operating cycle.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Only unplug the temperature sensor when the unit is switched off.</li> </ul>
Er 13	SAFE TEMP P (Probe) adjustable safety circuit fault.	<ul style="list-style-type: none"> <li>• Difference between the setpoint and actual values of the adjustable safety circuit for the minimum temperature monitoring.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Switch ON and set the SAFE TEMP P to another value. If this clears the fault when the unit is switched on again, you can return to using the previous value.</li> </ul>
Er 14	Temperature sensor 1, short-circuit	<ul style="list-style-type: none"> <li>• Short-circuit at temperature sensor plug.</li> <li>• Short-circuit in connecting cable or on temperature sensor.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Check the plug.</li> <li>• Substitute the temperature sensor.</li> </ul>
Er 15	Temperature sensor 2, short-circuit	<ul style="list-style-type: none"> <li>• Short-circuit at temperature sensor plug.</li> <li>• Short-circuit in connecting cable or on temperature sensor.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Check the plug.</li> <li>• Substitute the temperature sensor.</li> </ul>
Er 16	Temperature sensor 1 has exceeded the SAFE TEMP P (Probe)	<ul style="list-style-type: none"> <li>• SAFE TEMP P (Probe) is set at a temperature lower than the actual temperature on the temperature sensor 1.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Allow the medium to cool.</li> <li>• Set the SAFE TEMP P (Probe) at a higher temperature</li> </ul>
Er 17	Temperature sensor 2 has exceeded the SAFE TEMP P (Probe).	<ul style="list-style-type: none"> <li>• SAFE TEMP P (Probe) is set at a temperature lower than the actual temperature on the temperature sensor 2.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Allow the medium to cool.</li> <li>• Set the SAFE TEMP P (Probe) to a higher temperature</li> </ul>
Er 21	Fault during heating plate safety test.	<ul style="list-style-type: none"> <li>• Safety relay does not open.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Switch off the unit. Allow it to cool down and then switch on again.</li> </ul>
Er 22	Fault during heating plate safety test.	<ul style="list-style-type: none"> <li>• S_CHECK cannot generate a higher H_S_TEMP.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Switch off the unit. Allow it to cool down and then switch on again.</li> </ul>
Er 23	Fault in SAFE TEMP H (Hotplate) adjustable safety circuit	<ul style="list-style-type: none"> <li>• Difference between the setpoint and actual values of the adjustable safety circuit for the heating plate temperature monitoring</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Switch ON in operating mode A and set the SAFE TEMP P to another value. If this clears the fault when the unit is switched on again, you can return to using the previous value.</li> </ul>

Error code	Description	Cause	Effect	Solution
Er 24	Heating plate temperature has exceeded the SAFE TEMP H (Hotplate)-setting	<ul style="list-style-type: none"> <li>• SAFE TEMP H (Probe) is set a temperature lower than the heating plate actual temperature.</li> <li>• Breakage of the heating plate temperature control sensor</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Allow the heating plate to cool.</li> <li>• Set SAFE TEMP H (Hotplate) to a higher temperature.</li> </ul>
Er 25	Heating switching element monitoring	<ul style="list-style-type: none"> <li>• There is a short-circuit in a switching element (Triac) of the heating control circuit</li> <li>• Safety relay has shutdown the heating circuit.</li> <li>• The heating or the connection is broken</li> <li>• The heating plate safety temperature sensor is broken.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Switch off the unit. Allow it to cool down and then switch on again.</li> </ul>
Er 26	Difference between the internal heating plate control and safety temperature sensors is too large.	<ul style="list-style-type: none"> <li>• Irregular temperature distribution on the heating plate due to point heat dissipation</li> <li>• Fault in the control or safety temperature sensors.</li> </ul>	Heating OFF.	<ul style="list-style-type: none"> <li>• Switching off the unit. Allow in to cool down and then switch on again.</li> <li>• When using metal blocks on the heating plate surface, ensure regular heat dissipation</li> </ul>
Er 31	Fault in the heating switch element (Triac)		Heating OFF.	
Er 41	Stirring bar breakaway	<ul style="list-style-type: none"> <li>• Stirring bar has broken away on four consecutive occasions without any changes to the setpoint value.</li> </ul>	The setpoint speed is reduced to 300 rpm (min.). Heating OFF.	<ul style="list-style-type: none"> <li>• Select a lower speed setpoint.</li> <li>• Kleineres Rührstäbchen wählen</li> </ul>
WD (Watch Dog)	In remote operation (PC) with watchdog function 2 enabled: No communication between PC and RET control visc /C	<ul style="list-style-type: none"> <li>• PC does not send OUT_WDx@m data during the watchdog time</li> <li>• The connection/cable to the PC is broken</li> </ul>	The setpoint temperature is set to the WD safety temperature The setpoint speed is set to the WD safety speed.	<ul style="list-style-type: none"> <li>• Change the watchdog time.</li> <li>• Send data from the PC within the watchdog time (OUT_WD@m).</li> <li>• Check the connector cable and plug.</li> </ul>

## 16 Maintenance and cleaning

The IKAMAG RET control-visc /C is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

**When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.**

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices.

Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

Send in the case of service the equipment back in the original packaging zurück. Storage packing is not sufficient for the back dispatch. Use additionally a suitable transportation packing.

### 16.1 RET control-visc C Technical Enamel


The technical enamel makes the heating plate easier to care for and more resistant to acids and bases. Because of it, however, the heating plate is also more susceptible to extreme fluctuations in temperature and the force of impact. This can result in cracks forming or the coating flaking off.

Make certain that the bottom of the placing vessel is even,

clean and dry. The bottom of the placing vessel must not have any sharp grooves, sides or edges. Remove residues of bases and immediately. We recommend most strongly that you clean the heating plate regularly.

#### 16.1.1 Protective conductor RET control-visc C

On the lower surface of the heating plate is a metallicly bright place. It can be used as a test point for the protective conductor.

## 17 Associated standards and regulations

### Construction in accordance with the following safety standards

EN 61010-1	UL 3101-1
CAN/CSA C22.2 (1010-1)	EN 61010-2-10

### Construction in accordance with the following EMC standards

EN 61326-1

### Associated EU guidelines

EMC-guidelines:	89/336/EWG
Machine guidelines:	73/023/EWG

## 18 Accessories

- Pt 100.50** Pt100 Temperature measuring sensor stainless steel DIN EN 60751 Cl. A
- Pt 100.51** Pt100 Temperature measuring sensor glass-coated DIN EN 60751 Cl. A
- Pt 1000.50** two-way Pt1000 Temperature measuring sensor stainless steel DIN EN 60751 Cl. A
- Pt 1000.51** two-way Pt1000 Temperature measuring sensor class-coated DIN EN 60751 Cl. A
- H 15** Bath top stainless steel 1,5l
- H 16V** Support rod
- H 36** Holding rod
- H 44** Cross sleeve
- H 28** Bath top stainless steel 1l
- H 29** Oil bath top aluminium 1l
- H 30** Oil bath top aluminium 1,5l
- RS 1** set of stirring bars
- RSE** PTFE-stirring bar remover
- PC 5.1** Adapter                      **PC 1.2** Adapter
- PC 1.5** Cabel                        **PC 2.1** Cabel
- AK 2.1** Cabel                        **AK 2.2** Cabel

### labworldsoft

#### Recommended IKA-stirring bars

**TRIKA**-stirring bar PTFE, triangular  
lengths 25, 40 and 55 mm

**IKAFLO**N-stirring bar PTFE

ø 6mm      lengths 10 and 15 mm

ø 8mm      lengths 20, 25, 30, 40, 50 and 60 mm

## 19 Technical data

Design voltage:	<b>VAC</b>	220... 230 ±10%
or	<b>VAC</b>	115±10%
or	<b>VAC</b>	100±10%
Design frequency:	<b>Hz</b>	50/60
<b>Stirring</b>		
Motor:		EC-motor with ball bearings, without brushes
Max. power input of motor:	<b>W</b>	12
Max. power output of motor:	<b>W</b>	5
Speed range:	<b>rpm</b>	0; 50 at 1500
Setpoint speed setting resolution:	<b>rpm</b>	10
Actual speed display resolution:	<b>rpm</b>	1
Speed variation (no load (115/230V) at 1500 1/min, RT 25°C)	<b>rpm</b>	±30
Stirred quantity (H <sub>2</sub> O)	<b>l</b>	20
<b>Heating</b>		
Heating plate material:		Stainless steel
Heating plate material Vers. "C":		Technical enamel white
Ø of the heating plate:	<b>mm</b>	135
Heating power: (115/230V)	<b>W</b>	600 ±10%
(100V)	<b>W</b>	455 ±10%
Temperature range:	<b>°C</b>	RT ...340
Setpoint temperature setting resolution:	<b>K</b>	0,5 <100°C
	<b>K</b>	1 ≥100°C
Actual temperature -		
medium temperature:	<b>K</b>	0,1 <100°C
	<b>K</b>	1 ≥100°C
heat transfer temperature:	<b>K</b>	1
Heating plate adjustable safety circuit:	<b>°C</b>	50 - 350
Medium plate adjustable safety circuit:	<b>°C</b>	50 - 350
Safety circuit setting resolution:	<b>K</b>	1
Temperature sensor in medium:		Pt100 or 2xPt1000 DIN EN 60751 Kl. A
Temperature sensor -	<b>K</b>	<=±(0,15 + 0,002°ITl)
variation (Pt100 oder 2xPt1000, DIN EN 60751 Kl. A		

Temperature measuring range variation:	<b>K</b>	<=±0,3 <100°C
	<b>K</b>	<=±0,5 >100°C
Temperature measuring range + sensor variation:	<b>K</b>	<=±0,45 at 0°C
	<b>K</b>	<=±0,85 at 100°C
Heating plate temperature variation:		
RET control visc	<b>K</b>	±3 without vessel, heating plate centre at 100°C
RET control visc C	<b>K</b>	±6 without vessel, heating plate centre at 100°C
Control hysteresis with temperature sensor: (500ml water, 600ml glass beaker, stirring bar 30, 800 rpm, PT100.50, 50°C)	<b>K</b>	±0,2 (FUZZY-optimised)
Heating plate control hysteresis		
RET control visc	<b>K</b>	±3 without vessel, heating plate centre at 100°C
RET control visc C	<b>K</b>	±6 without vessel, heating plate centre at 100°C
Heating speed for 1 litre water in H15: (25°C at 99°C, 600W)	<b>min</b>	ca. 16
Analog output: 0...10V		
Temperature		1V/100K
Speed		1V/1000 rpm
Viscosity trend		100%=3,5V -100%=1,5V
Analog output voltage variation:	<b>V</b>	±0,05
Permitted ambient temperature:	<b>°C</b>	+5 at +40
Permitted relative humidity:	<b>%</b>	80
Permitted duration of operation:	<b>%</b>	100
Protection class acc. DIN 60529:		IP 42
Fuses:		F1/F2 6,3AT/250V Id-Nr 1408000 F3 1A (Trafo secondary 20Vac) F4 0,75A(Trafo secondary 10Vac)
Contamination level:		2
Überspannungskategorie:		II
Protection class:		1(protective earth)
Operation at a terrestrial altitude:	<b>m</b>	max. 2000 above sea level
dimensions:	<b>mm</b>	160 x 280 x 100
Weight:	<b>kg</b>	2,8

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
**20 List of spare parts RET cv /C**

Item	Designation
1	countersunk screw
5	spring nut
6	nut M7 x 0,75
7	screwed plug
22	operation knob
23	plastic cover
24	plastic cover
25	cap
27	cover
28	frame
35	heating plate
40	BLP power unit
43	EC - motor
44	bushing
45	washer
46	magnet
50	BLP - control
51	foil keyboard
55	washer
60	reference lable
3003	rocker switch
3004	bushing
3005	SUB-D - bushing plug
3006	main plug 120°C
3008	G-fuse
5001	ball bearing
5002	grooved ball bearing

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## **IKA® Works, Inc.**

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

North America

---

### **IKA® Works, Inc.**

2635 North Chase Pkwy SE  
Wilmington, NC 28405-7419  
USA  
TEL. +1 800 733-3037  
TEL. +1 910 452-7059  
FAX +1 910 452-7693  
E-mail: usa@ika.net

## **IKA® Works, (Asia) Sdn Bhd**

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

Asia - Australia

---

### **IKA® Works (Asia) Sdn Bhd**

No. 17 & 19, Jalan PJU 3/50  
Sunway Damansara Technology Park  
47810 Petaling Jaya  
Selangor, Malaysia  
TEL. +60 3 7804-3322  
FAX +60 3 7804-8940  
E-mail: sales@ika.com.my

## **IKA® Japan Y.K.**

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

Japan

---

### **IKA® Japan Y.K.**

293-1 Kobayashi-cho  
Yamato Koriyama Shi  
639-1026 Japan  
TEL. +81 74358-4611  
FAX +81 74358-4612  
E-mail: japan@ika.de

## **IKA® Works Guangzhou**

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

China

---

### **IKA® Works Guangzhou**

173-175 Friendship Road  
Guangzhou Economic & Technological  
Development Zone  
Guangzhou 510730, P.R.CHINA  
TEL. +86 20 8222-6772  
FAX +86 20 8222-6776  
E-mail: sales@ikagz.com.cn