Mastercycler[®] ep *realplex* *

Instruction Manual

* New! Modified Slip Cover
(light transmissive cover for optical openings in heated platen assembly, US patent pending).

eppendorf

Eppendorf Operating Manuals Feedback

We regularly assess our operating manuals. It is important that they meet your needs and expectations. As a user, your opinion is highly valued. Therefore you have the possibility to evaluate the quality of our operating manuals and to tell us your opinion. Please use our online questionnaire on www.eppendorf.com and receive ep-points as a small thankyou for your participation or send your replies to us either by fax or mail to:

Eppendorf AG
 Technical Writing / HH-UT
 D-22331 Hamburg

Fax: ++49 - 40 - 53 80 18 40

Which document do you refer to? Product: _____ Document ID: _____ (found on the second or last page of the manual, e.g. MixMate / B 5353 900.015-01/0506)	
Language of the manual:	
How often do you use this manual? regularly <input type="radio"/> , frequently <input type="radio"/> , rarely <input type="radio"/> , never <input type="radio"/>	
Please check, where applicable: ☺ : good ☺ : average ☹ : unsatisfactory ○ : don't know	
1) Is the manual readily understandable?	☺ ☹ ☹ ○
2) How do you find the structure of the manual?	☺ ☹ ☹ ○
3) How easily can you find information? (by table of contents, index etc.)	☺ ☹ ☹ ○
4) How useful are the figures?	☺ ☹ ☹ ○
5) Would you like to have:	
• more figures	☺ ☹ ☹ ○
• more text	☺ ☹ ☹ ○
• no change	☺ ☹ ☹ ○
6) How happy are you with the	
• print of the manual?	☺ ☹ ☹ ○
• binding of the manual?	☺ ☹ ☹ ○
• color of the manual?	☺ ☹ ☹ ○
7) Remarks and suggestions (e.g. missing or too much content) as well as statements, particularly for critical ratings in 1) to 6).	
Job title:	E-Mail: (for possible further inquiries only)

Certificate

Licensed Real-Time Thermal Cycler

This instrument, Serial No. _____, is an Authorized Thermal Cycler.

Its purchase price includes the up-front fee component of a license under United States Patent Nos. 4,683,195, 4,683,202 and 4,965,188, owned by Roche Molecular Systems, Inc., and under corresponding claims in patents outside the United States, owned by F. Hoffmann-La Roche Ltd., covering the Polymerase Chain Reaction ("PCR") process to practice the PCR process for internal research and development using this instrument. The running royalty component of that license may be purchased from Applied Biosystems or obtained by purchasing Authorized Reagents. This Instrument is also an Licensed Real-Time Thermal Cycler for use with applications licenses available from Applied Biosystems. Its use with Authorized Reagents also provides a limited PCR license in accordance with the label rights accompanying such reagents. Purchase of this product does not itself convey to the purchaser a complete license or right to perform the PCR process. Further information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California, 94404, USA.

No rights are conveyed expressly, by implication or estoppel to any patents on real-time methods, including but not limited to 5' nuclease assays, or to any patent claiming a reagent or kit.

Applied Biosystems does not guarantee the performance of this instrument.

This Real-Time Thermal Cycler or Real-Time Temperature Cycling Instrument is licensed for use solely in research and all other applied fields except human *in vitro* diagnostics under Applera's United States Patent No. 6,814,934 and corresponding claims in non-U.S. counterparts thereof. No right is conveyed expressly, by implication or by estoppel under any other patent claim, such as claims to apparatus, reagents, kits, or methods such as 5' nuclease methods. Further information on purchasing licenses, contact the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California, 94404, USA.

No part of this publication may be reproduced without the prior permission of the copyright owner Copyright© 2005
Eppendorf AG, Hamburg.

Mastercycler® is a registered trademark – US Reg. No. 2,273,536 – Canadian Reg. No. TMA510572

Covered by one or more of US Patent Nos. 6,210,958; 6,767,512 and 7,074,367

Proprietary particle protection: Patent pending.

Table of contents

1	Safety precautions.....	9
1.1	Sicherheitshinweise	11
1.2	Consignes de sécurité	13
1.3	Avvertenze di sicurezza	15
1.4	Instrucciones de seguridad	17
1.5	Avisos de segurança	19
1.6	Veiligheidsvoorschriften	21
1.7	Sikkerhedshenvisninger	23
1.8	Turvallisuusohjeet	25
1.9	Wskazówki bezpieczeństwa	27
2	Installation	29
2.1	Delivery package	29
2.2	Setting up the device	29
2.3	Functional units	30
2.3.1	Mastercycler ep <i>realplex</i>	30
2.4	Startup	32
2.4.1	Connection of <i>realplex</i> module to thermo module	32
2.4.2	Connection of Mastercycler ep <i>realplex</i> to a computer	32
2.4.3	Connection of other components	33
2.4.4	Installation of <i>realplex</i> software	34
2.4.4.1	Main installation procedure	34
2.4.4.2	Restore database	35
2.4.4.3	Start the <i>realplex</i> software	35
3	Operation	36
3.1	Login as administrator	36
3.2	Login	36
3.3	User Change	36
3.4	System configuration	38
3.4.1	User administration	38
3.4.1.1	Setting up new users	39
3.4.1.2	Editing of users	39
3.4.1.3	Deleting users	39
3.4.2	User levels	40
3.4.3	Cycler	40
3.4.4	<i>realplex</i> Module	43
3.4.5	System level	44
3.4.5.1	Properties	44
3.4.5.2	ErrorLog	44
3.4.5.3	UserLog	44
3.4.5.4	Software update	45
3.5	Assay management	46
3.5.1	Creating assays	46
3.5.2	Opening of assays	46
3.5.3	Saving assays	47
3.5.4	Saving templates	47
3.5.5	Working with templates	47
3.5.6	Copying assays	47
3.5.7	Exporting of assays	48
3.5.8	Importing of assays	48
3.5.9	Creating folders	49
3.5.10	Deletion of assays and folders	49
3.6	Probe administration	50
3.7	Dye administration	51
3.8	Calibration	53
3.8.1	Background calibration	53
3.8.2	Color calibration	55
3.8.2.1	Color Calibration using commercially available calibration kits	55
3.8.2.2	Calibration of customer-specific dyes	56

Table of contents

3.9	Database Tool.....	57
3.9.1	Automatically saving of current database	58
3.9.2	Saving of current database.....	59
3.9.3	Importing a database.....	60
3.9.4	Archiving a database	62
4	Programming (Assay Setup)	63
4.1	Creating a new assay	63
4.2	Creation of a plate layout.....	63
4.2.1	Determination of dye	63
4.2.2	Input of sample volume	64
4.2.3	Determination of probe type.....	64
4.2.4	Determination of background.....	64
4.2.5	Definition of sample type.....	64
4.2.5.1	Standards	65
4.2.5.2	Positive Control	66
4.2.5.3	Negative control	66
4.2.5.4	Unknown sample.....	67
4.2.6	Automatic creation of replicates and standard series	68
4.2.7	Sample types for relative quantification	69
4.2.8	Sample types for relative quantification assays	71
4.2.8.1	Set up of monoplex relative quantification assays	71
4.2.8.2	Set up of multiplex relative quantification assays	73
4.2.8.3	Examples	74
4.2.9	Sample types for melting point analysis.....	76
4.2.10	Sample types for endpoint determination	76
4.2.11	Sample types for a +/- assay.....	76
4.2.12	Sample types for gene identification	77
4.2.13	Copying and cutting of samples.....	78
4.2.14	Create Links	78
4.2.15	Deleting samples	78
4.2.16	Grouping of several samples	78
4.2.17	Determination of subassays	79
4.2.18	Editing of a plate layout of a completed assay.....	80
4.2.19	Color for Samples.....	81
4.3	Creation of a PCR program	82
4.3.1	Create a PCR program with a program template	82
4.3.2	Create a PCR program without a program template.....	82
4.3.3	Insertion and editing of program steps.....	83
4.3.3.1	Temp	83
4.3.3.2	Cycle (1 – 5 StepCycle).....	86
4.3.3.3	Hold	86
4.3.3.4	Pause.....	86
4.3.3.5	Sound	86
4.3.3.6	Melting Curve	86
4.3.3.7	Endpoint	87
4.3.4	Copying of program steps.....	87
4.3.5	Deletion of program steps	87
4.3.6	Determination of measuring point	88
4.3.7	Header-settings	88
5	Operation (Monitoring)	90
5.1	Sample loading.....	90
5.2	Block and Lid temperature control	90
5.3	Start assay.....	90
5.3.1	Quick start of an assay	91
5.4	Selection of individual samples	92
5.5	Selection of individual sample types	92
5.6	Cycler status.....	92

Table of contents

5.7	Interruption of an ongoing assay	93
5.8	Continue assay	93
5.9	Abort assay	93
5.10	Use of Mastercycler ep <i>realplex</i> for standard PCR	94
5.11	Use of Mastercycler ep <i>realplex</i> as Fluorimeter	94
5.12	Shut down	94
6	Analysis.....	95
6.1	General	95
6.1.1	Open Assay	95
6.1.1.1	Open assay during an ongoing assay	95
6.1.1.2	Open subassay	96
6.1.2	Selection of individual samples	96
6.1.3	Deactivation of samples	98
6.1.4	Selection of individual sample types	99
6.1.5	Selection of analysis module	100
6.1.6	Zooming charts.....	100
6.1.7	Modifying the scaling of charts	101
6.1.8	Restore default settings.....	102
6.1.9	Copy and save data.....	102
6.1.10	Copy and save charts.....	102
6.1.11	Copy and save results as TXT, CSV, and XLS data format.....	103
6.1.12	Copy and save results to special software tools for relative quantification	104
6.1.12.1	qBase.....	105
6.1.12.2	geNorm.....	105
6.1.12.3	REST 200x.....	105
6.1.12.4	GenEx	105
6.1.13	Compile reports	106
6.1.13.1	Work with report templates	107
6.1.14	Save analyses.....	108
6.1.15	Call up analyses.....	108
6.2	Quantification	108
6.2.1	Determination of threshold	109
6.2.2	Determination of baseline.....	110
6.2.3	Standard curve	111
6.2.4	Analysis.....	112
6.2.5	Multiplex analysis	112
6.3	Relative quantification	113
6.3.1	Performing a relative quantification analysis	113
6.3.2	Results table for monoplex assays	113
6.3.3	Results table for multiplex assays	114
6.4	Melting curve analysis	115
6.4.1	Determination of threshold	116
6.4.2	Analysis.....	117
6.5	Endpoint measurement	118
6.5.1	Standard curve	119
6.5.2	Analysis.....	120
6.6	+/- Assay	121
6.6.1	+/-Assay with endpoint determination as data source.....	121
6.6.1.1	Determination of threshold	121
6.6.2	+/-Assay with PCR data as data source	122
6.6.2.1	Determination of threshold	122
6.6.3	Analysis.....	123
6.6.4	Multiplex analysis	124
6.7	Gene identification.....	125
6.7.1	Determination of threshold	125
6.7.2	Analysis	125
6.8	Inverted data.....	127

Table of contents

7	Maintenance.....	128
7.1	General instructions on cleaning of device	128
7.2	Cleaning of thermoblock	129
7.3	Validation and calibration of thermoblock	129
8	Troubleshooting	130
8.1	General Errors.....	130
8.2	Errors during SYBR Green assays.....	132
8.3	Errors during qRT-PCR	132
8.4	RealSpector	133
9	Short instructions	134
9.1	Login	134
9.2	Setting up new users	134
9.3	Creating assays	135
9.4	Saving assays.....	135
9.5	Performing an assay	136
9.5.1	General	136
9.5.2	Plate Layout.....	136
9.5.3	PCR Program.....	136
9.5.4	Start Assay	137
9.5.5	Interruption of an ongoing assay	137
9.5.6	Continue assay	137
9.5.7	Abort assay.....	137
9.5.8	Analysis.....	137
10	Technical data	138
11	Ordering information	139
12	Ordering information North America	141

1 Safety precautions

- Repairs should only be carried out by an Eppendorf AG authorized specialist. Such authorization can only be obtained by participating in a certified course of training.

- **Transfer**

If the device is passed on to someone else, please include the instruction manual.

- **Disposal**

In case the product is to be disposed of, the relevant legal regulations are to be observed.

- **Information on the disposal of electrical and electronic devices in the European Community**

The disposal of electrical devices is regulated within the European Community by national regulations based on EU Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

According to these regulations, any devices supplied after 13.08.05 in the business-to-business sphere, to which this product is assigned, may no longer be disposed of in municipal or domestic waste. They are marked with the following symbol to indicate this.

As disposal regulations within the EU may vary from country to country, please contact your supplier if necessary.

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2 Installation

2.1 Delivery package

The delivery package includes the following parts:

- 1 mounting set for Mastercycler ep *realplex*
- 1 mains power cable
- 1 PC CAN bus cable
- 1 CAN bus cable
- 1 USB CAN adapter
- 1 Torx screw driver
- 1 Allen head key
- 1 CD with instruction manual and software
- 1 instruction manual in English (printed)
- 1 twin.tec PCR plate 96, skirted / semi-skirted
- 1 clip for rollflex cable
- 1 certificate of conformity

2.2 Setting up the device

To avoid damage caused by condensation, the device should only be switched on 12 hours after setup. This also applies to a change of location.

The Mastercycler ep *realplex* must stand properly on a stable work surface.

Make sure that the Mastercycler ep *realplex* is not placed close to instruments that could cause vibrations (e.g. centrifuges).

There must be enough space available to prevent blocking of the front and rear ventilation slots and to allow cooling air to pass beneath the device.

The distance between the rear ventilation slots and the wall should be min. 10 cm.

Space requirement:	Width:	26.0 cm
	Depth:	41.0 cm
	Height:	39.6 cm

Mains connection:	A safety plug socket for the thermo module.
	The <i>realplex</i> module is supplied via the mains connection of the thermo module.

Make sure that the ventilation slots remain freely accessible and the inlets and outlets are not blocked. There should be no paper or other objects under the device that might block ventilation.

The mains voltage must match the data on the identification plate.

The device must be connected to a grounded socket.

The device should be set up so it can be easily disconnected from the mains.

The delivery carton should be retained to return the device safely in the case of repair.

Note: Transport and storage must not exceed max. 50 °C.
The device should be set up to prevent direct exposure to sunlight.

2 Installation

2.3 Functional units

2.3.1 Mastercycler ep *realplex*


Figure 1: Front view,
realplex module closed

- 1 *realplex* module
- 2 Status display
- 3 Sealing clamp (locked)
- 4 Holding clip for control panel (not necessary for Mastercycler ep *realplex*)
- 5 Thermo module


Figure 2: Front view,
realplex module opened

- 1 Sealing clamp (opened)
- 2 *realplex* module
- 3 Thermoblock

After samples have been loaded onto equipment the *realplex* module is pulled to the front, over the inserted tubes or plate, closed with the clamp, and a program then started. When moving the sealing clamp into position, make sure you do not catch your fingers.

Do not reach beneath the lid when the *realplex* module is open.

2 Installation


Figure 3: Rear view

- 1 *realplex* module
- 2 Connection strip of *realplex* module (see Fig. 4)
- 3 Thermo module
- 4 Ventilation slots
- 5 Mains power switch 0 = off / I = on
- 6 Mains power socket
- 7 Connection strip of thermo module (see Fig. 5)


Figure 4: Connection strip of *realplex* module

- 1 RS232: serial interface
- 2 Address: rotary switch Address
- 3 CAN_out
- 4 CAN_in
- 5 Terminal resistance


Figure 5: Connection strip of thermo module

- 1 Printer: parallel interface
- 2 RS232: serial interface
- 3 Control panel: connection of control panel (not necessary for Mastercycler ep *realplex*)
- 4 CAN_out
- 5 CAN_in
- 6 Terminal resistance
- 7 Address: rotary switch Address


Note: Only devices complying with the standards EN/IEC 60950 and UL 60950A should be connected to the interfaces.

2 Installation

2.4 Startup


2.4.1 Connection of *realplex* module to thermo module

CAN_out

Connect the CAN bus connecting cable to the CAN_out connection of the thermo module.

CAN_in

Connect the free end of the CAN bus connecting cable to the CAN_in connection of the *realplex* module.


Term

ON

Set the Term switch at the rear of the *realplex* module to ON and that of the thermo module to OFF.

Address

Move "Address" on the *realplex* module using the rotary switch address at the rear of the *realplex* module to the setting 11. The arrowheads of both rotary switches are each set to 1. The "Address" on the thermo module is set to 02 (arrowhead of the left-hand rotary switch set to 0 and the arrowhead of the right-hand rotary switch set to 2).

Note: Make sure that the right addresses are allocated so that after switch-on of the devices they are detected automatically by the software and shown in the software display.

2.4.2 Connection of Mastercycler ep *realplex* to a computer

Computer hardware requirements:

AMD Athlon / Intel Pentium III @ 1.0 GHz, min. 512 MB RAM
 4 GB hard-disk storage (up to 1 GB required for software and database)
 USB 1.0/1.1 or USB 2.0 connector
 CD-ROM writer + software for writing CDs (integrated or external)
 CAN USB adapter, firmware 4.2 (standard accessories)


Requirement on operating system of computer:

Windows XP (Home Edition or higher)

Optional for a network environment:

Network card for connection to additional workstation computers
 Installation of a TCP/IP network

CAN_in

Connect the PC CAN bus to the CAN_in connection of the thermo module.

Connect the free end of the PC CAN bus cable to the USB CAN adapter.

Connect the USB CAN adapter to the USB interface of the computer.


Before switching on the device check that the correct mains voltage is available for the device. Compare the data on the identification plate at the left hand side of the thermo module.

Connect the mains power cable to the thermo module and then connect to the mains power supply (see Fig. 3 Rear view). The mains power supply of the *realplex* module is via the connected thermo module.

After the user has switched on the device and started the *realplex* software, the connected Mastercycler ep *realplex* is detected automatically by the software and shown on the status display.

The computer must not be disconnected from the Mastercycler ep *realplex* while running.

2 Installation

2.4.3 Connection of other components

Note: Only devices complying with the standards EN/IEC 60950 and UL 60950A should be connected to the interfaces.

2 Installation

2.4.4 Installation of *realplex* software

2.4.4.1 Main installation procedure

First close all programs that are currently running on your computer. Insert your USB CAN adapter into the USB port. A few seconds later a hardware installation wizard will be displayed.


Close Wizard with **Cancel**.

Note: Do not remove the USB CAN adapter until the installation is complete!

Attention: For correct installation of the *realplex* software computer administration rights are highly recommended!

Now insert the *realplex* CD-Rom, and start **setup.exe** in the CD-Rom directory (if Autorun is active, the setup will start automatically). Then the software will ask you for the components you wish to install:

Continue with **Next**. Several windows will pop up to the screen, accept the default settings with **Next**, **Yes**, **Install** and **Finish**, respectively. When the installation is complete, restart the computer.

Now you will find the new group Eppendorf in the startup menu, from where you can start the software you have installed. Additionally, the software automatically creates icons on the desktop of the computer.

2 Installation

2.4.4.2 Restore database

If the *realplex* software is installed the first time on startup, the software will recognize that no database is available and will automatically start the **Database Tool**. Just accept the default settings given on the sheet **Restore Database**, and press button **Start Restore**. When restore has finished, the *realplex* software can be started by clicking the button **Start realplex** in the upper right hand corner.

If *realplex* software has already been installed earlier, a restore of the database will not be necessary, because the existing database from the earlier installation can be used.

2.4.4.3 Start the *realplex* software

Switch on Mastercycler ep *realplex* and wait until the LED on the top switches to green. Now you are ready to start the *realplex* software.

You also may start the *realplex* software with the cyclor switched off or even with no cyclor hardware at all. In the latter case, an error message will be displayed. Ignore the message and continue.

Attention: Do not use any programs, which require major parts of the computer's capacity (e.g. CD Creator, Media Player) during usage of the *realplex* software!

Never run the computer in battery mode! Change the settings for the power supply of the computer to be sure that it will not switch to the standby mode during an ongoing assay!

3 Operation


3.1 Login as administrator

When the *realplex* software is opened for the first time, a dialogue **Login** (see 3.2) is displayed, prompting the Administrator named Eppendorf to enter his password. This is factory-set and can be changed in the dialogue **User Administration** if required (see section 3.4.1). Whenever the program is opened again, not only the Administrator but also all users he has set up on the system can then log in (see section 3.2).

3.2 Login

To work with the Mastercycler ep *realplex* users have to log in when opening the program. All users known to the system are listed in a selection box. Select the required user and enter the relevant password.

New users are set up by the Administrator (see section 3.4.1.1).

If the wrong password has been entered accidentally, the display shows **Login failed**. After confirmation with the function key **OK**, the Login procedure is repeated.

If the dialogue is ended using the function key **Cancel**, the main screen is displayed (see section 3.3). However, all function keys are inactive. To work in the program, the user has to open the dialogue **Login** again (see section 3.2) and repeat the Login procedure.

Note: The logged-in user now has full read and write access to his folders and assays and can start his assays on the connected Mastercycler ep *realplex*. The folders and assays of other users can be viewed. To start and edit the assays of other users it is necessary to copy them to his own folder.

3.3 User Change

To log in to the system when another user is already logged on, proceed as follows:

Logout is performed by selecting the appropriate icon.

If the icon **Login** is selected, the dialogue **User Login** is displayed. The user can then log in as described in section 3.2.

Alternatively these functions can also be accessed via the menu item **File**.

Only one user can ever be logged on. The login of a new user automatically results in logout of the existing user.

Current User: Beate

The current user is shown in the status display at the bottom of the screen.

The main screen appears after login as the Administrator (or user). The main screen is divided into Navigator, work area and status display. If individual program functions are selected in the navigation tree on the left, they will be displayed in the work area. The status display provides information about the status of the thermo module and the user currently logged on.

3 Operation

3

Operation

The screenshot shows the Mastercycler ep realplex software interface. The interface is divided into several sections:

- 1 Menu bar:** Located at the top, containing 'File', 'Edit', 'Commands', 'Plate Layout', 'Setup', and 'Help'.
- 2 Tool bar:** Located below the menu bar, containing various icons for file operations and assay control.
- 3 Navigator:** A tree view on the left side showing the software structure: 'Assay Setup - <->', 'Plate Layout', 'PCR Program', 'Monitoring', and 'Analysis - <->'. Below it is a 'Samples' table with columns for 'Position' and 'Name'.
- 4 Work area:** The central area displaying the 'Plate Layout' configuration. It includes settings for 'Filter + Dyes' (with filters for 520 nm, 550 nm, 580 nm, and 605 nm), 'Sample Vol.' (in µl), 'Probe', and 'Background' (set to 'twin.tec skirted 20 µl'). Below these settings is a 12-column by 8-row grid (rows A-H, columns 1-12) for defining the assay layout. To the right of the grid are 'Sample Info', 'Comment', and 'Hints' (displaying 'No dyes selected!') fields.
- 5 Status display:** A bar at the bottom showing 'Current User: EPPENDORF', 'Cycler Status: Idle', and 'Database: DATA_20080603.GDB'.

- 1 Menu bar
- 2 Tool bar
- 3 Navigator
- 4 Work area
- 5 Status display


The Navigator can be closed or opened by selection of the icon **Show/Hide Navigator**. Alternatively this function can also be accessed via the menu item **Setup**.

Note: All actions which are available for the respective program functions and system configurations can be accessed via the tool bar or the menu bar. Selected actions are additionally accessible via the right mouse button.

3 Operation


3.4 System configuration

In the menu bar the dialogue **realplex System configuration** can be opened under the menu item **Setup** by selecting **System Configuration**. This lists all users as well as the connected modules on the left. In addition, nodes for system settings and other functions are available. In the navigation tree different levels are indicated by symbols. The subordinate folders become visible if individual nodes are selected with the cursor. To close the nodes, reselect. If levels or subordinate folders are selected, the relevant window is displayed on the right of the dialogue.

3.4.1 User administration

If the node **Users** is selected, the tab sheet **User Administration** is displayed, listing all users. The logged-in Administrator can use this dialogue to set up new users and edit or delete users who are already set up.


The function keys are only accessible for a logged-in Administrator. If a user is logged-in without Administrator rights, these function keys are then inactive.

3 Operation

3.4.1.1 Setting up new users

New users can only be set up by the Administrator when he is logged in. To do so, the dialogue **realplex System configuration** is opened under the menu item **Setup**.

If the node **Users** is selected, the tab sheet **User Administration** can be seen on the right of the dialogue. The function key **New** is used to open another dialogue **New User**.

The name and password of the new user is entered in the relevant fields.

The individual access rights of a user are determined via the user group. There are two preset user groups.

Administrator: has access to all program and administration functions.

Standard User: have access to all program functions.

In addition, it is possible to enter a remark in the field Comment. After confirmation the dialogue is closed and the newly setup user displayed in the dialogue **realplex System configuration** under **Users**. Color-coding indicates the group to which a user belongs: administrators have red jackets and standard users green ones.

Fixed access rights are assigned to the two user groups and cannot be changed.

Note: If a user has forgotten his password, the logged-on Administrator can allocate a new password for the user (see section 3.4.1.2). The Administrator is not able to view user passwords.

3.4.1.2 Editing of users

In the dialogue **realplex System configuration** all data for an existing user can be edited by the Administrator. To do so, the user is first selected in the window **User Administration** and the dialogue **Edit User** then opened using the function key **Edit**. After the data have been changed, this is confirmed with **OK**.


3.4.1.3 Deleting users

Users can only be deleted by the Administrator. To do so, he must log in as described in section 3.1 and open the dialogue **realplex System configuration**. After a user has been selected in the window **User Administration**, he is deleted via the function key **Delete**.

3 Operation

3.4.2 User levels


If individual users are selected under **realplex System configuration**, user-specific information is shown on the right side of the dialogue. Besides the user's name and his user group, this tab sheet also includes the number of folders and programs set up by the user. The time of last access is also shown. A logged-in user can change the comment and password here.

To change the password inputs have to be made, firstly in the field **Password** and secondly in the field **Confirmation**. The entries are confirmed with the function key **OK** and the dialogue ended with **Close**. The user's password cannot be viewed by the Administrator.

3.4.3 Cyclcer

All information available for the connected Cyclcer can be viewed in the Cyclcer node.


3 Operation

The device characteristics are listed under **Properties**. Besides data relating to the firmware the last user to be logged on is specified here, as well as the last program run.

The tab sheets **Incubate** and **Lid** can be used to bring the thermoblock and the heated lid to an adjustable temperature. To do so, enter the required value in the field **New temperature** and confirm by selecting the check box **Incubate**. To stop tempering, deactivate the checkbox **Incubate**.

The documentation of the PCR programs can be viewed in the tab sheet **Protocol**.

The header contains information about the logged-in user who started the PCR program as well as the name and start time of the PCR program. This is followed by data relating to properties and header settings of the program. In addition, information about the individual program steps is displayed. The protocol can not be deleted.

Save

If this function key is selected, the protocol is saved as a TXT file.

Print


This function key can be used to print out the protocol to a connected printer.

3 Operation

3

Operation

Double-clicking on the Cyclor node, the node **System** is displayed with information about the block properties.


- Properties: The tab sheet **Properties** shows the block type and serial number of the thermo module as well as the data for the last validation and calibration of the thermoblock.
- Sensors: Display of current sensor temperatures.
- Temp. Validation: Testing of device temperature using a temperature sensor which is part of the temperature validation system (see Ordering information). The execution is described in the operating manual of the temperature validation system.
- Temp. Calibration: Adjustment of device temperature with the help of a temperature sensor (forms part of temperature validation system). Performance is described in the operating manual of the temperature validation system.

3 Operation


3.4.4 *realplex* Module

System-specific information for the *realplex* module is shown under the node **realplex Module**.

Besides information about the version of the Mastercycler ep *realplex* and the version of the device software, the serial numbers of certain assemblies are shown here. There is also data on the filter equipment and the photomultiplier voltage selected. In addition, **realplex Module** indicates which CAN bus address has to be set at the *realplex* module and whether the connection between the thermo and *realplex* module is via a CAN bus cable. **Cover State** shows the current position of the heated lid.

Double-clicking on the node **realplex Module** displays another node underneath. The directory **Module Log** documents all errors occurring during operation. The table can be saved as a TXT file using the function key **Save** or printed out to a connected printer with **Print**.


3 Operation

3

Operation

3.4.5 System level

The system node provides information about the current software version and includes UserLog and ErrorLog to monitor user activities and any errors.

Alternatively the dialogue **realplex System Info** can also be opened using the icon **Show System Info**.

3.4.5.1 Properties

Here the status displays of the thermo and *realplex* module can be viewed.

3.4.5.2 ErrorLog

ErrorLog documents all error messages. This data can be saved using the function key **Save** or printed out to a connected printer with **Print**. The ErrorLog can not be deleted.


3.4.5.3 UserLog

UserLog documents all actions carried out on the Mastercycler ep *realplex*. The system records all Login / Logout procedures and all actions which involve saving data, on a user-related basis. In addition, all program interactions such as start, pause or abort are documented. The UserLog can be saved or printed out to a connected printer. The UserLog can not be deleted. Whereas the tab sheet **UserLog** covers only the last approximately 5000 entries, all existing actions are listed in the **Archive**.


3 Operation

3.4.5.4 Software update

Double-clicking on the system node displays two other nodes. After selection of one of the nodes a software update can be carried out in the right-hand work area.

If the icon **Browse** has been selected, the dialogue **Load firmware file** is displayed.

After selection of a software version, it is loaded with **Open** in the field **Source Path**. The software update is carried out by selecting the function key **Start**.

Note: Updating the *realplex* module from firmware 2.xx is not possible by the way described here. The **Firmware Update Tool** from the start menu has to be used instead.

3 Operation

3.5 Assay management

The following functions are available in the tool bar to manage assays:

1 2 3 4 5 6

- 1 New Assay
- 2 Open Assay
- 3 Save Assay
- 4 Print Report
- 5 Import Assay from File
- 6 Export Assay to File

Note: Alternatively, these functions are also accessible via the menu item **File** or the right mouse button when the cursor is positioned in the navigation tree.

3.5.1 Creating assays


If the icon **New Assay** is selected, a blank plate layout and a standard PCR program are loaded in the work area.

3.5.2 Opening of assays

After selection of the icon **Open Assay**, the dialogue **Assays** is opened.

If a user or folder is highlighted, all assays stored there are shown on the right. In addition to the name, the date of creation and updating of the assay are specified in the panel **Assays**. Moreover, the table shows the status of the assay.

The assay marked **Setup** has not yet started. On opening it is loaded to **Assay Setup** as well as to **Analysis**. However, no data are available in the latter.

Assays which can be used as a template for the creation of new assays are marked **Template**. See section 3.5.4 and 3.5.5 for further information about the creation and usage of templates.

In the case of assays marked **Complete** data acquisition has already taken place. If this assay is opened, the data are displayed under **Analysis** in the navigation tree. In addition, the assay appears in **Assay Setup** unless another assay is currently running (see section 6.1.1.1).

If individual assays are selected, all analyses of this assay saved to date are shown in the lower window **Analysis**. The date of creation and the name of the user who performed analysis are also listed here.

Assays which have been aborted e.g. due to a power outage or system failure are marked **Aborted**. They may also be analyzed if sufficient data were acquired before the assay was aborted.

3 Operation



3.5.3 Saving assays

A logged on user can save assays under his node or in one of his folders. The editing functions for these programs are only accessible to the logged-in user. Other users are only able to view and copy these assays.

If the icon **Save Assay** is selected, the current assay is saved.

To save an assay under a name or to another folder, the relevant dialogue can be opened by selecting the function **Save as** under **File** in the menu.

After the inputting of the name, Assay is selected using the check box **Assay type** and confirmed with **OK**.

Note: Only the following characters are allowed for assay names:
alphanumeric characters, umlaut $\text{ä}\text{ö}\text{ü}$, blank, and $\text{.,_}\#()\text{=}$.

3.5.4 Saving templates

If an assay should be used more than once after creation, it is saved as a template. To do so, proceed as under **Saving assays**, selecting Template in the combo box **Assay type**.

3.5.5 Working with templates

Templates can be used for the simple creation of new assays. To do so, the template is loaded in **Assay Setup** using the function **Open Assay** and then saved under a new name and Assay as Assay type in the user's node.

3.5.6 Copying assays

A user can copy assays and templates of other users to his own directory in order to start or edit them in his own folders. To do so, **Open Assay** is used with the dialogue **Assays** to select and open the relevant assay in the folder of another user. The user then saves the assay under his own user node via **Save as** (see chapter 3.5.3).

The copied assay can now be processed or started directly on a connected Mastercycler ep *realplex*.


3 Operation

3.5.7 Exporting of assays

The export and import functions are used to exchange assays between several workstations. Here the assays are exported or imported with the plate layout, PCR program and measured data.

The icon **Export assay to file** can be used to save the current assay in ASY format to a computer at another workstation for further processing.


To save an assay the upper combo box is used to select a storage location, followed by entering the file name and selection of **realplex-Exchange (*.asy)** in the combo box **File type**.

3.5.8 Importing of assays

If the icon **Import assay from file** is selected, an assay in ASY format can be loaded to a blank Assay Setup. Therefore, a new assay has to be opened first.

To import an assay, the relevant storage location must be selected in the upper combo box followed by the format **realplex-Exchange (*.asy)** as the File type.

3 Operation

3.5.9 Creating folders

A logged-in user can create folders under his name (node) and file associated assays there so he can find them more easily later on. To create a new folder proceed as follows:


Open the dialogue **Assays** using the icon **Open Assay**.

New

After selecting his own user node, the user adds a new folder with the function key **New**. After selecting a folder, subfolders can be created in the same way.

Rename

The created folders can be renamed by highlighting and then selecting **Rename**. Also assays can be renamed by this function.

Folder and subfolders are arranged alphabetically.

Assays can now be copied or created in the generated folders (see section 3.5.4 and 3.5.2).

3.5.10 Deletion of assays and folders

Open the dialogue **Assays** using the icon **Open Assay**. The user can delete a selected assay from his own user node using the function key **Delete** and confirming with **OK**. The same procedure is used for folders.

Delete

After an assay or folder is selected, it can be deleted using the function key **Delete** and confirming with **Yes**.


3 Operation

3.6 Probe administration

Probe administration allows a probe type to be selected later on in the plate layout of an assay. This is solely for the purpose of documentation and has no influence on the PCR program. The measuring point necessary for data acquisition is determined on creation of the PCR program (see section 4.3.6).


If the icon **Probes** is selected, the relevant dialogue will open.

In this dialogue it is possible to manage the probes used which are available to every user under **Plate Layout**.

The name of the probe type is entered in the dialogue displayed. Confirm with **OK**.

The name of an existing probe can be amended in the dialogue displayed. Confirm with **OK**.

Confirm the deletion of a probe type using **Yes**.

3 Operation

3.7 Dye administration

If the icon **Dyes** is selected, the relevant dialogue is opened.

This dialogue is used for the management of the dyes which have to be specified for fluorescence measurement. Here it is the emission maximum of the respective dye that is important as the user has to select the appropriate filter (e.g. FAM and filter 520 nm (510 nm – 530 nm)). Please consult the documentation of the respective manufacturer for the emission maximum of the dye used.

Note: Before performing color calibration (see section 3.8.2) the relevant dyes have to be defined.

Name	Wavelength	Channel	Sensitivity	Calibrated by	Calibration date
FAM	520 nm	1	2444	Melanie	May/19/2005
SYBR	520 nm	1	3544	Melanie	May/20/2005
NED	580 nm	3	7632	Melanie	May/19/2005
JOE	550 nm	2	6433	Melanie	May/20/2005
ROX	605 nm	4	7951	Melanie	May/20/2005
VIC	550 nm	2	6543	Melanie	May/20/2005
TAMRA	580 nm	3	7732	Melanie	May/20/2005

Buttons: Edit/View..., Delete, New Dye..., Close

Name

Name of dye (max. 12 characters).

Wavelength

Wavelength of emission maximum of dye in nm.

Channel

Channel for detection of the dye. This is automatically preset by the software based on the wavelength specified.

Sensitivity

Sensitivity level of the photo multiplier specified through dye calibration. This consists of a four-digit number code, with a number being specified for each of the four filters (1 = low and 9 = high sensitivity level). This number code has its seeds in the color calibration and is preset by the software. The sequence of numbers corresponds to the emission filters 520, 550, 580, and 605 nm (from left to right).

Calibrated by

User who performed the last dye calibration.

Calibration date

Date of last dye calibration.


3 Operation

3

Operation

New Dye

In the dialogue displayed all data for a new dye can be entered.

This dialogue is used to enter the name and the filter in which the new dye can be detected. The filter must match the emission wavelength of the new dye.

Edit/View...

In the dialogue displayed the comment for an existing dye can be amended. Confirm with **OK**.

Delete

Confirm the deletion of a dye using **Yes**.

Attention: If the filter of the dye is changed, the calibration data of this dye will be lost.

3 Operation

3.8 Calibration

For fluorescence measurement with the Mastercycler ep *realplex* it is necessary to perform a background and a color calibration. Calibration is carried out with the same temperatures as for subsequent measurement.

Note: To ensure that there is no liquid on the sealing cover and no air bubbles in the liquid, the calibration plates have to be centrifuged prior to calibration. Make sure that no liquid has evaporated; otherwise use a new calibration plate. Every well of a PCR plate has to be filled with the same volume of liquid. If PCR tubes are used, every position of the thermoblock has to be loaded. Make sure that the instrument has been warmed up for at least 15 minutes before starting a calibration.

Attention: The software automatically determines the necessary sensitivity of the photomultiplier for every dye and filter. No other inputs are required from the user.

Background calibration must always be performed before color calibration!

3.8.1 Background calibration

If the menu item **Setup** is selected, the dialogue **Background calibration** can be opened, listing all background calibrations stored. This list also includes the date of calibration and the user who performed calibration.

Plate	Default	Sensitivity	Calibration date	Calibration by
twin.tec skirted 10 µl		5466	Apr/29/2005	Beate
twin.tec skirted 20 µl	X	3277	Apr/29/2005	Beate
twin.tec skirted 50 µl		4365	Apr/29/2005	Beate

Plate

Name of plate used (max. 32 characters).

Default

The background calibration, which is preset as the default setting is marked with a cross.

Sensitivity

Sensitivity level of the photomultiplier specified through dye calibration. This consists of a four-digit number code, with a number being specified for each of the four filters (1 = low and 9 = high sensitivity level). This number code is preset by the software.

Calibrated at

Date on which background calibration was performed.

Calibrated by

User who performed background calibration.

New

A new background calibration process is performed in the dialogue displayed.

3 Operation

- After the dialogue **New** is opened, a name is entered for the background plate.
- The block and lid temperature are entered in the field **Incubation** in the dialogue on the right. The default settings are 60 °C for the block and 105 °C for the lid temperature.
- After insertion of the background plate the calibration is performed by the function key **Calibrate**. Confirm with **OK**.
- After heating and the first reading the software will prompt to turn the plate. Follow the software prompt to turn the plate and confirm with **OK**.
- If the background calibration is successfully performed, the plate layout shows 96 green positions of the background plate. The prompt **Calibration successfully finished** can be closed by confirming with **OK**.

If the background calibration failed in some positions, the calibration data can be discarded by selecting Cancel. In this case the former calibration data for this consumable will be used. By selecting Proceed, the background calibration data will be excepted anyhow.

View/Run...

In the dialogue displayed all temperature data for an existing background calibration can be changed, and the calibration process then repeated. The data of the previous background calibration will be overwritten.

Default

If the function **default** is selected, this calibration will be the default setting for all assays. This is shown with a cross in the list of background calibrations.

Delete

Confirm deletion of a background calibration using **Yes**.

3 Operation

One reason for a failed background calibration could be a fluorescence contamination in some wells of the thermoblock. For cleaning of the thermoblock see 7.2.

3.8.2 Color calibration


During color calibration the different sensitivity levels of the 96 light paths are adjusted.

In addition, color calibration is the prerequisite for the performance of multiplex assays, which are based on the detection of more than one dye. The emission spectrum of these fluorescence dyes is generally so wide that the light emitted can be detected in more than one filter. Color calibration results in calibration data which the software can use to calculate this crosstalk when performing multiplex assays.

Note: For color calibration it is generally possible to use all tubes and sealing covers which are suitable for real-time PCR. However, it should be remembered that background calibration which is used during color calibration must also be performed with the same consumables.

3.8.2.1 Color Calibration using commercially available calibration kits

Note: Only calibration kits in 96-well format are suitable for color calibration of the Mastercycler ep *realplex*. These kits must include one plate per dye (e. g. ABI Prism 7000 SDS 96-well Spectral Calibration Kit). If a commercially available calibration kit is used for color calibration, background calibration must be carried out using the background plate included in the kit. Selection of the dyes required for calibration depends on the filters contained in the Mastercycler ep *realplex*.

- Open the dialogue **Color Calibration** in the menu under **Setup**.
- After the dialogue **Color Calibration** has been opened, the combo box is used to select the background plate included in the kit. The name of the dye is also selected here.
- The block and lid temperature are entered in the field **Incubation** in the dialogue on the right. The default settings are 60 °C for the block and 105 °C for the lid temperature.
- After insertion of the color calibration plate the calibration is performed by the function key **Calibrate**. Confirm with **OK**.
- After heating and the first reading the software will prompt to turn the plate. Turn the plate and confirm with **OK**.
- If the color calibration is successfully performed, the plate layout shows 96 green positions of the color calibration plate. The prompt **Calibration successfully finished** can be closed by confirming with **OK**.

If the color calibration failed in some positions, the calibration data can be discarded by selecting **Cancel**. In this case the former calibration data for this dye will be used. Selecting **Proceed** the color calibration data will be excepted anyhow.

This procedure is repeated for every dye.

After calibration, a sensitivity code is displayed in the panel **Calibration**. This code consists of a numerical combination, with every digit specifying the sensitivity level of the photomultiplier used to measure the respective dye in the relevant filter. The sensitivity levels for the dyes of the calibration kit are automatically preset by the software.

3 Operation

3

Operation

3.8.2.2 Calibration of customer-specific dyes

For each dye a concentration has to be prepared according to the recommended concentrations shown in the table below. For each concentration a volume of 50 µl has to be pipetted in each well of a 96-well plate. After sealing the plate a color calibration has to be performed for each dye concentration (see chapter 3.8.2.1).

Note: Use the same consumables during background and color calibration. A background calibration has to be done first.

Main filter (nm), which is used during calibration of customer-specific dye	Recommended dye concentration during color calibration in nM
520	200
550	300 – 500
580	1000 – 2000
605	2000

Note: The emission maximum of the dye must be within the main filter.

The sensitivity code should preferably not contain 1 in one of the positions. The code for the main filter should be between 2 and 8. If the sensitivity code contains 9 in the main filter, it is highly recommended to repeat the color calibration with a higher dye concentration.

3 Operation

3.9 Database Tool

All data, analyses, settings and user structures are automatically stored in a database. It is recommended to create backups of the database at regular intervals.


The software **Database Tool** can be used for the management of databases.

Note: Database and instrument are correlated by the serial number of the *realplex* module. If a database is connected to an improper instrument, the instrument has to be calibrated first.

Double-clicking on the icon **Database Tool** displays the dialogue **realplex Database Administration**.

The upper panel **Selected Database** shows the current database and its status.

The tab sheet **Databases** shows a list of all databases and archives which are available. In parentheses, the size of the database is shown. There are two different types of databases:

- The actually used database in which all data are currently saved. The name of the database is composed of the term **DATA** and the date when this database was started. In the panel **Selected Database** this database is marked in green and as open.
- The archived databases have a read-only status and include previous data. They are named with **Archive** and the date of creation. In the panel **Selected Database** this database is marked in red and as closed. Archives can be deleted by the function key **Delete Archive**.

If one of the listed databases is selected it will be shown in the panel Selected Databases and will be connected to the *realplex* software.


If the function key **Start realplex** is selected the *realplex* software will be opened directly.

3 Operation

3.9.1 Automatically saving of current database

After opening the *realplex* software a backup is saved automatically on a regular basis to the directory C:\Program Files\Eppendorf\realplex\Security. The settings for the automatic backup can be accessed in the tab sheet **Configuration**.

Note: Only users with administrator rights have access to the configurations (see chapter 3.9.4).


Note: If more than 10 automatic backups have been created, the oldest one will be deleted automatically. Therefore, creation of customer specific backups on a regular basis is recommended. These manual backups should be stored on an external disc or server.

3 Operation

3.9.2 Saving of current database

Saved copies can be made of the current database using the tab sheet **Backup Database**. This backup contains all data and settings including assays, templates, calibration data and user structure.


Note: It is strongly recommended to save the current database before installing a new *realplex* software version!

The panel shows the database of which a backup is to be made and which is selected in the sheet **Databases**. In general, it is possible to create backups of databases and archives.

The name of the backup is created automatically by the software. Depending on the source the name is composed of backup_DATA or backup_ARCHIVE and the date of creation.

... If the icon **Browse** is selected, the dialogue **Set backup file name** is displayed.

The current database is given a name and loaded to the panel **Backup** in GBK format using the function key **Save**.

Note: For faultless recognition of the database by the *realplex* software it is necessary to precede each name with backup_DATA or backup_ARCHIVE depending on the source!


3 Operation

Start Backup

If the function key **Start Backup** is selected, a backup is created in the directory C:\Program Files\Eppendorf\realplex\Backup. The process is documented in the tab sheet **Log**.


The backup is saved automatically to the directory C:\Program Files\Eppendorf\realplex\Backup. The path can be changed in the tab sheet **Configuration**. If the computer is part of a network it is recommended to select a directory on a server.

Note: Only users with administrator rights have access to the configurations (see chapter 3.9.4).

3.9.3 Importing a database


To use a backup of an older database as the active database in the *realplex* software, the current database must be replaced by a backup. To do so, first select the backup in the tab sheet **Restore Database**.

3 Operation

If the icon **Browse** is selected, the dialogue **Open backup file** is displayed.

A database in GBK format, is loaded to the tab sheet **Restore Database** using the function key **Open**.

If this function key **Start Restore** is selected, the database is loaded from the directory C:\Program Files\realplex\Backup.

The process is documented in the tab sheet **Log**.

Attention: All data including assays, settings, calibration data and user structure are overwritten by this database! It is recommended to create a backup of the currently used database before!
For using a restored database as actually used main database make sure that all calibration data are still valid.

Note: It is also possible to make a restore of the database with a backup of an archive. After this backup has been loaded to the tab sheet **Restore Database** the check box **Create Main Database from Archive** can be activated and the restore can be performed. It is recommended to create a backup of the currently used database before!


All users with Administrator rights have access rights to database management. After selection of the user name in the combo box, the password is then entered and confirmed with **OK** (see section 3.9.4). The password is the same as for the *realplex* software.

3 Operation


3.9.4 Archiving a database

If the database exceeds approximately 250 MB the user will be asked to archive the database after opening the *realplex* software. After archiving the user can continue with a new database including all settings concerning templates, calibration data and user structure.

All users with Administrator rights have access rights for this part of the database management. After selection of the user name in the combo box, the password is then entered and confirmed with **OK**. The password is the same as for the *realplex* software.

After login the tab sheet **Create Archive** is accessible.

The panel shows the database of which an archive is to be made and which is selected in the sheet **Databases**.

The name of the archive is created automatically by the software and it is composed of backup_ARCHIVE and the date of creation.

The backup of the archive is saved automatically to the directory C:\Program Files\Eppendorf\realplex\Backup. The path can be changed in the tab sheet **Configuration**.

In the panel **New Database** the name of the new active database is displayed. The name is created automatically by the software and has the date of creation.

Start Archiving

If the function key **Start Archiving** is selected, an archive is created in the directory C:\Program Files\Eppendorf\realplex\Backup.

The process is documented in the tab sheet **Log**.

4 Programming (Assay Setup)

The software of the Mastercycler ep *realplex* basically consists of two modules **Assay Setup** and **Analysis**. The module **Assay Setup** is divided into the sub-items **Plate Layout**, **PCR Program** and **Monitoring**. The first two sub-items are used to edit assays, while under **Monitoring** an ongoing assay can be monitored in real-time. The module **Analysis** is described in detail in section 6.


Note: Parallel operation of both modules with two different assays is possible. While one assay is running, the analysis of a completed assay can be performed (section 6).

4.1 Creating a new assay

A new assay can be created under the menu item **File** with **New Assay**. This is first displayed without a name in the window on the left of the screen in the navigation tree with the subordinate items **Plate Layout**, **PCR Program** and **Monitoring**.

They are displayed in the work area on the right if individual sub-items are selected.

Alternatively, the different modules can also be loaded to the work area using the functions of the menu item **Commands**.

A blank plate layout is available in the work area as the default setting after opening of the software.

4.2 Creation of a plate layout

To create a plate layout **Plate Layout** is selected in the navigation tree on the left. On the right the 96 positions of a PCR plate can be seen. To edit this, the dyes to be detected have to be specified first.

To aid the user, instructions on further procedure are shown in the field **Hints** on the bottom right next to the work area.

To create plate layouts, the following functions are available in the toolbar:


- | | |
|--------------------|--|
| 1 Unknown | 7 Delete Samples |
| 2 Standard | 8 Group as Replicates |
| 3 Positive Control | 9 Ungroup Replicates (inactive) |
| 4 Negative Control | 10 Create Subassay |
| 5 Create Links | 11 Delete Subassay (inactive) |
| 6 Clear Plate | 12 Edit mode for Plate Layout (inactive) |
| | 13 Color for Samples |

Note: Alternatively, these functions can also be accessed via the menu item **Plate Layout** or the right mouse button.

4.2.1 Determination of dye

Here the dye is selected in the panel **Filter + Dyes** with the relevant filter using the combo box. If the required dye is not included in the selection, it must be first added to the list **Dyes** under the menu item **Setup** (see 3.6). Calibration of the dye must be performed (see 3.8.2).

If the *realplex* software is connected to a Mastercycler ep *realplex*² filter 580 nm and filter 605 nm are not shown in the panel.

Note: The Mastercycler ep *realplex* does not need an internal reference dye to give accurate results.

4 Programming (Assay Setup)

4

4.2.2 Input of sample volume

In the top right panel of the work area the volume of the samples can be entered under **Sample Vol.** Although this input is not necessary, it is recommended for full documentation of the assay.

4.2.3 Determination of probe type

In the panel on the top right in the Plate Layout the probe type can be selected using a combo box. If the required probe type is not included in the selection, it must first be added in the list **Probes** under **Setup** (see 3.7).

Note: Definition of the probe type is not absolutely necessary as it has no influence on the measuring point in the PCR program. This input is however recommended for full documentation of the assay.

4.2.4 Determination of background

In the panel on the top right in Plate Layout the relevant plate type for the background is selected using a combo box. This has to be the plate type in which the assay is also to be performed. If the required plate type is not included in the selection, background calibration first has to be performed with this consumable (see 3.8.1).

4.2.5 Definition of sample type

The relevant dialogues are opened for editing by highlighting one or more positions in Plate Layout and subsequent selection of an icon in the toolbar for the different sample types. All functions used to edit the plate layout are also accessible via the right mouse button or the menu item **Plate Layout**.

After the inputs have been made, the samples appear in the left-hand window with data relating to the type, name of the sample and position in the table **Sample**. Data for a sample selected in Plate Layout are listed in the right-hand window **Sample Info**.

Sample Info	
Position	A1
Sample Name	Sample 1
Type	Unknown
Target 1	Actin
Dye	SYBR
Target 4	Reference Dye
Dye	ROX

4 Programming (Assay Setup)

4.2.5.1 Standards

If the icon **Standard** is selected, the relevant dialogue is opened.

The 'New Standard' dialog box contains the following elements:

- Position:** A1 (Callout 1)
- Name:** A text input field (Callout 2)
- Target gene(s):** A table with columns: Target, Dye, Amount, Unit.
 - Target 1: Checked, FAM, 1, Copies (Callout 3)
 - Target 2: Unchecked, empty, empty, empty (Callout 4)
 - Target 3: Unchecked, empty, empty, empty (Callout 5)
 - Target 4: Unchecked, empty, empty, empty (Callout 6)
- Use of target gene(s) as:**
 - Internal Positive Control (IPC): Unchecked, with a text input field (Callout 8)
 - Dye selection: A dropdown menu (Callout 9)
 - Symbol for sample type: A dropdown menu (Callout 10)
- Buttons:** OK, Cancel, Auto Series -->

- | | |
|--|---|
| 1 Position of sample in the plate layout | 6 Unit for standard concentration |
| 2 Name of sample | 7 Concentration of standard |
| 3 Check box for selection of target gene | 8 Check box for selection of an internal positive control (IPC) |
| 4 Name of target gene | 9 Combo box for selection of dye for IPC |
| 5 Dye used to detect target gene | 10 Symbol for sample type |

Here data are entered for the name of the standard and the target gene to be examined. It is also necessary to specify the amount or concentration of the standard in the relevant fields. Exponential values are entered for example as 1E 8. The unit **copies** and amount = 1 are preset as the default setting. Standard series can be created with the function key **Auto Series** (see 4.2.6).

Note: The number of target genes which can be edited in this dialogue depends on the number of dyes selected in Plate Layout.

4 Programming (Assay Setup)

4

4.2.5.2 Positive Control

If the icon **Positive Control** is selected, the relevant dialogue is opened.

Target 1	Dye 1	Amount 1	Unit 1	Use as
<input checked="" type="checkbox"/>	FAM	0	Copies	<input type="checkbox"/> Allelic Ctrl.
<input type="checkbox"/>				<input type="checkbox"/> Allelic Ctrl.
<input type="checkbox"/>				<input type="checkbox"/> Allelic Ctrl.
<input type="checkbox"/>				<input type="checkbox"/> Allelic Ctrl.

Besides data for the name of the control and the target gene(s), it is possible here to define the positive controls as alleles for the analysis module **Gene identification**. To do so, it is necessary to edit a positive control as allele (see 4.2.12). Each allele is automatically assigned a color in this dialogue for subsequent evaluation under **Analysis**.

It is also possible to enter values in the field **Amount**. However this is not absolutely necessary for subsequent analysis.

4.2.5.3 Negative control

If the icon **Negative Control** is selected, the relevant dialogue is opened.

Only data for the name of the negative control and the target gene(s) are entered here.

4 Programming (Assay Setup)

4.2.5.4 Unknown sample

If the icon **Unknown** is selected, the relevant dialogue is opened.

Besides data for the name of the unknown sample and target gene(s), replicates can also be edited here using the function key **Auto Series** (see 4.2.6).

At the bottom this dialogue includes the panels **Use of target gene(s) as** and **Relative quantification / gene expression**:

Use of target gene(s) as: It is also possible to set target genes either as an internal positive control or as a housekeeping gene by selecting the relevant function and then specifying the respective dye using a combo box. The target gene for which an internal positive control or the housekeeping gene has been determined for analysis is then highlighted in yellow. You can only use one dye for every assay / subsassay.

4 Programming (Assay Setup)

4

4.2.6 Automatic creation of replicates and standard series

With the sample types **Standard** and **Unknown** it is possible to create replicates automatically. When editing standards it is also possible to create standard series automatically. To do so, select the area (at least two positions) in the plate layout in which standard series and/or replicates are to be created.

Auto Series -->

To simplify the creation of standard series and replicates using the dialogue **New Standard**, it has been extended to include the option **Auto Series**.

Note: This function is only activated when more than one position is selected in the plate layout. The positions have to be vacant and contiguous.

The extended dialogue indicates how many positions in the plate layout have been selected and is divided into the areas **Replicates** and **Dilution**.

Replicates: Here it is possible to enter the number of replicates in the spin edit. In addition, this item can be used to define the pattern forming the basis for arrangement of the replicates on the plate. The related replicates can be arranged next to each other either in horizontal or vertical alignment. The direction in which different replicates follow each other can be in rows or columns. The pattern is displayed in the window **Preview** for the purpose of clarification.

Dilution: The spin edits are used to specify the dilution factor and the number of dilution steps. With the combo box **Order** it can also be specified whether the dilution series is to be ascending or descending. In the field **Preview calculated amounts** a dye can be selected via a combo box and the final concentration displayed according to the initial concentration defined for the target gene.

Note: This area is inactive in the extended dialogue **Unknowns**.

After confirmation of the dialogue with **OK** the same number of groups are created in the selected area as the number of parameters. When creating replicates with unknown samples the entire area selected is filled.

Sample names and all other associated data are allocated automatically. With standard series the concentration data for the individual standards are also displayed automatically in the plate layout.

4 Programming (Assay Setup)

4.2.7 Sample types for quantification

To quantify the initial amount of nucleic acid in a sample a standard curve is created from samples with a known initial amount of nucleic acid. Standards are entered into the plate layout followed by unknown samples and negative controls.

Example of programming a standard series:

The aim is to create a standard series for 10 – 100000 copies of a target gene in triplicates, with the triplicates being arranged next to each other in rows. The highest number of copies is to be at positions A1 – C1, and the dilution series is to be in descending order in steps of 10 in horizontal alignment.

- Select area A1 – E3 in the plate layout.
- Open the dialogue **Standard** and expand the dialogue by selecting the function key **Auto Series**.
- Input parameters.
- Confirm input with **OK**.

Result in plate layout:

	1	2	3	4
A	GAPDH... 3: 100000	GAPDH... 3: 100000	GAPDH... 3: 100000	
B	GAPDH... 3: 10000	GAPDH... 3: 10000	GAPDH... 3: 10000	
C	GAPDH... 3: 1000	GAPDH... 3: 1000	GAPDH... 3: 1000	
D	GAPDH... 3: 100	GAPDH... 3: 100	GAPDH... 3: 100	
E	GAPDH... 3: 10	GAPDH... 3: 10	GAPDH... 3: 10	
F				

4 Programming (Assay Setup)

4

Programming (Assay Setup)

Example of programming replicates with unknown samples:

The aim is to define the target gene in triplicates for an unknown sample, with the triplicates being arranged next to each other at positions A1 – C1.

- Select area A1 – C1 in the plate layout.
- Open dialogue **Unknown** including additional item **Auto Series**.
- Input parameters.
- Confirm with **OK**.

New Unknown

Position: **A1**

Name: Sample 001_#/N/#R

Target gene(s)

Target	Dye	Amount	Unit
<input checked="" type="checkbox"/> Bcr2	FAM		
<input type="checkbox"/>			
<input type="checkbox"/>			
<input type="checkbox"/>			

Use of target gene(s) as

Internal Positive Control (IPC)

Housekeeping Gene

Relative quantification / gene expression

Set sample as calibrator

Well with housekeeping gene (Monoplexing): <none>

Automatic Series

Number of selected rows: 3
Number of selected columns: 1

Replicates: Number of Replicates: 3

Define Replicate Pattern:
Alignment: vertical
Direction: by Rows

Preview:

Buttons: ? | OK | Cancel

Result in plate layout:

	1	2
A	Sample ...	
B	Sample ...	
C	Sample ...	
D		

4 Programming (Assay Setup)

4.2.8 Sample types for relative quantification

Relative quantification is used for comparison of the expression levels of one or more target genes between two samples (e.g. from a cell culture subjected to different treatments or agents). The expression level(s) of the target gene(s) are normalized to the expression level of a non-regulated housekeeping gene (reference gene, endogenous control) in both samples. The normalized expression levels of one sample are then related to those of the other sample (designated as **Calibrator**) by calculating the relative expression levels according to the $\Delta\Delta C_T$ method (see chapter 6.3).

For relative quantification with the $\Delta\Delta C_T$ method no external standards are used. Therefore all samples are created as sample type **Unknown**. For each target gene a complete set of samples has to be set up according to the scheme shown below. Otherwise the Relative Quantification analysis module will not be available.

PCR Reaction type	Description	Setup
Monoplex	The target gene and the housekeeping gene are detected in separate reactions. All reactions must contain the same dye (e.g. SYBR Green or sequence specific probes for the target gene and housekeeping gene all containing FAM as reporter dye).	1 target gene reaction or replicate group linked to 1 separate housekeeping gene reaction or replicate group and 1 target gene reaction or replicate group linked to 1 separate housekeeping gene reaction or replicate group both defined as Calibrator .
Multiplex	The target gene(s) and its corresponding housekeeping gene are detected in the same reaction. Sequence specific probes containing different dyes for target gene(s) and housekeeping gene have to be used (e.g. in duplex reactions FAM for the target gene and JOE for the housekeeping gene).	1 reaction or replicate group for target gene and housekeeping gene and 1 reaction or replicate group for target gene and housekeeping gene defined as Calibrator .

Note: For each calibrator to be used, a separate subassay containing a full setup of target gene(s), housekeeping gene and calibrator has to be defined (see chapter 4.2.16).

4.2.8.1 Set up of monoplex relative quantification assays

In **Assay Setup / Plate Layout** select the dye, background and probe type. Enter the reaction volume.

Select the position(s) to use for target gene reactions in the plate layout and open the **Unknown** dialog. Enter the sample and target gene names. Make sure that the checkbox for the corresponding dye is selected. Confirm with **OK**. If applicable, group the samples as replicates.

4 Programming (Assay Setup)

4

Programming (Assay Setup)

Select the position(s) to use for housekeeping gene reactions in the plate layout and open the **Unknown** dialog. Enter the sample and housekeeping gene names. Make sure that the checkbox for the corresponding dye is selected. Activate the checkbox **Housekeeping Gene**. Confirm with **OK** and group the samples as replicates.


Repeat steps 2 and 3 for the calibrator samples.

Open the **Create Links** dialog by clicking the button . In the mini plate layout select the positions containing the target gene reactions and while holding the [↑ Shift] key select the position containing the corresponding housekeeping gene reactions. Click the **Create Link** button. Select the positions containing the target gene and the positions containing the housekeeping gene of the calibrator sample. Click the **Create Links** button.

Note: If a position is selected in the Create Links dialog, all positions currently linked to that position are highlighted automatically.

4 Programming (Assay Setup)

Select the calibrator group in the mini plate layout and click the **Set Calibrator** button. Confirm and close the dialog by pressing **OK**.

Calibrator positions will be marked with a light blue bar.


Note: If more than one replicate group for the target gene and / or housekeeping gene is present in the calibrator set, a confirmation dialog concerning the Relative Quantification calculation will be shown. In this case, the arithmetic average C_T values of all target genes and / or housekeeping genes in the calibrator will be used in the calculation of the relative expression levels. To accept and proceed click the **OK** button.

Arrange all corresponding samples in a subassay before adding additional samples to the plate layout.

4.2.8.2 Set up of multiplex relative quantification assays


In **Assay Setup / Plate Layout** select the dyes, background and probe type. Enter the reaction volume.


Select the position(s) to use for the combined target gene(s) and housekeeping gene reactions in the plate layout and open the **Unknown** dialog. Enter the sample name and the name(s) of the target gene(s). Make sure that the checkbox for the corresponding dyes is selected. Activate the checkbox **Housekeeping Gene** and select the dye to use for the detection of the housekeeping gene in the combo box. Confirm with **OK**. If applicable, group the samples as replicates.

Repeat step 2 for the calibrator samples.

4 Programming (Assay Setup)

Open the **Create Links** dialog by clicking the button . In the mini plate layout select the positions containing the calibrator sample. Click the **Set Calibrator** button. Confirm and close the dialog by pressing **OK**.

Calibrator positions will be marked with a light blue bar.

Note: If more than one replicate group for the target gene and / or housekeeping gene is present in the calibrator set, a confirmation dialog concerning the Relative Quantification calculation will be shown. In this case, the arithmetic average C_T values of all target genes and / or housekeeping genes in the calibrator will be used in the calculation of the relative expression levels. To accept and proceed click the **OK** button.

Arrange all corresponding samples in a subassay before adding additional samples to the plate layout.

4.2.8.3 Examples

Example 1: Expression of 1 target gene is compared between samples from 3 different experiments and a control sample. The control sample is used as calibrator. Triplicate monoplex reactions employing SYBR Green are used.

Plate Layout <->

Filter + Dyes
 Filter 520 nm SYBR Filter 550 nm Filter 580 nm Filter 605 nm

Sample Vol. 25 μ l
 Probe SybrGreen
 Background twintec ss HSF

	1	2	3	4	5	6	7	8	9	10	11	12
A	Exp1	Exp1	Exp1	Exp2	Exp2	Exp2	Exp3	Exp3	Exp3	Ctrl	Ctrl	Ctrl
B	Exp1 1: HKG	Exp1 1: HKG	Exp1 1: HKG	Exp2 1: HKG	Exp2 1: HKG	Exp2 1: HKG	Exp3 1: HKG	Exp3 1: HKG	Exp3 1: HKG	Ctrl 1: HKG	Ctrl 1: HKG	Ctrl 1: HKG
C	samples for experiment 1			samples for experiment 2			samples for experiment 3			calibrator samples		
D												
E												

Sample Info

Position **B10**
 Sample Name **Ctrl**
 Type **Unknown Calibrator**
 Target 1 **Hkg1**
 Dye **SYBR HKG**

4 Programming (Assay Setup)

Example 2a: Expression of 3 target genes is compared between an experiment sample and a control sample. The control sample is used as calibrator. Triplicate monoplex reactions employing SYBR Green are used.

Plate Layout <->

Filter + Dyes

Filter 520 nm SYBR Filter 550 nm

Filter 580 nm Filter 605 nm

Sample Vol. 25 μ l

Probe SybrGreen

Background twintec ss HSF

	1	2	3	4	5	6	7	8	9	10	11	12
A	C_Trgt1	C_Trgt1	C_Trgt1		C_Trgt2	C_Trgt2	C_Trgt2		C_Trgt3	C_Trgt3	C_Trgt3	
B	C_Trgt1 1: HKG	C_Trgt1 1: HKG	C_Trgt1 1: HKG		C_Trgt2 1: HKG	C_Trgt2 1: HKG	C_Trgt2 1: HKG		C_Trgt3 1: HKG	C_Trgt3 1: HKG	C_Trgt3 1: HKG	
C	C_Ctrl1	C_Ctrl1	C_Ctrl1		C_Ctrl2	C_Ctrl2	C_Ctrl2		C_Ctrl3	C_Ctrl3	C_Ctrl3	
D	C_Ctrl1 1: HKG	C_Ctrl1 1: HKG	C_Ctrl1 1: HKG		C_Ctrl2 1: HKG	C_Ctrl2 1: HKG	C_Ctrl2 1: HKG		C_Ctrl3 1: HKG	C_Ctrl3 1: HKG	C_Ctrl3 1: HKG	
E	subassay for target gene 1				subassay for target gene 2				subassay for target gene 3			

Sample Info

Position **A9**

Sample Name **C_Trgt3**

Type **Unknown**

Target HKG **B9**

Target 1 **Trgt3**

Dye **SYBR**

Comment

4 Programming (Assay Setup)

4

Example 2b: Expression of 3 target genes is compared between an experiment sample and a control sample is used as calibrator. Triplicate duplex reactions employing sequence specific probes labeled with FAM for the target genes and JOE for the housekeeping gene.

Plate Layout <->

Filter + Dyes

Filter 520 nm FAM Filter 550 nm JOE

Filter 580 nm Filter 605 nm

Sample Vol. 25 µl

Probe TaqMan

Background twintec ss HSF

	1	2	3	4	5	6	7	8	9	10	11	12
A	D_Trgt1 2:HKG	D_Trgt1 2:HKG	D_Trgt1 2:HKG		D_Trgt2 2:HKG	D_Trgt2 2:HKG	D_Trgt2 2:HKG		D_Trgt3 2:HKG	D_Trgt3 2:HKG	D_Trgt3 2:HKG	
B	D_Ctrl1 2:HKG	D_Ctrl1 2:HKG	D_Ctrl1 2:HKG		D_Ctrl2 2:HKG	D_Ctrl2 2:HKG	D_Ctrl2 2:HKG		D_Ctrl3 2:HKG	D_Ctrl3 2:HKG	D_Ctrl3 2:HKG	
C	subassay for target gene 1				subassay for target gene 2				subassay for target gene 3			
D												
E												

Sample Info

Position B9

Sample Name D_Ctrl3

Type Unknown
Calibrator

Target 1 Trgt3
Dye FAM

Target 2 Hgk1
Dye JOE
HKG

Comment

4.2.9 Sample types for melting point analysis

No special aspects have to be taken into account for melting curve analysis.

4.2.10 Sample types for endpoint determination

To perform endpoint determination standards, negative controls and unknown samples need to be ascertained. See section 4.2.5 for further details.

4.2.11 Sample types for a +/- assay

To carry out a +/- assay positive and negative controls have to be defined in addition to the unknown samples.

4 Programming (Assay Setup)

4.2.12 Sample types for gene identification

To identify genes and alleles in unknown samples, a positive control has to be defined for every gene or allele. This is shown in the following example.

Allele 1

- Select area B2 – B4 in the plate layout.
- Open dialog **Positive control**.
- Select both alleles in the panel **Target gene(s)**.
- The concentration of the positive control can be specified if known.
- Selection of the genes as Allelic Control using the relevant check box.
- Confirm with **OK**.
- Grouping of the samples at positions B2 – B4.

Allele 2

- Select area C2 – C4 in the plate layout.
- Open dialog **Positive control**.
- Select both genes in the panel **Target gene(s)**.
- The concentration of the positive control can be specified if known.
- Selection of the genes as Allelic Control using the relevant check box.
- Confirm with **OK**.
- Grouping of the samples at positions C2 – C4.

Allele 1

Allele 2

After the samples have been edited, the plate layout looks as follows:

	1	2	3	4	5
A					
B		Allele 1 1: 1000	Allele 1 1: 1000	Allele 1 1: 1000	
C		Allele 2 2: 1000	Allele 2 2: 1000	Allele 2 2: 1000	
D					

4 Programming (Assay Setup)

4.2.13 Copying and cutting of samples

The commands **Copy**, **Cut** and **Paste** required for this procedure can be accessed either via the right mouse button or the menu item **Plate Layout**.

To copy samples of the plate layout, these are first selected and then copied to the clipboard using the **Copy** command. Then a position in the plate layout is selected and the sample type inserted using the **Paste** command.

The **Cut** command can be used to cut samples and then insert them at another location in the plate layout using the **Paste** command.

4.2.14 Create Links

Create Links is used for special settings in a relative quantification assay:

- To define samples as calibrators.
- To link samples (target and housekeeping gene).

For details see section 4.2.8.

4.2.15 Deleting samples

To delete one or more samples they are first selected and then deleted by selecting the icon **Delete Samples**.

If the icon **Clear Plate** is selected, all samples of a plate layout are deleted.

Deleting samples in assays, which are already started or completed is possible (see section 4.2.18). The software automatically saves all fluorescence data concerning a certain position in the plate layout. Only entries like name, amount etc. are deleted.

4.2.16 Grouping of several samples

In addition to the option described in section 4.2.6, it is possible to create replicates by grouping together samples already present in the plate layout. The prerequisite for grouping is that the sample type and all other data correspond. The sample names in a group do not have to be identical.

If several samples are highlighted in the plate layout and the icon **Group** then selected, these samples are grouped together to permit the automatic statistical analysis of replicates.

The grouped samples are marked in the plate layout by highlighting in color on the bottom right. The color can be changed by reselecting the icon **Group as replicates**.

	1	2	3	4
A				
B		NTC	NTC	NTC

Grouping can be undone using the icon **Ungroup replicates**.

4 Programming (Assay Setup)

4.2.17 Determination of subassays

To perform various assays on a PCR plate, the individual assays are marked as subassays to permit separate analysis by the software at a later date. This means it is possible to carry out multiple assays in parallel, whereby each assay may involve different samples with different target genes and dyes. An example is shown in the plate layout, where two different target genes are to be determined in two different standard series and different unknown samples.

The prerequisite for this is that both assays include the same PCR program and that acquisition of the fluorescence data is carried out in the same temperature step.


If all samples of a subassays are highlighted in the plate layout and the icon **Create Subassay** selected, these samples are marked by a frame.

Note: Alternatively blank positions in the plate layout may initially also be selected and defined as a subassay. This may be followed by determination of the individual sample types of the assay.

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B		TNF_1E... 1: 1E9	TNF_1E... 1: 1E9	TNF_1E... 1: 1E9	Mr Jones		HIV_1E11/1 2: 1E11	HIV_1E11... 2: 1E11	HIV_1E11... 2: 1E11	Sample 1		
C		TNF_1E... 1: 1E8	TNF_1E... 1: 1E8	TNF_1E... 1: 1E8	Mr Jones		HIV_1E1... 2: 1E10	HIV_1E1... 2: 1E10	HIV_1E1... 2: 1E10	Sample 1		
D		TNF_1E... 1: 1E7	TNF_1E... 1: 1E7	TNF_1E... 1: 1E7	Mr Jones		HIV_1E9/1 2: 1E9	HIV_1E9/2 2: 1E9	HIV_1E9/3 2: 1E9	Sample 1		
E		TNF_1E... 1: 1E6	TNF_1E... 1: 1E6	TNF_1E... 1: 1E6	Mrs Li		HIV_1E8/1 2: 1E8	HIV_1E8/2 2: 1E8	HIV_1E8/3 2: 1E8	Sample 2		
F		TNF_10... 1: 100000	TNF_10... 1: 100000	TNF_10... 1: 100000	Mrs Li		HIV_1E7/1 2: 1E7	HIV_1E7/2 2: 1E7	HIV_1E7/3 2: 1E7	Sample 2		
G		NTC	NTC	NTC	Mrs Li		NTC	NTC	NTC	Sample 2		
H												


The display of the samples of a subassay during performance or analysis (see section 6.1.1.2) is performed separately and data is saved under separate subitems in the navigation tree.

4 Programming (Assay Setup)

The individual subassays and their position in the plate layout are shown under the node **Monitoring** and **Analysis**. If individual subassays are selected, the related data are loaded to the work field.

The name of a subassay is generated automatically and includes the well positions of the subassay (e.g. Subassay B2-G5). The subassay can be renamed by right mouse clicking the name field in the navigation tree. 32 characters are allowed in length. For specification of the characters which are permitted see section 3.5.3.

Grouping can be undone using the icon **Delete Subassay**. To do so the samples first have to be selected.

Note: The prerequisite for the performance of two different assays in one PCR plate is that both assays are performed with the same PCR program and that acquisition of the fluorescence data is carried out in the same temperature step.

4.2.18 Editing of a plate layout of a completed assay

After the PCR program of an assay has been started or the assay is completed, the plate layout is locked. To make subsequent changes, the plate layout has to be set to **Edit Mode**.

The blocked plate layout of a completed assay can be reactivated by selecting the icon **Edit Mode** for plate layout. Once editing has been performed, the plate layout is deactivated again by reselecting the icon.

All subsequent changes in the plate layout are documented in the UserLog.

EditMode allows changing and deletion of initially defined samples using the tools described in chapters 4.2.5. to 4.2.17. Deleted samples will be marked with a gray background color indicating that fluorescence data for these positions had been or is being collected during the run. Deletion of samples in Edit Mode excludes this data from analysis but the data itself is not deleted from the assay. The data can be reaccessed by redefining samples at these positions.

Note: Note that no data is acquired for positions where no samples were set up before starting the assay. By using Edit Mode already present analysis settings are reset to default.

4 Programming (Assay Setup)

4.2.19 Color for Samples

If the icon **color for samples** is selected, the color of the sample and the corresponding curves can be changed.

Different basic colors are available. User-specific colors can be defined via the color field. For allocating the same color to different samples, select the corresponding positions before the dialogue **color** is opened.

4 Programming (Assay Setup)

4.3 Creation of a PCR program

The following functions are available in the toolbar to edit the PCR program:


- 1 Insert Step
- 2 Edit Step
- 3 Delete Step
- 4 Clear PCR Program
- 5 Show Gradient
- 6 Set Measuring Point

Note: Alternatively these functions can also be accessed via the menu item **PCR Program** or the right mouse button.

4.3.1 Create a PCR program with a program template

If the sub-item **PCR Program** is selected in the navigation tree, a program template, which consists of a 3 Step Cycle, will be displayed in the work field for a new assay.

To edit the program, individual steps are selected and the icon **Edit Step** then used to call up the dialogue in which the data for temperature, time and number of cycles can be amended. Alternatively the dialogue **Edit Step** can be opened by double-clicking on a step. In addition, it is possible to use this dialogue to make all other settings for increments, gradient function and heating and cooling rates (see section 4.3.3.1).

Alternatively the data for temperature, time and number of cycles can be amended directly in the PCR program by selecting and then adding the change.

Note: The properties of the heated lid and the thermoblock can be amended in **Program Header** (see section 4.3.7).

4.3.2 Create a PCR program without a program template

To create a program without a template, the icon **Clear PCR program** can be used to delete the entire program and a new program created by inserting program steps.

4 Programming (Assay Setup)


4.3.3 Insertion and editing of program steps

The cursor is used to select the area of the PCR program in which the new program step is to be inserted. Steps are always inserted to the left of the highlighting.

If the icon **Insert Step** is selected, the relevant dialogue is displayed and can be used to insert various program steps.

If the program step is selected, it will be inserted in the PCR program to the left of the highlighting.

After the program step has been inserted, the data for temperature, time and number of cycles can be directly amended in the PCR program. To do so, they are selected and the change then made.

If other adjustment options are to be amended for program steps, the icon **Edit Step** can be used to call up the dialogue in which all other settings and changes can be made.

4.3.3.1 Temp

A temperature step can be inserted at any point in the program. The main parameters of the temperature step, temperature and holding time, are shown in the work area of the PCR program and can be directly edited there.

TEMP	nominal temperature in degrees Celsius
Permissible values:	4.0 °C to 99.0 °C
Input steps:	0.1 °C
Time:	holding time in minutes and seconds
Permissible values:	00:01 to 99:59 mm:ss
Input steps:	1 s

Additional functions of the temperature step can also be accessed via this window.

Temp incr. – Temperature increment:

Increase / decrease in the temperature of the programmed temperature step to the next cycle by this amount.

Default setting:	0.0 °C
Permissible values:	0.0 to 10.0 °C
Input steps:	0.1 °C
Temperature increase:	+
Temperature decrease:	-


4 Programming (Assay Setup)

The increment is marked by a star in the step of the PCR program: $\frac{55\text{ °C}^*}{00:30}$

Gradient:

When using the gradient function, every column of the thermoblock is set to a different temperature. To program a temperature gradient activate the **Gradient** check box and add the minimum and the maximum temperature into the **From** and **To** fields, respectively. Activating **Gradient** disables the **Temp Incr** and **Time Incr** fields.

Default setting:	deactivated
from temp:	temperature of left block side
Permissible values:	30.0 °C to 99.0 °C
Input steps:	0.1 °C
to temp:	temperature of right block side
Permissible values:	30.0 °C to 99.0 °C
Input steps:	0.1 °C


If the icon **Show Gradient** is selected, a dialogue is displayed, showing the temperatures of the individual columns.

A thermoblock can be selected in the panel Parameters. In addition, the temperatures of the gradient can be entered here.

4 Programming (Assay Setup)

4

The gradient is marked in the step of the PCR program by additional underlining.

55 °C
00:30

The maximum spread of the gradient depends on the thermo module used (see Technical data). If the gradient function is switched off, there is homogeneous temperature distribution over the entire block.

Time incr. – Time increment:

Increase / decrease in the holding time of the programmed temperature step to the next cycle by this amount.

Default setting: 00:00
 Permissible values: 00:00 to 01:00
 Input steps: 1 s
 Time increase: +
 Time decrease: -

The increment is marked by a star in the step of the PCR program.

55 °C
00:15*

Ramp – Heating and cooling rate:

Speed at which the set temperature of the programmed step is achieved.

Default setting: 100 %
 Permissible values: 1 to 75 % and 100 %
 Input steps: 1 %

For procedures involving slow heating or cooling, different values can be entered for the ramp or the bar in the graph moved using the cursor. The current value will then be displayed in the input field. The heating and cooling rates depend upon the thermo module used (see Technical data).

4 Programming (Assay Setup)

4.3.3.2 Cycle (1 – 5 Step Cycle)

The Cycle command contains the number of temperature steps (1 – 5) in one cycle and the number of repetitions used for execution of the temperature commands.

Default setting: 40
 Permissible values: 1 to 99
 Input steps: 1

4.3.3.3 Hold

The program holds the temperature at the entered value until the PCR program is ended. To proceed press Start.

Default setting: 4.0 °C
 Permissible values: 4.0 °C to 99.0 °C
 Input steps: 0.1 °C

4.3.3.4 Pause

The program is stopped, and the temperature remains at the value of the last temperature command. The program continues after the key **Start program** is pressed.


4.3.3.5 Sound

It is possible to emit an acoustic signal, e.g. to announce a programmed pause or the end of the program.

Default setting: 9
 Permissible values: 1 to 9 (number of signals)
 Input steps: 1


4.3.3.6 Melting Curve

Selection of the item **Melting Curve** in the dialogue **Insert Step** automatically inserts four program steps. Besides the actual step of the melting curve, this includes an initial denaturation step as well as a temperature step each for the initial and final temperature.


To amend the parameters for determination of the melting point, the step of the melting curve is selected, followed by the icon **Edit step**.

4 Programming (Assay Setup)

- Ramp time: Input of the time in which the melting curve is to be created.
Permissible values: 5 – 30 min
Input steps: 1 s
Fluorescence data are measured continuously.
- Resolution: Alternatively the measuring rate for melting curve can be entered.
Permissible values: 0.1 – 2.0 °C
Input steps: 0.1 °C
Fluorescence data are measured after heating of x °C is finished.

In the PCR program the melting curve is marked by a symbol.

4.3.3.7 Endpoint

After the dialogue **Insert Step** has been opened choose the entry **Endpoint**. The data for temperature and duration can be easily edited by clicking at the values in the graph.

- TEMP nominal temperature in degrees Celsius
Permissible values: 4.0 to 99.0 °C
Input steps: 0.1 °C
- Time: holding time in minutes and seconds
Permissible values: 01:00 to 99:59 mm:ss
Input steps: 1 s

4.3.4 Copying of program steps

The commands **Copy** and **Paste** required for this procedure can be accessed either via the right mouse button or the menu item **PCR program**. To copy a program step, this is first selected and then copied to the clipboard using the **Copy** command. Then a position in the PCR program is selected and the program step inserted there using the **Paste** command.

Here the program step is always inserted on the left next to the highlighting.

4.3.5 Deletion of program steps

Individual temperature steps are deleted by first being selected and then deleted by selecting the icon **Delete Step**.

4 Programming (Assay Setup)

4.3.6 Determination of measuring point

In general, data acquisition can be carried out during PCR program, endpoint assay or melting point analysis. Each item can be programmed once in an assay.

The measuring point is determined in a temperature step by highlighting a step and selecting **Set measuring point** in the **Temp Step** dialogue.

Attention:

The measuring point is defined by the probe that is used. It is not set up by the software automatically!

Note:

Data acquisition can only be carried out in one temperature step within the PCR cycle. The PCR program can be composed of repeated PCR cycles. Data acquisition can be carried out in each PCR cycle, but the software based evaluation is only possible in the first PCR cycle.

Attention:

When setting measuring points it should be ensured that the duration of the temperature steps is sufficient. With one dye per filter rotor 8 seconds are required for data acquisition while the measuring time is extended to 16 seconds for two dyes per filter rotor.

The measuring point is marked in the step of the PCR program by a relevant symbol.

55 °C 
00:30

To remove the measuring point the relevant temperature step is first selected and then the check box **Set measuring point** deactivated in the dialogue **Edit step** or the icon **Set measuring point** selected in the toolbar.

4.3.7 Header settings

The header settings form the program head in which the controller strategy and the behavior of the heated lid are each set on a program-related basis.

Program Header			
Lid Temp	<input type="text" value="105"/> °C	<input checked="" type="checkbox"/> TSP Heated Lid	<input type="checkbox"/> Switch off Lid at low Block Temp
Temp. Mode	<input type="text" value="Fast"/> ▼	<input type="checkbox"/> Impulse	<input type="checkbox"/> Simulate Mastercycler gradient

Lid Temp

The heated lid temperature of 105 °C is preset and cannot be changed.

Temp.Mode

This command is used to specify the type of regulation for the block temperature mode.

Standard Default setting / block temperature mode for standard applications. This mode should be used for sample volumes in the range from 20 µl to 50 µl. Additionally, this setting will be recommended if the mode **fast** results in poor amplification results.

Fast The temperature at the thermoblock is measured and the nominal temperature then set. This mode should be used for Fast-PCR and low sample volumes (< 20 µl) as well as for templates with low GC-content. The program run time is decreased in comparison to the mode **standard**.

TSP Heated Lid

The PCR program is only started after the lid temperature has been achieved. If this function is not selected, the PCR begins directly after the program start. For endpoint assays deactivation of the combo box is not possible.

4 Programming (Assay Setup)

Impulse

For devices with silver block (thermo module Mastercycler ep gradient S) it is possible to accelerate the first heating process using the Impulse function. If another block is selected or if the function **Simulate Mastercycler gradient** is selected, this function is not active.

Switch off lid at low block temp

Depending on the input the heated lid is switched off at low block temperatures. If block temperatures are held below 15 °C for more than 5 minutes, the heated lid is switched off automatically.

If this function is not selected, the lid heating is only switched off at the end of a program.

Note: If the last command in a program is a Hold step, the heated lid remains switched on until the key is operated. The prerequisite for this is that the temperature of the Hold step is <15 °C and the function **Switch off lid at low block temp** is deactivated.

Simulate Mastercycler gradient

PCR programs adopted from other Cyclers may have to be re-optimized as the temperature characteristics vary and may influence the PCR results.

Simulate Mastercycler gradient permits simulation of the temperature characteristics of the Mastercycler gradient to guarantee the simple transfer of established programs. This setting allows the temperature control behavior of a Mastercycler gradient to be simulated and thus to establish protocols from the Mastercycler gradient in general without further optimization on the Mastercycler ep *realplex*.

To do so, the PCR program of the Mastercyclers gradient is entered in Assay Setup under the PCR program and the check box **Simulate Mastercycler gradient** then activated in the program header.

After being saved the program can be started on the Mastercycler ep *realplex* without re-optimization of the parameters.

5 Operation (Monitoring)

5.1 Sample loading

The Mastercycler ep *realplex* can be equipped with up to 96 PCR tubes or one 96-well PCR plate. To ensure uniform pressure of the heated lid when using tubes it is recommended positioning a tube in each corner of the thermoblock. Other tubes should be distributed evenly in the block.

Please remember that with real time PCR analyses both the tubes and the sealing systems selected must be suitable for fluorescence detection (see Ordering information).

For optimum temperature transition it should be ensured that the tubes are properly positioned in the thermoblock and do not wobble. Tubes that are not sufficiently heat resistant (approx. to 120 °C) should not be used. This generally also applies to covers used to seal PCR plates.

Note: Different tube types should not be used together as they will then not be sealed properly by the heated lid. In addition, no suitable calibration data have been stored. When using individual 0.2 ml PCR tubes, tubes should also be positioned in each corner of the thermoblock.

Attention: PCR plates should only be labeled at the side as labeling may also be fluorescent!

5.2 Block and Lid temperature control

For incubation experiments or to accelerate the start of a method, both the thermo module and the heated lid of the *realplex* module can be subjected to temperature control. To do so, the required block or lid temperature is selected under the cycler node in the dialogue **realplex System configuration** (see section 3.7).

Note: At the start of a program the settings of the respective program have priority. During an ongoing program separate block temperature control of the thermoblock cannot be switched on.

5.3 Start assay


After the created assay has been saved and the sample tubes inserted, the heated lid is slid forwards. The sealing clamp is then moved down and the program started.

Note: Make sure that the instrument has been warmed up for at least 15 minutes before starting an assay.


If the icon **Start current cycler program** is selected, the assay is started. Alternatively this function is also available under the menu item **Command**.

After the assay has started, the dialogue **Start PCR** is displayed and can be used to select the correct background plate with a combo box.

5 Operation (Monitoring)

If a background plate has already been specified in the plate layout, this is automatically shown in the combo box. It is possible to select a background plate other than the one already specified in the plate layout. If this is done, it is documented in the report (see section 6.1.10).

The data acquisition can be observed in the mode **Monitoring**.

The fluorescence intensity is shown here as a function of the number of cycles and time. Underneath the temperature profile is shown in the course of time.

If the check box **Legend** is selected, a legend is inserted on the right of the fluorescence profile.

With the acquisition of measured data for multiplex assays it is also possible to choose between individual dyes in the panel **Show Dyes** by selecting the relevant check box.

5.3.1 Quick start of an assay

To quick start an assay, select the dyes and the background to be used and enter the sample volume. Then select all positions in the plate layout and define them as **Unknowns**. All fields of the dialogue may be left blank. Enter the PCR program and start the run. While the assay is running or at any later time before analysis, the actual plate layout can be entered and unused samples deleted using **Edit Mode** (see 4.2.18).

5 Operation (Monitoring)

5.4 Selection of individual samples

All samples are listed in the navigation tree under **Samples**.

This table shows the names of the individual samples and their positions in the plate layout. In addition, the samples are marked with the symbols of the sample types. To facilitate the assignment of individual fluorescence curves to the samples, the color of the fluorescence curves is also shown in the table.

Position	Name
- A9	NTC
! B9	1
! C9	10
! D9	100
! E9	1000
! F9	10000
! G9	100000
- A10	NTC
! B10	1
! C10	10
! D10	100
! E10	1000
! F10	10000
! G10	100000

If one or more samples are selected, the related fluorescence curves are shown in the work area. To select several samples in one row, first the upper sample is highlighted and then the lowest sample by simultaneously pressing the <Shift> key. If several samples which are not positioned next to each other are to be selected, the <Ctrl> key has to be pressed during selection.

5.5 Selection of individual sample types

To show or hide individual sample types they are activated or deactivated in the window **Show Sample Type**.

Show Sample Type	
<input checked="" type="checkbox"/>	? Unknowns
<input checked="" type="checkbox"/>	! Standards
<input checked="" type="checkbox"/>	+ Pos. Controls
<input checked="" type="checkbox"/>	- Neg. Controls

5.6 Cyclor Status

Directly after the start all information for the ongoing PCR program is shown under **Monitoring** in the panel **Cyclor Status**. The progress of an ongoing program can be monitored in the status display.

Cyclor Status		
Step 3	Step 55.0°C	Start 12:42:06
Cycle 5/40	Block 55.0°C	End 13:48:11
Time 00:00:13	Lid 105°C	Remaining Run Time 00:57:22


The current cycle is shown in the left-hand block of the panel, in addition to the current temperature step and the time still remaining for this temperature step.

The middle block of the panel **Cyclor Status** shows the nominal temperature of the current program step and the current temperatures of the block and heated lid.

The right-hand block provides information about the start and end time of the ongoing PCR program and the time period for which this program will last.

5 Operation (Monitoring)

By changing from **Monitoring** to **PCR Program** it can also be seen in the PCR program which is the current temperature step. This is shown in red.

5.7 Interruption of an ongoing assay

The assay can be interrupted, e.g. to add a reagent. The current temperature of the block and the lid remain unchanged.

The icon **Pause current cycler program** can be used to interrupt the ongoing assay.

Cycler Status: Pause The status display then shows **Pause**.

The program can be continued at any time (see section 5.8) or aborted (see section 5.9).

5.8 Continue assay

If the assay is to be continued, the heated lid has to be closed first.

A program that has been stopped can be continued at any time with the icon **Start current cycler program**.

Cycler Status: Running After the assay has been continued, the status display then shows **Running**.

5.9 Abort assay

An assay is interrupted when the icon **Pause current cycler program** is selected.

Afterwards the assay is aborted by selecting the icon **Abort current cycler program**.

Cycler Status: Idle The information **Idle** is shown in the status display.

5 Operation (Monitoring)

5.10 Use of Mastercycler ep *realplex* for standard PCR

The Mastercycler ep *realplex* can also be used solely as a thermo cycler. To do so, a new assay is opened (see section 3.10.2) and only a PCR program created in the Assay-Setup (see section 3.5.1). Make sure that no measuring point is set in the case of a standard PCR. After creation of the PCR program it is saved and started.

Note: The plate layout remains blank, and no dye should be selected as otherwise there will be an error message on start.

5.11 Use of Mastercycler ep *realplex* as Fluorimeter

The Mastercycler ep *realplex* can also be used for fluorescence measurement independently of a PCR. In this case solely endpoint determination is carried out (see section 4.2.10).

Here it must be remembered that background calibration is also necessary for such determination. In addition, color calibration must be performed for every dye to be detected.

Note: On calibration the temperature of the block and heated lid must correspond to the conditions of subsequent fluorescence measurement.

5.12 Shut down

To shut down the Mastercycler ep *realplex* the ongoing assay is first ended and saved. The software is then closed, computer shut down and the device switched off.

Saving of assays is only necessary if the data were evaluated, the raw data will be saved automatically.

Note: The *realplex* module should be kept in a closed position when the Mastercycler ep *realplex* is in off state to prevent contamination of the block.

6 Analysis

An assay is analyzed by first being opened with the command **Open Assay** and then selecting the node **Analysis** in the navigation tree. If the assay finished shortly before and no other assay was opened afterwards, the data will be automatically accessible under **Analysis**.


Note: While one assay is running, it is possible to analyze an assay that has already been completed under **Analysis** (see section 6.1.1.1).

Information about the plate layout and PCR program can be viewed under the relevant node of the navigation tree when the relevant assay also appears under **Assay Setup**. If another assay is running at the same time, access to these data is not possible.

6.1 General

6.1.1 Open Assay


Open an assay using **Open Assay**.

After the assay is opened, the raw data are displayed automatically under **Analysis** in the right-hand work field. These data are the fluorescence values acquired during detection which were set against the factors from the background and color calibration. Display of the raw data differs depending on the analysis module and is described in the individual sections for evaluation of the different analysis modules.

6.1.1.1 Open assay during an ongoing assay

During data acquisition of an ongoing assay it is possible to analyze an assay that has already been completed. To do so it is selected in the dialog **Assays**. Another dialog appears, stating that the assay is for evaluation only as another assay is currently running.


After confirmation with **OK** the assay is displayed under **Analysis**.

6 Analysis

6


Analysis

Note: During data acquisition of an ongoing assay it is not possible to view information on the plate layout or PCR program of the assay already completed. Sample types and their positions can be seen in the Mini Plate Layout under Analysis.

6.1.1.2 Open subassay

If subassays have been defined in the plate layout, they are shown under **Analysis** as individual nodes. The designation **Subassay B2-G5** indicates the position in the plate layout (for renaming a subassay see section 4.2.18). If a subassay is selected, only the data associated with the subassay are loaded to the work field where they can be analyzed.

All analyses of the sub-assays are saved under the assay.

6.1.2 Selection of individual samples

All samples are listed in the Navigator under **Samples**.

Position	Name
A9	NTC
B9	1
C9	10
D9	100
E9	1000
F9	10000
G9	100000
A10	NTC
B10	1
C10	10
D10	100
E10	1000
F10	10000
G10	100000


This table shows the names of the individual samples and their positions in the plate layout. In addition, the samples are marked with the symbols for the sample types. The color is also shown in the table to simplify the assignment of individual fluorescence curves to the samples. In the individual analysis modules this table is extended into a table of results.

If one or more samples are selected, the associated fluorescence curves are shown in the work area. To select several samples in one row, first highlight the top sample and then the bottom sample simultaneously pressing the <Shift> key. If several samples which are not positioned next to each other are to be selected, the <Ctrl> key has to be pressed during selection.

The sequence of the samples can be changed by double-clicking on the heading of individual columns in the sample table.


6 Analysis

Alternatively, samples can be selected using the Mini Plate Layout:

If the Mini Plate Layout is not visible it can be activated using the **Show / Hide Mini Plate Layout** button.

In the default view mode the check box **Show All** in the right upper corner of the results table is activated. Whereas in the work area only the selected samples are presented as curves or columns, all samples which have been defined in the plate layout are listed in the sample table.


6 Analysis

6

Analysis


By deactivating the check box **Show All** another mode of sample presentation is chosen. In this case the results table shows only the samples selected in the **MiniPlate Layout**. If an assay contains a lot of samples, the sample table will be more clearly arranged by this way.

Note: A Report shows the currently selected view mode.

6.1.3 Deactivation of samples

To exclude samples from analysis they can be deactivated. To do so one or more samples are highlighted in the sample table under **Analysis** and the function **Deactivate Sample(s)** selected using the right mouse button (alternatively: menu item "Commends").

The deactivated samples are marked in the sample table by a cross through the sample type symbol.

6 Analysis

Position	Name
- A2	NTC
! B2	1
 C2	10
! D2	100
! E2	1000
! F2	10000
! G2	100000
- A3	NTC
! B3	1
! C3	10
! D3	100
! E3	1000
! F3	10000
! G3	100000


In addition, the deactivated sample is marked with a cross in the plate layout.

	1	2	3
A		NTC	NTC
B		1 3:1	1 3:1
C		10 3:10	10 3:10
D		100 3:100	100 3:100

The analysis is then updated automatically. In order to reactivate the deactivated samples, the sample in the sample table is highlighted again and the function **Activate Sample(s)** selected using the right mouse button.

6.1.4 Selection of individual sample types

To show or hide the curves of individual sample types in the work field, they are activated or deactivated in the window **Mini Plate Layout** using the appropriate check box.

This selection does not apply to the sample table, nor does it affect the analysis – unlike the function **Deactivate Sample(s)** or **Activate Sample(s)**.


6 Analysis

6

Analysis

6.1.5 Selection of analysis module

In the top panel different analysis modules can be selected using the combo box **Type of Application**. These data are then displayed in the work field according to the module. In addition, a legend can be inserted by selecting **Legend** on the right.

Note: If the required analysis module is not available in the selection, the settings necessary for this analysis have not been made in the plate layout.

Subsequent alteration of the plate layout is possible in the mode **Edit mode for plate layout** (see section 4.2.19).

If analysis modules are not available, this may also be due to the fact that settings necessary for this analysis are not included in the PCR program. In this case subsequent alteration is not possible.

With multiplex assays, the combo box **Dye** can also be used to select individual dyes. It is possible to display the data selectively for one dye or all of the dyes at the same time via **all dyes**.

Special settings for analysis (e.g. setting of threshold) are not possible in the view **all dyes**.

6.1.6 Zooming charts

The following icons are available to magnify and shrink the data in the work field:

The zoom tools are deactivated using **Unselect Zoom Functions**.

Parts of the charts can be magnified using **Zoom In**.

Charts can be shrunk by selecting **Zoom Out**.

If the icon **Unzoom** is selected, altered charts will be displayed according to the default.

To magnify or shrink charts, first select the relevant icon and then click on the chart with the tool then available.

This process is repeated until the chart attains the required size. To deactivate the tool select the icon **Unselect Zoom Functions**.


6 Analysis

6

Analysis

6.1.7 Modifying the scaling of charts

By default, all charts generated in the **Analysis** section are auto scaled. To set the axis scalings manually, right click inside the chart area and select **Scaling** from the context menu.

In the **Chart Scaling** dialog all charts from the modules that have been included in the analysis are listed in the tree. Select the radio button **Manual** for the chart to be rescaled and enter the minimum and maximum values for the X- and Y-axis. By clicking the **View** button, the modified chart can be previewed.

Note: The View button is inactive if the selected chart belongs to an analysis module other than the currently active.

6 Analysis

6

Analysis


6.1.8 Restore default settings

Restore Defaults

If this function is selected, the defaults will be restored for all parameters altered.

6.1.9 Copy and save data

For further analysis in other software programs, copy and save of fluorescence data is possible. Fluorescence data (raw data, PCR data, and melting curve data) are exported as indicated in the software (taking into account e.g. curve smoothing and base line correction).

To copy and save fluorescence data of samples the corresponding entry in the popup menu of the chart has to be selected. The software exports the data which are presented in the chart, i.e. only the samples selected in the sample table. Data of deactivated samples are not included. The data are displayed line-by-line (A1 -> A12).

The command **Copy Data** can be used to first copy the fluorescence data to the clipboard and then insert it in other documents.

If the command **Save Data As...** is selected, the file selection dialogue will be shown. The following data formats can be chosen: Text files (*.txt), comma separated (*.csv), and Microsoft Excel 97 (*.xls).

Note: Copying data of melting curves with a variety of data points per well could result in difficulties if pasted to Microsoft Excel or other spreadsheet applications due to a limitation in maximum column numbers.

6.1.10 Copy and save charts

To save and copy individual charts first select with the cursor and then press the right mouse button. A popup menu will appear in the chart (see figure section 6.1.9).

If the command **Save Chart As Metafile Picture** is selected, the chart can be saved in EMF format.


The chart can be opened in Windows Explorer by double-clicking on it. In the program <Windows Picture and Fax Viewer> the chart can then be saved in other image formats or transferred to the image editing program <Windows Paint>.

The command **Copy Chart To Clipboard** can be used to first copy the chart to the clipboard and then insert it in other documents.


6 Analysis

6.1.11 Copy and save results as TXT, CSV, and XLS data format

Relevant commands in the menu item **Edit** can be used to copy and save results from the table **Samples/Analysis** in the navigation tree. Alternatively, the table Sample/Analysis can be selected with the cursor and a popup menu then opened with the right-hand mouse button.

The selection dialogue enables user-specific settings like data format as well as sample and analysis selection.

Note: Deactivated samples will not be exported by the software.

The following settings can be selected:

- Data format: The appropriate data format can be chosen by selecting the corresponding entry in the combo box.
- Sample selection: The selection can be carried out by activating the radio buttons for **Selected Samples** or **All Samples**. By choosing **Selected Samples**, only the samples selected in the results table are exported. **All Samples** refers to all samples present in the plate layout. The check box **Across (Sub-) Assays** will only be of interest, if Subassays are in use.

6 Analysis

6

Analysis

- Selection of Analysis data (only selectable for Text format): Depending on what is selected all given columns of the results table or only a limited number of columns will be exported. For deselection of columns, the corresponding check box has to be deactivated. **Include Analysis Parameters** gives information concerning the analysis settings, e.g. threshold and baseline settings.

After the command **Copy Results** is selected the required settings have to be made in the selection dialogue. After confirming the dialogue by **OK**, the corresponding data is saved to the clipboard and can be inserted in another document. **Save Results As...** can be used to save the data from the table Sample/ Analysis as TXT, CSV, and XLS data format. For this purpose choose the format TXT first and close the dialog by **OK**. File directory, file name, and file type have to be chosen in the file selection dialogue **Save data as...**


6.1.12 Copy and save results to special software tools for relative quantification

There are several specific software tools for calculation of gene expression and related relative quantification assays. The *realplex* application enables export of data to the most commonly used of these software tools (qBase, geNorm, REST and GenEx).

Note: For further information concerning the software tools please see corresponding literature and internet websites. In this manual, only the procedure of exporting data to these tools is described.

qBase	http://medgen.ugent.be/qBase/
geNorm	http://medgen.ugent.be/genorm/
REST	Relative Expression Software Tool
	http://www.gene-quantification.com
GenEx	http://www.multid.se

Access to **Copy Results** and **Save Results As...** is given as described in section 6.1.11. In the selection dialogue the software tool of interest can be selected in the combo box.

The entries qBase, geNorm, REST 200x and GenEx will only be present, if the relevant assay contains PCR data. Melting curve and endpoint data can not be used. Because results and sample information have to be arranged in a format specific to each of these software tools, the analysis data selection is deactivated.

6 Analysis

The following options are available:

- Software tool
- **All Samples or Selected Samples**
- In the case of All Samples the check box **Across (Sub-) Assays** has to be activated if necessary.

File names of exported data are created based on the assay name and the corresponding ending of the file type. The software checks the exported data sets automatically with regard to completeness and consistency in sample and target name. If necessary, replicate sample numbers are removed automatically. If the requirements for the corresponding software tool are not fulfilled, one of the following error messages will be displayed.

- "The plate layout contains replicate groups with inconsistent sample names and targets. To resolve this, make sure that all samples within a replicate group have the same sample names and targets."
- "The plate layout contains samples with undefined sample names and/or target. To resolve this, enter sample names and target names for all samples."
- "The plate layout contains samples with duplicate combinations of sample name and target, which are not grouped as replicates. To resolve this, either group these samples as replicates or rename the samples."

In this case inconsistencies have to be corrected first. Afterwards **Copy / Save Results** can be repeated.

6.1.12.1 qBase

Export of real time PCR results to qBase is offered for qBase Version 1.3.5 and higher. Only samples of type Unknown, Standard, and Negative Control are exported in this format. In the All Dyes view, the format qBase can not be selected. If results of multiplexed assays are to be transferred to qBase, export the data of each dye to an extra file and import them in qBase to a single experiment data set. Output of data occurs line-by-line, sorted by well position.

Note: When importing assay data in qBase select "qBase" in the import dialog, since the "Eppendorf" import is not yet implemented in the current version of qBase.

Available file formats for qBase are: XLS (recommended), TXT, and CSV.

6.1.12.2 geNorm

Export of real time PCR results to geNorm is offered for geNorm Version 3.4 and higher. Only samples of type Unknown are exported in this format. Output of data based on C_T values is computed as described in the geNorm documentation and is arranged in a matrix of sample names and target names.

Available file formats for geNorm are: XLS (recommended), TXT, and CSV.

6.1.12.3 REST 200x

Export of real time PCR results to REST is offered for REST 200x Version 1.9.9 and higher. Only samples of type Unknown and Standard are exported in this format. For REST, the following information is exported from the Plate Layout: target names, housekeeping genes, calibrators. REST automatically calculates amplification efficiency for a target gene if a corresponding standard curve is present.

For REST 200x the only file format available is REST.

6.1.12.4 GenEx

Export of real time PCR results to GenEx is offered for GenEx Version 4.2. and higher. Only samples of type Unknown are exported in this format. Amounts of Mean Amounts (for replicate groups) are exported if a corresponding standard curve is present for each target to be exported. Otherwise C_T values or Mean C_T values are exported. Data are arranged as a matrix of sample names and target names.

Available file formats for GenEx are: XLS (recommended), TXT, and CSV.

6 Analysis

6


Analysis

6.1.13 Compile reports

Reports can be used to make extensive compilations of experimental parameters and results which go beyond the standard screen display. They are output as a preview on screen, a printout to a printer or as a printable file in PDF format. After the function **Report** is selected under the menu item **File**, the dialog **Report Configuration** appears. This function is only active when in the module **Analysis**.

Alternatively, the dialog **Report Configuration** can also be opened using the icon **Print Report**.

By checking individual check boxes in the panel **Data** the user can select the data to be listed in the report. This includes general document information and the data from the assay setup. This is followed by tab sheets for raw data and the analyses, which have been saved for the assay or currently created. By selecting the appropriate check box the charts of the fluorescence data and the analyses can be inserted in the report as well as the associated legends. In addition, the user can select the dyes for which data are to appear in the report. To exclude complete analyses and all raw data from the report the top check box is deactivated in the relevant tab sheet. In the default setting all data of an assay are included in the report.

Note: It is only possible to select the analyses already performed by the user with the data record concerned under **Analysis** or previously saved by him. The function keys **Preview**, **Print** and **Export** are only accessible if at least one printer driver is installed on the computer!

With the function **Page Header & Footer**, items such as company names and logos can be inserted in the header line of the report.

Edit

Alternatively, the dialog **Edit Page Header & Footer** can also be opened with the function key **Edit**.

The name of the company or institute can be entered in the field **Company Name**.


Clear Logos already present in the field **Company Icon for Report** can be deleted.

Load... If this function is selected, logos are loaded to the field **Company Icon for Report**.
Note: Only icons in BMP format can be inserted (size: max. 250 x 50 Pixel).

6.1.13.1 Work with report templates

User-specific report templates can be saved in the panel **Report Template** to simplify the preparation of reports.

Save After amending the report this function key is used to open the dialog **Save template**. The name of the report is then entered and the input confirmed with **OK**.

To use a report template to prepare a report, a template is selected after opening the dialog **Report Configuration** using the combo box **Select Template** in the panel **Report Templates**.

Delete Report Templates can be deleted with this function key.


6 Analysis

6

Analysis

6.1.14 Save analyses

If the acquired data of an assay are analyzed, the user is automatically prompted to save the settings on opening another assay or closing the assay.


This also applies when changes have been made to an existing analysis.

6.1.15 Call up analyses

When opening an assay for which analyses have already been saved, the raw data first appear in the work field in the mode **Analysis**. If the relevant analysis is selected in the combo box **Type of Application**, the analysis appears with all saved settings.

6.2 Quantification

After completing an assay or opening a completed assay the raw data are loaded to the work area under **Analysis**. Here the fluorescence intensity is displayed as a function of the number of cycles and time. Underneath the temperature profile is shown in the course of time.


6 Analysis

6

Analysis

If **Quantification** is selected in the combo box **Type of Application**, the fluorescence curves are automatically shown in logarithmic display. If the data are to be presented in logarithmic display, the check box **log scale** must be activated.


The header now also shows the field **Data Source**, indicating which data source was used as the basis for the analysis displayed. Basically there are three different data sources: endpoint determination, PCR data and melting curve analysis. The data sources which are available for analysis are determined by the PCR program in Assay Setup.

Note: Different data are only available for selection under **Data Source** with the analysis modules Gene Identification and +/- Assay.

6.2.1 Determination of threshold

The threshold is used to specify C_T values of samples. Here the C_T value is the cycle in which the fluorescence signal intersects with the threshold. The threshold is determined using the field on the right-hand side, with the data points of the exponential phase of the amplification curves being used to calculate the C_T values.

Three different methods for specifying the threshold are available in the panel **Threshold**:


The radio button **Noiseband** can be used for automatic determination of the threshold by the software. Here the threshold is specified so that it is significantly (10 times the standard deviation) above the noise of the baseline. **Noiseband** is preselected as the default setting by the software.

6 Analysis

6

Analysis

The function **CalQplex** can be used for automatically analysing single curves (Eppendorf algorithm). Fluorescence data are analysed based on a continuous growth model in combination with a baseline correction. The parameters of a perfect PCR curve are aligned to the experimental fluorescence curve of a given sample. The **CalQplex** method results in a reference cycle number for each fluorescence curve. This method does not rely on a threshold level. The reference cycle numbers are referred to as **Reference Points**. They correspond to the C_T values of the conventional methods and will be therefore displayed in the column C_T values. In default view mode the Reference Points are not shown in the amplification plots. To change this, the item **Show Reference Points CalQplex** has to be chosen in the right mouse popup menu. Reference Points are marked as crosses.


With the function **Best R²** the software locates the threshold in the area in which all curves of the standards have a maximum value for correlation. If no standards have been defined, this function is inactive. At least three standard concentrations are necessary for this function to be active.

If the radio button **Manually** is selected, the location of the threshold is specified manually by inputting a fluorescence value. Alternatively, the threshold can also be directly moved in the chart using the left mouse button.

6.2.2 Determination of baseline

The panel **Baseline Settings** is used to specify the parameters for calculation of the baseline.

Automatic calculation of the baseline is selected as the default setting. In this case the optimum baseline is calculated for every sample individually.

To adjust the baseline the radio button **Manual Baseline** is selected, now making the previously inactive **Cycle Range** accessible. In the relevant fields the user can enter the minimum and maximum cycle between which the baseline is to be calculated for all curves. Here the area between cycle 3 and 15 is given as the default setting.


In addition, the drift in the baseline can be corrected in this panel by activating the check box **Drift Correction**. This correction will affect all curves equally.

Note: All settings made in this panel affect the curve shapes and the C_T values.

6 Analysis


6.2.3 Standard curve

The software automatically uses the C_T values of the standards to calculate a standard curve which shows the C_T values as a function of the amount of nucleic acid used. The unit of the standard concentrations was determined when editing the standards in the plate layout.


Slope	Increase in standard curve
Y-Intercept	Point at which the standard curve intersects with the Y-axis (based on amount = 1 [unit])
Efficiency	Efficiency of the PCR $E = 10^{[-1/\text{slope}] - 1}$ The value for E should be approx. 1
R ²	Correlation coefficient of standard curve

The measured values of the standards are marked with red squares (see above) and of the unknown samples with black crosses. If individual standards are selected in the sample table, the associated measured value is indicated by a yellow square in the standard curve.

If individual unknown samples are selected, the calculated C_T value and associated concentration of the unknown sample are indicated by red lines in the standard curve.

6 Analysis

6

Analysis

6.2.4 Analysis

The standard curve allows the user to conclude the initial amount of nucleic acid of the unknown samples from a specific C_T value. The results are calculated automatically by the software and can be found in the table **Sample/Analysis**.

Sample/Analysis								
Pos	Name	Ct FAM ...	Ct Mean ...	Ct Dev. ...	Amount ...	Amount Mean ...	Amount Dev. ...	Target ...
! A1	Std 1	11.69	11.62	0.07	1.00E+08			Lambda
! A2	Std 1	11.59	11.62	0.07	1.00E+08			Lambda
! A3	Std 1	11.57	11.62	0.07	1.00E+08			Lambda
! A4	Std 1	11.64	11.62	0.07	1.00E+08			Lambda
! A5	Std 1	11.60	11.62	0.07	1.00E+08			Lambda
! A6	Std 1	11.58	11.62	0.07	1.00E+08			Lambda
! A7	Std 1	11.74	11.62	0.07	1.00E+08			Lambda
! A8	Std 1	11.56	11.62	0.07	1.00E+08			Lambda
? A9	Sample A	11.72	11.79	0.08	8.21E+07	7.88E+07	4.15E+06	Lambda
? A...	Sample A	11.90	11.79	0.08	7.34E+07	7.88E+07	4.15E+06	Lambda
? A...	Sample A	11.81	11.79	0.08	7.78E+07	7.88E+07	4.15E+06	Lambda
? A...	Sample A	11.72	11.79	0.08	8.19E+07	7.88E+07	4.15E+06	Lambda
! B1	Std 2	14.99	15.07	0.12	1.00E+07			Lambda
! B2	Std 2	15.17	15.07	0.12	1.00E+07			Lambda
! B3	Std 2	15.11	15.07	0.12	1.00E+07			Lambda
! B4	Std 2	15.29	15.07	0.12	1.00E+07			Lambda
! B5	Std 2	15.13	15.07	0.12	1.00E+07			Lambda

Pos

Specifies the position of the sample in the plate layout. For simple assignment the samples are also marked with the symbols for the respective sample type.

Name

Name of sample

Ct

C_T value

If an internal reference dye is used in a concentration that gives relative fluorescence signals < 100, the entry **Ref.dye too low** will appear. In this case, no curve is shown in the amplification curve.

IPC

Internal positive control (only present if selected in plate layout)

Amount

Initial amount of nucleic acid, specified by the user with these standards. With unknown samples and negative controls this is calculated by the software based on the standard curves.

Mean

Mean value for C_T value and amount, respectively, when determining unknown samples in replicates. **Note:** Only included when replicates defined in plate layout.

Deviation

Standard deviation for C_T value and amount, respectively, when determining unknown samples in replicates. **Note:** Only included when replicates defined in plate layout.

Target Gene

Target gene (name determined in plate layout when editing sample).

Note: The sequence of the samples can be changed by double-clicking on the heading of individual columns in the sample table.

6.2.5 Multiplex analysis

If different genes are quantified in parallel in a multiplex analysis, a separate standard curve must be created for every gene. To do so the relevant dye is selected in the top panel using the combo box **Dye**.

If **All Dyes** is selected in the combo box **Dye**, the curve shapes for all dyes are displayed. In this mode it is not possible to display the standard curves at the same time.

6 Analysis

6.3 Relative quantification

The Relative Quantification analysis module is used to compare expression levels of a target gene among differently treated samples. The expression levels are calculated using the $\Delta\Delta C_T$ method (Livak and Schmittgen. 2001, *Methods* 25(4): 402–408). A prerequisite to yield valid relative expression values is that the efficiency is comparable *and* close to 1 in the PCR systems of both target gene and housekeeping gene.

The expression levels R are calculated as follows:

$$R = 2^{-\Delta\Delta C_T} = 2^{-(\Delta C_T \text{ sample} - \Delta C_T \text{ calibrator})}$$

$$= 2^{-[(C_T \text{ sample} - C_T \text{ housekeeping gene}) - (C_T \text{ calibrator} - C_T \text{ housekeeping gene})]}$$

The assignment of target and housekeeping gene pairs and calibrators is made through the **Create Links** dialog. The relative quantification analysis option is only available if a suitable set of samples has been defined in the plate layout (see chapter 4.2.9) and the PCR run is completed.

6.3.1 Performing a relative quantification analysis

1. In the navigator window move to the **Analysis** section and select **Relative Quantification** from the **Type of Application** combo box.
2. Select a **Dye** from the combo box and mark the samples to be included in the **Sample / Analysis** list. The amplification curves are then displayed in the graph window. Adjust the **Threshold** and **Baseline Settings** according to chapters 6.2.1 and 6.2.2.
3. For multiplex assays, repeat step 2 for all additional dyes used.

The resulting expression levels are displayed in the **Sample / Analysis** list.

6.3.2 Results table for monoplex assays

For every sample individual C_T values and mean C_T and standard deviation of replicates are displayed. The mean expression levels and the range of individual expression levels are listed in columns **Expr. Level** and **Range** for positions where a non-calibrator target gene is defined. The expression level of the calibrator target gene is 1 per definition and is not displayed. Expression level values correspond to a fold expression of the target gene compared to the calibrator sample.

Sample/Analysis							
Pos	Name	Ct FAM...	Ct Mean FAM...	Ct Dev. FAM...	Expr. Level FAM...	Range FAM...	
? B5	Exp_B	24.78	24.65	0.14	2.1	1.8 - 2.5	
? B6	Exp_B	24.74	24.65	0.14	2.1	1.8 - 2.5	
? B7	Exp_B	24.47	24.65	0.14	2.1	1.8 - 2.5	
? B8	Exp_B	24.60	24.65	0.14	2.1	1.8 - 2.5	
? B9	Exp_C	22.56	22.47	0.11	10	8.8 - 11	
? B...	Exp_C	22.33	22.47	0.11	10	8.8 - 11	
? B...	Exp_C	22.55	22.47	0.11	10	8.8 - 11	
? B...	Exp_C	22.45	22.47	0.11	10	8.8 - 11	
? C1	Control	25.56	25.80	0.19			
? C2	Control	26.02	25.80	0.19			
? C3	Control	25.79	25.80	0.19			
? C4	Control	25.82	25.80	0.19			
? C5	Exp_A	26.01	25.85	0.17	0.93	0.80 - 1.1	
? C6	Exp_A	25.76	25.85	0.17	0.93	0.80 - 1.1	
? C7	Exp_A	25.97	25.85	0.17	0.93	0.80 - 1.1	

Pos

Specifies the position of the sample in the plate layout. For simple assignment the samples are also marked with the symbols for the respective sample type. The sample defined as the calibrator in the plate layout is marked with a C.

Name

Name of sample.

6 Analysis

Ct FAM

C_T value for target gene (here: FAM-marked probe).
Mean value and standard deviation for C_T value are shown, when determining unknown samples in replicates.

Expr. Level...

Expression level of sample calculated according to $\Delta\Delta C_T$ method.

Range

Range of Expression level according to $\Delta\Delta C_T$ method.

6.3.3 Results table for multiplex assays

If **All Dyes** is selected from the **Dyes** combo box, individual C_T values and mean C_T and standard deviation of replicates for all dyes used are displayed. The mean expression levels and the range of individual expression levels are listed in columns **Expr. Level** and **Range** for all dyes used for detection of a non-calibrator target genes is defined. The expression level of the calibrator target gene is 1 per definition and is not displayed.


Sample/Analysis									
Pos	Name	Ct FAM...	Ct Mean FAM...	Ct Dev. FAM...	Expr. ...	Range F...	Ct JOE...	Ct Mean JOE...	Ct ...
? A1	Control	25.92	25.10	0.74			28.84	27.87	0.80
? A2	Control	24.33	25.10	0.74			27.36	27.87	0.80
? A3	Control	25.50	25.10	0.74			28.18	27.87	0.80
? A4	Control	24.65	25.10	0.74			27.09	27.87	0.80
? A5	Exp_A	25.00	25.53	0.40	0.92	0.53 - 1.6	27.39	28.17	0.87
? A6	Exp_A	25.44	25.53	0.40	0.92	0.53 - 1.6	28.71	28.17	0.87
? A7	Exp_A	25.83	25.53	0.40	0.92	0.53 - 1.6	27.46	28.17	0.87
? A8	Exp_A	25.85	25.53	0.40	0.92	0.53 - 1.6	29.10	28.17	0.87
? A9	Exp_B	24.57	24.78	0.17	1.8	1.1 - 2.7	28.36	28.37	0.52
? A...	Exp_B	24.75	24.78	0.17	1.8	1.1 - 2.7	29.00	28.37	0.52
? A...	Exp_B	24.96	24.78	0.17	1.8	1.1 - 2.7	27.72	28.37	0.52
? A...	Exp_B	24.86	24.78	0.17	1.8	1.1 - 2.7	28.40	28.37	0.52
? B1	Exp_D	22.92	22.85	0.14	6.5	4.2 - 10	28.48	28.32	0.77

6 Analysis


6.4 Melting curve analysis

Melting curve analysis is generally performed after a PCR by continuously measuring the fluorescence intensity of the reaction samples during a slow increase or fall in temperature. This process can be used to determine the melting temperature (T_m) of a PCR product (e.g. marked with SYBR-Green I) or when using specific probes, the melting temperature of the resulting hybrids.

After completing an assay or opening a completed assay the raw data are loaded to the work area under **Analysis**. Here the fluorescence intensity is displayed as a function of time. Underneath the temperature profile is shown in the course of time.

To analyze the melting points the analysis module **Melting Curve** first has to be selected in the combo box **Type of Application**. The fall in the fluorescence intensity is shown in the chart below as a function of temperature. At the melting point the measured signal curves show a maximum in the fall in fluorescence intensity.


6 Analysis

6

Analysis

To allow the maximum intensity changes to be identified more easily the melting curves are shown as peak curves. These peak curves are the negative first derivative of the melting curves as a function of temperature ($-dl/dT$). The unit of the y-axis is given in %, with 100 % representing the highest peak of the assay.


The location of the peak maximum corresponds to the melting temperature of the product. This is determined by the software and marked by vertical lines. A vertical line will only appear if a single sample is selected.

6.4.1 Determination of threshold

A threshold is determined for the melting curve analysis. Here it should be borne in mind that it is the peaks above the threshold that are analyzed, not those below. This is shown as a horizontal line in red and can be moved using the spin edit to the right of the chart. The default setting is 33 %.

Alternatively, the threshold can also be changed by manually moving the red line.

Melting points of individual samples can be checked manually using the green vertical line. To do so it is moved to the point of the peak maximum, and the melting point can then be read off on the right under the threshold input.

6 Analysis

6.4.2 Analysis

The melting points of the samples determined by the software can be found in the table **Sample/Analysis**.

Sample/Analysis		Show All		
Pos	Name	No. Tm...	Tm x (°C) SYBR	Tm y (°C) SYBR
- A1	NTC	1	84.0	
- A2	NTC	1	83.8	
- A3	NTC	1	83.7	
? A4	500	1	83.8	
? A5	500	1	83.9	
? A6	500	1	83.9	
? A7	1000	1	83.9	
? A8	1000	1	83.9	
? A9	1000	1	83.9	
? A...	2000	1	83.8	
? A...	2000	1	83.9	
? A...	2000	1	84.0	

Pos

Specifies the position of the sample in the plate layout. For simple assignment the samples are also marked with the symbols for the respective sample type.

Name

Name of sample

No. Tm

Indicates whether and how many melting points were identified.

Tm(s) (°C)

Melting point of sample determined by the software.

Tm x (°C)

Melting point no.1 of sample determined by the software.
In the case of two melting points, the lower temperature is listed here.

Tm y (°C)

Melting point no.2 of sample determined by the software.
In the case of two melting points, the higher temperature is listed here.
If only one melting point is present, there will be no entry in this column.

Mean

Mean value for Tm(s) when determining unknown samples in replicates.
Note: Only included when replicates were defined in the plate layout and only one melting point is determined per sample and dye.

Deviation

Standard deviation when determining unknown samples in replicates.
Note: Only included when replicates were defined in the plate layout and only one melting point is determined per sample and dye.


The software is capable of determining max. two melting points per sample. If a curve has more than two peaks above the threshold, no entry will be made in the sample table. In this case analysis can be carried out manually.

6 Analysis

6.5 Endpoint measurement


Endpoint measurement is performed by including the step **Endpoint** in the PCR program. Here data acquisition can take place after the creation of PCR data, but also independently of data acquisition during a PCR (see section 5.12).

After completing an assay or opening a completed assay the raw data are loaded to the work area under **Analysis**. Here the fluorescence intensity is displayed as a function of the time. Underneath the temperature profile is shown in the course of time.

If **Endpoint** is selected in the combo box **Type of Application**, the fluorescence data of the endpoint measurement are shown as bars in all allocated positions of the PCR plate or block when using individual tubes. The scale range of the height of the bars is related to the maximum value ascertained on endpoint measurement. This takes up the entire height. Negative values are not shown here but can be found in the sample table. Replicates defined in the plate layout as a group are indicated by colored triangles (see section 4.2.16).

6 Analysis


If samples have been defined as negative controls in the plate layout, the fluorescence values of the negative controls can be deducted from all other measured values by activating the check box **Neg. Control Subtraction** in the panel **General**. With more than one negative control their mean value is subtracted. In the presence of more than one group of negative controls, the mean of all of them is used.

All measured values can be found in the table **Samples/Analysis**.

6.5.1 Standard curve

If standards have been defined in the plate layout, a standard curve is calculated by the software. This shows the fluorescence intensity as a function of the amount used. The unit is determined when editing the standards in the plate layout.

- Slope Increase in standard curve
- Y-Intercept Point at which the standard curve intersects with the Y-axis
- R² Correlation coefficient of standard curve

The measured values of the standards are marked with red squares and of the unknown samples with black crosses.

6 Analysis

6

Analysis

6.5.2 Analysis

The standard curves are used to calculate the amount for every measured fluorescence value of an unknown sample or negative control. They can be found in the table **Samples/Analysis**.

Sample/Analysis						
Pos	Name	Endpoint SYBR	Amount SYBR [ng/μl]	Mean SYBR	Dev. SY...	
H2	Control	115	5.85E-02	-1.02E-02	9.72E-02	
H3	Control	74	-7.89E-02	-1.02E-02	9.72E-02	
F2	Sample B	206	0.364	0.431	9.48E-02	
F3	Sample B	246	0.498	0.431	9.48E-02	
B2	Sample C	2322	7.46	7.53	0.107	
B3	Sample C	2367	7.61	7.53	0.107	
A2	Std 1	3085	10.0			
A3	Std 1	3181	10.0			
C2	Std 2	346	1.00			
C3	Std 2	388	1.00			
D2	Std 3	260	0.500			
---	---	---	---	---	---	

Pos

Specifies the position of the sample in the plate layout. For simple assignment the samples are also marked with the symbols for the respective sample type.

Name

Name of sample

Endpoint

Fluorescence intensity of sample measured during endpoint determination.

Amount

Initial amount of nucleic acid, specified by the user with these standards. With unknown samples and negative controls this is calculated by the software based on the standard curves.

Mean

Mean value when determining unknown samples in replicates.

Note: Only active when replicates were defined in plate layout.

Deviation

Standard deviation when determining unknown samples in replicates.

Note: Only active when replicates were defined in plate layout.


6 Analysis

6.6 +/- Assay

Data from endpoint measurement or a PCR can be analyzed in the analysis module **+/- Assay**. When the assay is opened, raw data for a PCR run and/or for endpoint determination will appear accordingly in the right-hand work field.



6.6.1 +/- Assay with endpoint determination as data source

If **+/- Assay** is selected using the combo box **Type of Application** followed by the data source **Endpoint Data** with the combo box **Data Source**, the fluorescence data of this endpoint measurement are shown as bars in all allocated positions of a PCR plate. Here grouped replicates are indicated by colored triangles (see section 4.2.16).

6.6.1.1 Determination of threshold

The determination of a threshold defines when a sample should be seen as positive (+) or as negative (-). Here the samples below the threshold are categorized as negative and those above it as positive. To do so the type of threshold first has to be specified in the panel **General**. The threshold is shown as a red line in all occupied positions of the PCR plate.


6 Analysis

6

Analysis

If the radio button **Threshold** is selected, a fixed threshold can be specified in the panel **Threshold Settings**. This can be carried out either by entering a numerical value or by moving the track bar. The greatest possible threshold corresponds to the maximum value determined on endpoint measurement.


The location of the threshold is marked by a red line.

Note: Every change in the threshold has a direct influence on whether a sample is categorized as negative or positive.

If the radio button **Threshold Range** is selected, it is possible to enter a fluorescence range for the threshold in the panel **Threshold Settings**. Here too either a numerical value can be directly entered or the left track bar moved. If samples show a fluorescence within the range, they will be automatically categorized as inconclusive.

The threshold range is marked in grey.

Note: If the fluorescence range of the threshold is changed, this may lead to a change in the categorization as negative, positive and inconclusive results.

If samples have been defined as negative controls in the plate layout, the fluorescence intensity of the negative controls can be deducted from all other measured values by activating the check box **Neg. Control Subtraction**. If more than one negative control has been edited in the plate layout, the mean value is subtracted.

6.6.2 +/- Assay with PCR data as data source

If **+/- Assay** is selected using the combo box Type of Application followed by **PCR Data** under **Data Source**, the fluorescence curves will appear in the right-hand work field.

6.6.2.1 Determination of threshold

The determination of a threshold defines when a sample should be seen as positive (+) or as negative (-). Here the samples below the threshold are categorized as negative and those above it as positive. There are no inconclusive results. See section 6.2.1 for further details for the determination of thresholds.


6 Analysis

6.6.3 Analysis


View

The function **View** in the top panel can be used to display the results in the work area.

The results are shown as different-colored bars of the same height in the allocated positions of the PCR plate.

The legend is shown in the panel **+/- result**. Besides positive and negative results inconclusive results are also displayed. They are indicated by lying within the **Threshold range** and being marked as **inconclusive**. If no **Threshold range** but a **Threshold** was specified, a result marked **inconclusive** is not possible.

In addition, the results are shown in the table **Samples/Analysis**. The descriptions **positive**, **negative** and **inconclusive** are shown in the column **+/- result [Dye]** beside the fluorescence data.

Sample/Analysis			
Pos	Name	+/- Result	Endpoint
?	F7 sample 8	positive	16913
?	E7 sample 7	positive	11881
?	D7 sample 6	inconclusi...	8214
?	C7 sample 5	negative	5401
?	F6 sample 4	positive	12453
?	E6 sample 3	positive	12638
?	D6 sample 2	inconclusi...	8614
?	C6 sample 1	negative	5306
-	B6 Neg. Control	negative	1214
-	B7 Neg. Control	negative	1300
+	F4 Pos. Control	positive	16058
+	F5 Pos. Control	positive	17138

Evaluation

If this function is selected, the measured values of endpoint measurement or of the PCR are loaded to the work area again.

6 Analysis

6

Analysis


6.6.4 Multiplex analysis

If different genes are determined in parallel in a multiplex analysis, separate analysis must be carried out for every target gene. To do so the relevant dye is selected in the top panel using the combo box **Dye**.


All results can be displayed if **All Dyes** is selected in the combo box **Dye**.

With endpoint determination in the mode **Evaluation** the fluorescence values of the different dyes are shown as blue bars for every sample. The assignment of the bars to the target genes analyzed can be found in the legend.

Note: If PCR Data is selected in the combo box **Data Source**, the relevant curve shapes will appear in the mode **Evaluation**.

If you switch to the mode **View** the results of analysis will be shown as different-colored bars for every sample. The assignment of the bars to the target genes analyzed can be found in the legend. All results are listed in the table **Samples/Analysis**.

6 Analysis

6.7 Gene Identification

Gene identification can be based either on a PCR or endpoint determination as the data source. After the analysis module **Gene Identification** has been selected in the header as **Type of Application**, the data from the source forming the basis for gene identification will appear in the work area. If more than one data source is available, they can be selected using the combo box **Data Source**. For defining allelic controls under plate layout see 4.2.13.

6.7.1 Determination of threshold

Before actual analysis of the gene identification can be performed, the relevant dye must be selected for every allele in the top panel using the combo box **Dye**. Subsequent determination of the threshold then depends on whether a PCR or endpoint determination is used as the data source. The process can be found in the relevant sections for quantification (6.2.1) or +/- Assay (6.6.6.1).


6.7.2 Analysis

Once the threshold has been determined separately for every allele, all dyes are selected in the top panel using the combo box **Dye**.


View

The function **View** in the top panel can be used to display the results in the work area.

The results are shown as different-colored bars in the allocated positions of the PCR plate. If one allele is detected, it is shown as a single bar. If both alleles are detected in a sample, two different colored bars are seen.

The panel **Allele settings** shows the legend indicating the color-coded assignment of the alleles automatically set by the software. The dye used to detect the relevant allele is also indicated.

The results are listed in the table **Sample/Analysis**.

6 Analysis

6

Analysis

Sample/Analysis				
Pos	Name	Ct FAM	Genelident Result ...	Ct VIC
+ C5	Allel 1	18.55	Allel 1	
+ D5	Allel 2		Allel 2	19.88
? C3	Sample 001	32.10	Allel 1, Allel 2	33.12
? C4	Sample 001	31.82	Allel 1, Allel 2	33.25
? D3	Sample 002	31.86	Allel 1, Allel 2	33.45
? D4	Sample 002	32.14	Allel 1, Allel 2	32.59

Pos

Specifies the position of the sample in the plate layout. For simple assignment the samples are also marked with the symbols for the respective sample type.

Name

Name of sample

Ct FAM

C_T value for allele 1 (here: FAM-marked probe).

Genelident Result...

The alleles detected for the sample are given under **Genelident Result**.

Ct VIC

C_T value for allele 2 (here: VIC-marked probe).

Results for the unknown samples are only shown in the column **Genelident Result** if a C_T value has been determined for the positive controls for the alleles or the values are above the threshold, in the case of endpoint determination.

Note: Positive controls for alleles are always given as positive in the column **Genelident Result**, including where no C_T value was determined for them or the values are below the threshold, in the case of endpoint determination.



Evaluation

If this function is selected, the measured values of endpoint measurement or of the PCR are loaded to the work area again.

6 Analysis

6.8 Inverted data

There are real-time PCR chemistries which result in decreasing fluorescence signals the more PCR product is generated. Curves with decreasing fluorescence levels are shown in the Raw Data View. For analyzing these data in the different analysis modules of the realplex software, activate the check box **Inverted Data**. Subsequently, the curves are inverted and threshold and baseline adjustment can be done. Analysis can be performed as described in chapter 6.2 – 6.7.

Note: Analysis of endpoint data will not be possible for assays showing decreased fluorescence levels. Activation of the check box **Inverted Data** comes along with an entry in the document information of the report (see chapter 6.1.13).

7 Maintenance

7

Maintenance

The device may only be opened by a qualified technician. The warranty for the device will become null and void if there is unauthorized intervention in the device!

7.1 General instructions on cleaning of device

When maintenance and cleaning are carried out, the device must be switched off and the plug removed beforehand.

The thermoblock, the inner side of the heated lid and micro test tubes / PCR plates can reach temperatures of over 50 °C very quickly. There is a risk of burns!

The heated lid should be kept closed until temperatures of approx. 30 °C or less are reached.


Figure 7: Cleaning position

- 1 Head lock arm (locked)
- 2 Sealing clamp
- 3 *realplex* module (tilted)
- 4 Rollflex cable
- 5 Thermoblock

For cleaning the device, the heated lid can be opened. To do so the *realplex* module must be moved to the rear position. The heated lid can then be put in a vertical position and then locked in position using the head lock arm (see Figure 7).

While the *realplex* module is lifted up, the sealing clamp (see Figure 7) of the lid should not be lowered. This could damage the photomultiplier!

When moving down the heated lid, it should be held with one hand and the head lock arm pulled forwards with the other hand. The heated lid is then lowered carefully.

Do not insert your fingers between the lid and the housing of the device when opening or closing the lid or when opening to the cleaning position as they may get caught.

It should be borne in mind that when cleaning the device, the sharp edges of the heated lid may cause injury e.g. rollflex cable and guide rail.

The device should be protected from organic solvents and aggressive chemicals. Make sure that no liquids enter the device.

The operator is responsible for carrying out appropriate decontamination if hazardous material is split on or into the device.

The surface of the Mastercycler ep *realplex* may be wiped down with a damp cloth or a mild lab cleaner.

Note: Cleaning of the heated lid must not be performed. For further information, please call the Eppendorf Service.

7 Maintenance

Before cleaning or decontamination methods other than those specified by the manufacturer are used, the user should check with the manufacturer that the intended method will not damage the device.

Note: Cleaning and decontamination may also be necessary as precautionary measures if the device or accessories are to be sent in for maintenance, repair or passed on. The manufacturer supplies a form which the operator can use to confirm that the device has been subjected to such treatment.

Note: When switched off the *realplex* module should always be kept closed to prevent contamination of the block.

7.2 Cleaning of thermoblock

Before cleaning of thermoblock turn off the instrument and let instrument cool for 10 – 15 minutes.

- Clean contaminated positions using a cotton swab soaked with 70 % ethanol (stock solution: Ethanol p.a. (>99.8 %) undenatured). Turn the cotton swab carefully in the well.
- Dry the position with a new cotton swab.
- Repeating this process 3 to 5 times with new clean cotton swabs.
- Switch on the device and perform another background calibration (see chapter 3.8.1). Wells should be completely dry before using instrument.

If there are still positions marked as failed, repeat the procedure with new cotton swabs. In the case that the contamination can not be removed by this procedure, please call the Eppendorf Service.

7.3 Validation and calibration of thermoblock

Validation and calibration of the thermoblock is performed according to the instruction manual of the Temperature-Validation System for the Mastercycler ep.

8 Troubleshooting

Errors in the programming or handling of the device can be corrected by eliminating the cause shown in the display.

Technical faults may be caused by malfunctions, e.g. power failure or power fluctuations. In most cases it is sufficient to briefly switch the device off and back on after approx. 10 seconds. If the fault re-occurs, please note down the error code (e.g. 0x0B04) shown in the display and contact Eppendorf Service.

8.1 General Errors

Error	Cause	Solution
Background calibration was not successful	Contamination in thermoblock	Decontamination Chapter 7.2
	Contamination of background plate	Use a new background plate
	Contamination of <i>realplex</i> modul	
Color calibration was not successful	Evaporation in wells	Use a new color calibration plate
	Wrong background calibration was used during color calibration	Repeat with the right background calibration
	Contamination of color calibration plate	Use a new color calibration plate
	Wrong dye concentrations used or dye bleached	Use a new color calibration plate
No amplification can be detected	Component of reaction mix is missing	Check the reaction mix
	Reagents were handled wrong	<ul style="list-style-type: none"> • Store reagents at $-20\text{ }^{\circ}\text{C}$ • Protect the reagents from light • Avoid repeated thaw and freeze cycles
	Air bubbles	Centrifuge at 300 x g
	Reaction mix is not at the bottom of wells	Centrifuge at 300 x g
	Data acquisition has occurred in wrong PCR step or at wrong temperature	<ul style="list-style-type: none"> • Check PCR Program • Check measuring point
	Data acquisition in wrong positions	Check plate layout
	Wrong pipetting	<ul style="list-style-type: none"> • Check sample volume in wells • Use epMotion
	Non-optimized reagents	Titration of reagents (MgCl_2)
	Inhibition by contamination from nucleic acid purification	<ul style="list-style-type: none"> • Use kits for ultra pure nucleic acid purification • Do not use more than 10 μl DNA in a total volume of 20 μl • Check PCR system with positive control
	Large amplicons may reduce the PCR efficiency	Use short amplicons for real-time PCR
	GC-rich primer or probe	<ul style="list-style-type: none"> • Add DMSO to reaction mix • Redesign primer and probe
	Primer or probe does not bind	Check primer and probe design
	Logarithmic phase is not reached until end of PCR	Low template concentration
PCR system not optimized		<ul style="list-style-type: none"> • Titration of reagents (e.g. MgCl_2) • Check primer and probe design
Not enough cycles		Perform PCR with more cycles

8 Troubleshooting

Error	Cause	Solution
Amplification curve too flat	Low PCR efficiency	<ul style="list-style-type: none"> • Titration of reagents (e.g. MgCl₂) • Check primer and probe design • Check purity of template
Amplification curve decrease at the end of PCR	Due to the large amount of products at the end of PCR the re-annealing of these products prevents annealing of probes	The so-called "hook effect" has no influence on the analysis results
Negative control is positive	Primer-Dimer in SYBR Green assays	<ul style="list-style-type: none"> • Check primer design • PCR setup at 4 °C • Use Hotstart enzymes • Use higher amounts of template • Increase annealing temperature • Optimize MgCl₂ concentration
	Contamination	<ul style="list-style-type: none"> • Strict separation of nucleic acid purification, PCR setup and analysis • Use to avoid contamination from already finished assays
Poor reproducibility and accuracy	Wrong pipetting	<ul style="list-style-type: none"> • Use master mixes • Use epMotion or electronic pipettes • Pipette large volumes • Reduce amount of pipette steps
	Wrong pipetting	<ul style="list-style-type: none"> • Check sample volume in wells • Use epMotion
	Insufficient mixing	Vortex master mixes and samples for 5 s
	Air bubbles	Centrifuge at 300 x g
	Drops on sealing option	Centrifuge at 300 x g
	Drops on well	Centrifuge at 300 x g
	Contamination with flourogenic material	<ul style="list-style-type: none"> • Use powder-free gloves • Do not write on wells
	Contamination of <i>realplex</i> module	<ul style="list-style-type: none"> • Check with black plate • Call Eppendorf Service
	Low template concentration	unspecific binding of DNA to consumables can be avoided by the use of glycogen
Fluorescence signal too low	Probe or SYBR concentration too low	Check PCR setup
	Reagents were handled wrong	<ul style="list-style-type: none"> • Store reagents at – 20 °C • Protect the reagents from light • Avoid repeated thaw and freeze cycles
	PCR efficiency too low	<ul style="list-style-type: none"> • Titration of MgCl₂ • Check primer and probe design • Final concentration of primer and probe should be 300 – 900 nM • Optimize annealing and denaturing temperature with gradient function • Check PCR program
Fluorescence signal too high	Probe or SYBR concentration too high	Check PCR setup
	Signal enhancement is wrong	Call Eppendorf Service

8 Troubleshooting

8.2 Errors during SYBR Green assays

Error	Cause	Solution
Signals in stationary phase differ	High amount of bound dye due to a high initial template concentration	not use more than 50 ng genomic DNA or 100 ng RNA in a total volume of 20 μ l
	Bleaching effect	<ul style="list-style-type: none"> • Protect SYBR from light • Avoid repeated thaw and freeze cycles
Negative control is positive	Primer-Dimer in SYBR Green assays	<ul style="list-style-type: none"> • Check primer design • PCR setup at 4 $^{\circ}$C • Use Hotstart enzymes • Use higher amounts of template • Increase annealing temperature • Optimize $MgCl_2$ concentration • Perform melting curve analysis to detect primer dimers
No amplification	DNA is single stranded due to the high temperature during data acquisition	Check PCR program
Standard curve not OK	PCR setup does not match plate layout	Check plate layout and PCR setup


8.3 Errors during qRT-PCR

Error	Cause	Solution
Low fluorescence signals	Contamination with RNase	<ul style="list-style-type: none"> • Use powder-free gloves and RNase free H_2O • Store RNA at $-70^{\circ}C$ • Check integrity of RNA
	RT-PCR is inhibited	Use ultra pure RNA
	$MgCl_2$ concentration is not optimal	Titration of $MgCl_2$, concentration should be 3 – 8 mM
	RNA concentration not optimal	Titration of RNA, amount should be 1 – 500 pg
	cDNA concentration not optimal (not valid for two-step PCR)	Titration of cDNA
Amplification curve is arcuate during logarithmic phase	Temperature too high during reverse transcription	Reduce temperature for reverse transcription
	Primer not optimal	<ul style="list-style-type: none"> • Use pure primer • Check design of primers
Additional peak during melting curve analysis	Contamination with genomic DNA	<ul style="list-style-type: none"> • Use mRNA • Check purity of RNA
	Temperature too low during reverse transcription	Temperature during reverse transcription should be higher than $55^{\circ}C$
Peak in negative control during melting curve analysis	DNA contamination	Check all reagents


8 Troubleshooting

8.4 RealSpector

In the case of a potential instrument malfunction it will be helpful to give as much information as possible to Eppendorf Service. For this purpose a special software tool was created called **RealSpector**. The RealSpector application is available in the start menu (Programs / Eppendorf / RealSpector). The following start screen appears:

The necessary measurement readings are recorded automatically after pressing **OK**. The panels **realplex No.** and **Problem Description** may be filled with relevant information, if available. It does not matter which plate is inserted in the thermoblock. All relevant data are stored automatically by the software and will be available in the directory C:\Program Files\Eppendorf\realplex\info as zipped data file (name: "serial number of *realplex* module".zip).

Note: The *realplex* software application has to be shut down prior to a RealSpector run. Cancellation of the measurement process is not made possible.


In very rare cases hardware and/or software communication between *realplex* and thermo module as well as the computer is not possible at all. Then, serial number and problem description have to be entered manually by the user. An attachment of a document to the zipped data file can also be made.

Please send the zipped data file to your Eppendorf Service Organization.

9 Short instructions

9.1 Login


After opening the *realplex* software the Login is displayed.

Select the required user and enter the relevant password.

9.2 Setting up new users

New users can only be set up by the Administrator when he is logged in. To do so, the dialogue **realplex System configuration** is opened under the menu item **Setup**.

The function key **New** is used to set up a new user. Already existing users can be changed or deleted using the functions **Edit** and **Delete**.

9 Short instructions


9.3 Creating assays

If the icon **New Assay** is selected, a blank plate layout and a standard PCR program are loaded in the work area.

After selection of the icon **Open Assay**, the dialogue **Assays** is opened.

The assay marked **Setup** has not yet started. On opening it is loaded to **Assay Setup** as well as to **Analysis**. However, no data are available in the latter.

Assays which can be used as a template for the creation of new assays are marked **Template**. See section 3.6.4 and 3.6.5 for further information about the creation and usage of templates.

In the case of assays marked **Complete** data acquisition has already taken place. If this assay is opened, the data are displayed under **Analysis**. In addition, the assay appears in **Assay Setup** unless another assay is currently underway (see section 6.1.1).

If individual assays are selected, all analyses of this assay saved to date are shown in the lower window **Analysis**. The date of creation and the name of the user who performed analysis are also listed here.

Assays which have been aborted e.g. due to a power outage or system failure are marked **Aborted**. They may also be analyzed if sufficient data were acquired before the assay was aborted.

9.4 Saving assays

A logged on user can save assays under his node or in one of his folders. The editing functions for these programs are only accessible to the logged-in user. Other users are only able to view and copy these assays.


If the icon **Save Assay** is selected, the current assay is saved.

9 Short instructions

9.5 Performing an assay

9.5.1 General


The software of the Mastercycler ep *realplex* basically consists of two modules **Assay Setup** and **Analysis**. The module **Assay Setup** is divided into the sub-items **Plate Layout**, **PCR Program** and **Monitoring** which are needed for creating and running an assay.

They are displayed in the work area on the right if individual sub-items are selected. Alternatively, the different modules can also be loaded to the work area using the functions of the menu item **Commands**.

9.5.2 Plate Layout

The following settings are necessary depending on the assay type:

- Selection of the dyes in **Filters + Dyes**.
- Probe and Background can be selected in the appropriate combo boxes in the right panel. Moreover, the sample volume can be entered.
- Enter the sample types. To do so mark the position in the plate layout and use the following icons:

If the icon **Standard** is selected, the relevant dialogue is opened. This can be extended by the option **Auto Series** for simple creation of standard series and replicates.

The dialogue **New Positive Control** can be opened with the relevant icon.

If the icon **Negative Control** is selected, the relevant dialogue is opened.

If the icon **Unknown** is selected, the relevant dialogue is opened. A simple creation of replicates is possible by extending the dialogue by choosing **Auto Series**.

Note: If there are no dyes, probes or backgrounds available in the relevant combo boxes they have to be entered or created (see Chapter 3.6, 3.7 and 3.8).

9.5.3 PCR Program

After selecting the item **PCR Program** in the navigator a 3StepCycle is displayed in the work area. This PCR program can be modified with the following icons:

If the icon **Insert Step** is selected, the relevant dialogue is displayed and can be used to insert various program steps.

If other adjustment options are to be amended for program steps, the icon **Edit Step** can be used to call up the dialogue in which all other settings and changes can be made.

Individual temperature steps are deleted by first being highlighted and then deleted by selecting the icon **Delete Step**.

The measuring point is determined in a temperature step by highlighting a step and selecting **Set measuring point**.

9 Short instructions

9.5.4 Start Assay

After the created assay has been saved and the sample tubes inserted, the heated lid is slid forwards. The sealing clamp is then moved down and the program started.

If the icon **Start Current Cyclers Program** is selected, the assay is started. Alternatively this function is also available under the menu item **Command**.

After the assay has started, the dialogue **Start PCR** is displayed and can be used to select the background plate with a combo box.

9.5.5 Interruption of an ongoing assay

The assay can be interrupted, e.g. to add a reagent. The current temperature of the block and the lid remain unchanged.

The icon **Pause Current Cyclers Program** can be used to interrupt the ongoing assay.

Cycler Status: Pause The status display then shows **Pause**.

The program can be continued at any time (see section 9.5.6) or aborted (see section 9.5.7).

9.5.6 Continue assay

If the assay is to be continued, the heated lid has to be closed first.

A program that has been stopped can be continued at any time with the icon **Start Current Cyclers Program**.

Cycler Status: Running After the assay has been continued, the status display then shows **Running**.

9.5.7 Abort assay

An assay is interrupted when the icon **Pause Current Cyclers Program** is selected.

Afterwards the assay is aborted by selecting the icon **Abort Current Cyclers Program**.

Cycler Status: Idle The information **Idle** is shown in the status display.

9.5.8 Analysis

In the top panel different analysis modules can be selected using the combo box **Type of Application**. These data are then displayed in the work field according to the module.

The screenshot shows a software interface with two dropdown menus. The first dropdown menu is labeled 'Type of Application' and has 'Raw Data' selected. The second dropdown menu is labeled 'Dye' and has 'FAM' selected. To the right of the 'Dye' dropdown menu is a checkbox labeled 'Legend' which is currently unchecked.

With multiplex assays, the combo box **Dye** can also be used to select individual dyes. It is possible to display the data selectively for one dye or all of the dyes at the same time via **all dyes**.

After selection of an application type and the threshold setting the results are displayed in the table **Samples / Analysis**. Further information concerning the threshold settings can be found in the chapter for each application type (see chapter 6).

10 Technical data

10

Technical data

Mastercycler ep *realplex*

Excitation source	96 LEDs (470 nm)
Emission filters	520 ± 10 nm / 550 ± 10 nm (<i>realplex</i> ²) 520 ± 10 nm / 550 ± 10 nm / 580 ± 5 nm / 605 ± 15 nm (<i>realplex</i> ⁴)
Photomultiplier	1 Channel photomultiplier (<i>realplex</i> ²) 2 Channel photomultiplier (<i>realplex</i> ⁴)
Sample capacity	96 x 0.2 ml PCR tubes or one 96 PCR plate (unskirted, semi-skirted, skirted – as per SBS standard)
Temperature control range of block	4 °C to 99 °C
Range of gradient	1 °C to 20 °C (thermo module Mastercycler ep) 1 °C to 24 °C (thermo module Mastercycler ep S)
Temperature control range of gradient	30 °C to 99 °C
Temperature control range of lid	37 °C to 105 °C
Temperature of lid	105 °C
Block homogeneity	35 °C ± 0.3 °C 90 °C ± 0.4 °C
Control accuracy	± 0.2 °C
Heating rate*	approx. 4 °C/s (thermo module Mastercycler ep) approx. 6 °C/s (thermo module Mastercycler ep S)
Cooling rate*	approx. 3 °C/s (thermo module Mastercycler ep) approx. 4.5 °C/s (thermo module Mastercycler ep S)
Dimensions (W x D x H)	26 cm x 41 cm x 39.6 cm
Total weight	24 kg
Weight of thermo module	17 kg
Weight of <i>realplex</i> module	7 kg
Mains power connection	100 V – 130 V, 50 – 60 Hz, 200 V – 240 V, 50 – 60 Hz (Mains power deviation not greater than ±10 % of nominal mains power)
Power consumption	800 W
Over voltage category	II (IEC 610 10-1)
Contamination level	2
Protection class	I
IP protection	20
Ambient conditions (only indoors)	15 °C to 35 °C, 70 % rel. humidity, to 2000 m sl
Storage conditions	-20 °C to 50 °C, 85 % rel. humidity
Transport conditions	-20 °C to 50 °C, 85 % rel. humidity

* measured at block

Technical specifications subject to change!

Mains power connection supply and power consumption of *realplex* module via connected thermo module. The device has been awarded the CE mark.

11 Ordering information

Please only use the original accessories recommended by ourselves. The use of spare parts or disposables other than those recommended by ourselves may impair the function and safety of the device. Any warranty and liability for damage thus caused will be excluded.

Description	Order no.
Mastercycler ep <i>realplex</i> ² , 100 V – 130 V	950020202
Mastercycler ep <i>realplex</i> ²ⁱ , 100 V – 130 V	6300 000.906
Mastercycler ep <i>realplex</i> ² , 200 V – 240 V	6300 000.507
Mastercycler ep <i>realplex</i> ²ⁱ , 200 V – 240 V	6300 000.957
Mastercycler ep <i>realplex</i> ² S, 100 V – 130 V	950020211
Mastercycler ep <i>realplex</i> ² S, 200 V – 240 V	6300 000.604
Mastercycler ep <i>realplex</i> ⁴ , 100 V – 130 V	950020300
Mastercycler ep <i>realplex</i> ⁴ⁱ , 100 V – 130 V	6302 000.903
Mastercycler ep <i>realplex</i> ⁴ , 200 V – 240 V	6302 000.504
Mastercycler ep <i>realplex</i> ⁴ⁱ , 200 V – 240 V	6302 000.954
Mastercycler ep <i>realplex</i> ⁴ S, 100 V – 130 V	950020318
Mastercycler ep <i>realplex</i> ⁴ S, 200 V – 240 V	6302 000.601
<i>realplex</i> ² , 100 V – 130 V	950020105
<i>realplex</i> ² , 200 V – 240 V	6300 000.000
<i>realplex</i> ⁴ , 100 V – 130 V	950020113
<i>realplex</i> ⁴ , 200 V – 240 V	6302 000.008
Accessories	
Operating instructions for Mastercycler ep <i>realplex</i> (english)	6302 900.378
Temperature-Validation System	5331 222.005
Filter 520 nm	6302 101.009
Filter 550 nm	6302 103.001
Filter 580 nm	6302 104.008
Filter 605 nm	6302 105.004
PC CAN Bus cable and USB/CAN adapter	6302 070.405
CAN Bus cable	6302 070.456
Software <i>realplex</i> 2.2	6302 861.500
Software <i>routine-realplex</i> 2.3	6302 862.506
Consumables	
twin.tec PCR Plate 96, semi-skirted, colorless (25 pieces)	0030 128.575
twin.tec PCR Plate 96, semi-skirted, yellow (25 pieces)	0030 128.583
twin.tec PCR Plate 96, semi-skirted, red (25 pieces)	0030 128.613
twin.tec PCR Plate 96, semi-skirted, green (25 pieces)	0030 128.591
twin.tec PCR Plate 96, semi-skirted, blue (25 pieces)	0030 128.605
twin.tec PCR Plate 96, skirted, colorless (25 pieces)	0030 128.648
twin.tec PCR Plate 96, skirted, yellow (25 pieces)	0030 128.656
twin.tec PCR Plate 96, skirted, red (25 pieces)	0030 128.680
twin.tec PCR Plate 96, skirted, green (25 pieces)	0030 128.664
twin.tec PCR Plate 96, skirted, blue (25 pieces)	0030 128.672
twin.tec <i>real-time</i> PCR Plate 96, semi-skirted, blue (25 pieces)	0030 132.530
twin.tec <i>real-time</i> PCR Plate 96, semi-skirted, white (25 pieces)	0030 132.548
twin.tec <i>real-time</i> PCR Plate 96, semi-skirted, black (25 pieces)	0030 132.556

11 Ordering information

11

Ordering information

Description	Order no.
twin.tec <i>real-time</i> PCR Plate 96, skirted, blue (25 pieces)	0030 132.505
twin.tec <i>real-time</i> PCR Plate 96, skirted, white (25 pieces)	0030 132.513
twin.tec <i>real-time</i> PCR Plate 96, skirted, black (25 pieces)	0030 132.521
Heat Sealing Film (10 x 10 pieces)	0030 127.650
Recommended aids	
Heat Sealer 230 V / 50 Hz	5390 000.024
PCR-Cooler 0.2 ml pink	3881 000.023
PCR-Cooler 0.2 ml blue	3881 000.031
PCR-Rack, sorted according to color (10 pieces)	0030 127.455
Liquid Handling Workstation [®] epMotion 5070 basic device incl. control panel, software, optical sensor, waste box, MMC and reader, operating instructions, 50/60 Hz, 200 – 240 V	5070 000.000

12 Ordering information North America

Please only use the original accessories recommended by ourselves. The use of spare parts or disposables other than those recommended by ourselves may impair the function and safety of the device. Any warranty and liability for damage thus caused will be excluded.

Description	Order no.
Mastercycler ep <i>realplex</i> ² , 100 V – 130 V	950020202
Mastercycler ep <i>realplex</i> ² S, 100 V – 130 V	950020211
Mastercycler ep <i>realplex</i> ⁴ , 100 V – 130 V	950020300
Mastercycler ep <i>realplex</i> ⁴ S, 100 V – 130 V	950020318
<i>realplex</i> ² , 100 V – 130 V	950020105
<i>realplex</i> ⁴ , 100 V – 130 V	950020113
Accessories	
Operating instructions for Mastercycler ep <i>realplex</i> (english)	on request
Temperature-Validation System	950008008
Filter 520 nm	on request
Filter 550 nm	on request
Filter 580 nm	on request
Filter 605 nm	on request
Software <i>realplex</i> 2.2	950019010
Software <i>routine-realplex</i> 2.3	950019015
Consumables	
twin.tec PCR Plate 96, semi-skirted, colorless (25 pieces)	951020303
twin.tec PCR Plate 96, semi-skirted, yellow (25 pieces)	951020320
twin.tec PCR Plate 96, semi-skirted, red (25 pieces)	951020389
twin.tec PCR Plate 96, semi-skirted, green (25 pieces)	951020346
twin.tec PCR Plate 96, semi-skirted, blue (25 pieces)	951020362
twin.tec PCR Plate 96, skirted, colorless (25 pieces)	951020401
twin.tec PCR Plate 96, skirted, yellow (25 pieces)	951020427
twin.tec PCR Plate 96, skirted, red (25 pieces)	951020486
twin.tec PCR Plate 96, skirted, green (25 pieces)	951020443
twin.tec PCR Plate 96, skirted, blue (25 pieces)	951020460
twin.tec <i>real-time</i> PCR Plate 96, semi-skirted, blue (25 pieces)	951022043
twin.tec <i>real-time</i> PCR Plate 96, semi-skirted, white (25 pieces)	951022055
twin.tec <i>real-time</i> PCR Plate 96, semi-skirted, black (25 pieces)	951022067
twin.tec <i>real-time</i> PCR Plate 96, skirted, blue (25 pieces)	951022003
twin.tec <i>real-time</i> PCR Plate 96, skirted, white (25 pieces)	951022015
twin.tec <i>real-time</i> PCR Plate 96, skirted, black (25 pieces)	951022027
Heat Sealing Film (10 x 10 pieces)	951020060
Recommended aids	
Heat Sealer 230 V / 50 Hz	951023078
PCR-Cooler 0.2 ml pink	022510541
PCR-Cooler 0.2 ml blue	022510525
PCR-Rack, sorted according to color (10 pieces)	951010065
Liquid Handling Workstation® epMotion 5070	
basic device incl. control panel, software, optical sensor, waste box, MMC and reader, operating instructions, 50 / 60 Hz, 200 – 240 V	960000005

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The product named below fulfills the relevant fundamental requirements of the EC directives and standards listed. In the case of unauthorized modifications to the product or an unintended use this declaration becomes invalid.

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21.05.2007

Hamburg, Date:

Projektmanagement, Project Management:

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Vaudaux-Eppendorf AG
Eppendorf UK Limited
Eppendorf North America, Inc. · One Cantiague Road · P.O. Box 1019 · Westbury, N.Y. 11590-0207 · USA