Laborota 4000 / 4001 efficient

Laborota 4010 / 4011 digital

Laborota 4002 / 4003 control

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Instruction manual, p. 71

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Manual de instrucciones, p. 213

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CE

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Datum _____

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Table of contents

1	About this document	75
1.1	Release notes	75
1.2	About this manual.....	75
1.3	Reference documents	75
1.4	Icons and symbols.....	75
2	Basic safety instructions.....	77
2.1	General safety instructions.....	77
2.2	Intended use.....	77
2.3	Improper use	78
2.4	Use in potentially explosive environments	78
2.5	Responsibilities of the operator	78
2.6	Responsibilities of the operating personnel	78
2.7	Qualifications of the operating personnel.....	78
2.8	Safety-conscious work	78
2.9	Safety features of the instrument	79
2.10	Labels on instrument.....	79
3	Description of instrument.....	80
4	Setup and initial use	81
4.1	Scope of delivery.....	81
4.2	Transport	82
4.3	Setting up the instrument base	82
4.4	Initial operation	83
4.4.1	Attaching the heating bath.....	83
4.4.2	Filling the heating bath	83
4.4.3	Moving the heating bath	84
4.4.4	Connecting the instrument base	84
4.4.5	Using the lift.....	85
4.5	Installing glassware.....	86
4.5.1	Installing the condenser holder (G3–G6)	87
4.5.2	Installing the vapour tube	88
4.5.3	Installing the condenser	90

Table of contents

4.5.4	Mounting the vertical condenser (G3–G6) to the condenser holder	92
4.5.5	Installing the evaporating flask.....	93
4.5.6	Adjusting the angle of the evaporating flask	93
4.5.7	Adjusting the immersion depth of the evaporating flask	94
4.5.8	Releasing the evaporating flask from the vapour tube.....	95
4.5.9	Installing the inlet pipe.....	96
4.5.10	Installing the inlet hose.....	96
4.5.11	Installing the receiving flask	97
4.5.12	Connecting cooling water (not G5).....	97
4.5.13	Connecting the vacuum.....	97
4.6	Adding distillation material	102
4.7	Manual aeration.....	102
4.7.1	Installing/connecting accessories.....	102
5	Using the Laborota 4000 / 4001 efficient; 4010 / 4011 digital.....	106
5.1	Setting the rotating speed	106
5.2	Setting the heating bath temperature.....	107
5.3	Determining boiling temperature (Laborota 4010 / 4011, optional)	108
6	Using the Laborota 4002 / 4003 control	109
6.1	Setting the rotating speed	109
6.2	Setting the heating bath temperature.....	110
6.3	Determining the boiling temperature (with optional boiling temperature sensor)	111
6.4	Controlling the vacuum	111
6.4.2	p const mode	112
6.4.3	T auto mode	113
6.4.4	p auto mode	114
6.4.5	Interrupting evacuation.....	115
6.4.6	Applying the current value as the setpoint (p const mode only)	115
6.5	Time functions	115
6.5.1	Duration	115
6.5.2	Timer	116
6.6	Managing distillation parameters	116
6.6.1	Storing data records	117
6.6.2	Recalling data records.....	117
6.7	Ramp function	117

6.7.1	Programming a pressure ramp	117
6.7.2	Recalling a pressure ramp	118
6.7.3	Saving a pressure ramp	118
6.8	Auto Start Stop function	119
7	Faults and troubleshooting	120
8	Maintenance, cleaning, service	122
8.1	Maintenance	122
8.1.1	Calibrating the pressure sensor	122
8.1.2	Resetting the temperature limiter	123
8.1.3	Adjusting vacuum seal (glassware G6)	123
8.2	Replacing fuses	124
8.3	Cleaning	124
8.4	Service	125
9	Disassembly, storage, disposal	126
9.1	Disassembly	126
9.1.1	Removing glassware sets	126
9.1.2	Removing the Woulff bottle	126
9.1.3	Disconnecting cooling water/vacuum	126
9.1.4	Removing the heating bath	127
9.2	Storage	127
9.3	Disposal	127
10	Accessories, spare parts	128
10.1	Glass parts	128
10.2	Accessories	133
11	Appendix	134
11.1	Specifications	134
11.2	VAC senso T specifications	135
11.3	Solvent data	136
11.4	EU Declaration of Conformity	138
11.5	Warranty	138
11.6	Declaration on absence of hazards	139

1 About this document

1.1 Release notes

Version	Modification
1.02	10/2008

1.2 About this manual

This manual uses symbols and highlighting to help you find information quickly. Please read the notes on this topic in the next section.

Please read the safety instructions and warnings in this manual with special care. For safety instructions, please see Chapter 2; warnings can be found in the introductions of the individual chapters and before instructions related to actions.

Heidolph Instruments GmbH & Co KG has the copyright to all text and images in this manual.

1.3 Reference documents

For information on the Rotavac valve control and Rotavac vario control vacuum pumps, the Vac control automatic and the vacuum switch box, please see the operating manual 01-005-004-25.

1.4 Icons and symbols

Warnings

Symbol	Description
DANGER 	Immediate danger that will result in death or serious injury if not avoided.
WARNING 	Possible danger that may result in death or serious injury if not avoided.
CAUTION 	Possible danger that may result in minor injury or material damage if not avoided.
CAUTION	Possible danger that may result in material damage if not avoided.

Table 1-1: Warning indicators

About this document

Warnings have the following structure (here an example of the DANGER level):

DANGER

The type and source of the danger can be found here.

The possible consequences of not taking action to avoid the danger can be found here.

→ The measure to avoid the danger can be found here.

Symbols and highlighting

Symbol	Use	Explanation
✓	Precondition	This must be fulfilled before you follow an instruction.
→	Instruction	Description of required action
1. 2.	Instruction, multiple steps	Actions must be performed in the stated order. Deviating from the order can result in damage to the instrument or accidents.
Outcome	Result	The result of the previous action is shown indented.
• -	Bullet points, two levels	Lists
(see chapter 2, page 77)	Cross reference	References to images, tables, other chapters or other manuals
Example	Switch designation	Highlighting of designations of switches or buttons
	Note	Important information for understanding the instrument or optimising operation

Table 1-2: Symbols and highlighting

Images Position numbers in images are shown in parentheses in the text. If an image number is not stated, the position number refers to the image or chart directly above the text. When referring to a different image or chart, the image number will be specified, e.g. (Fig. 4-1 (11)).

2 Basic safety instructions

This state-of-the-art rotary evaporator has been designed to comply with recognised safety regulations. Nevertheless, certain hazards must be considered when installing, using and maintaining the instrument.

→ Follow safety instructions and warnings.

The basic safety instructions in this chapter are supplemented by specific warnings throughout this manual. These warnings explain precisely how to act in order to protect yourself or other persons from injury and property from damage.

This manual is a component of the Laborota 4000 /4001 efficient, Laborota 4010 / 4011 digital, Laborota 4002 / 4003 control rotary evaporators.

→ Keep this manual in a convenient location.

→ Be sure to give this manual to subsequent owners of the instrument.

2.1 General safety instructions

The rotary evaporator may only be used

- in technically sound condition
- according to the intended use
- if the operator has the necessary awareness of safety and potential hazards
- if the instructions in this manual are observed

Malfunctions – especially those that impair the instrument's safety – must be corrected without delay.

2.2 Intended use

Laborota 4000 /4001 efficient, Laborota 4010 / 4011 digital, Laborota 4002 / 4003 are rotary evaporators for:

- evaporating solutions
- powder drying
- distillation
- synthesis
- cleanup of chemicals
- Soxhlet extraction

Any application other than those stated above shall be deemed inappropriate.

The manufacturer shall bear no liability for damage arising from such applications.

The risk shall be borne solely by the operator. Intended use includes following the instructions in this manual and observing the terms of inspection and maintenance (see chapter 8.1, page 122).

2.3 Improper use

- Do not apply excess pressure to the instrument.
- Do not use the instrument in a corrosive atmosphere.

2.4 Use in potentially explosive environments

- Do not use the instrument in potentially explosive environments. The instrument is not explosion-proof.

2.5 Responsibilities of the operator

- Keep the instrument in technically sound condition.
- Ensure that the instrument is only used by qualified personnel.
- Ensure that users have been instructed in safe and responsible laboratory work practices.
- Ensure that the rotary evaporator is set up in a suitable location.

2.6 Responsibilities of the operating personnel

- Ensure that the distillation material can be evaporated without danger and that the residues are not explosive.

2.7 Qualifications of the operating personnel

The rotary evaporator is designed for use by qualified chemical and technical specialists. It may only be used by persons that have been instructed in its proper use.

2.8 Safety-conscious work

Heidolph Instruments will not be held liable for injury and/or damage arising from improper use.

- Please observe the following regulations:
 - laboratory regulations
 - safety regulations
 - ordinance on hazardous substances
 - any other generally applicable regulations related to safety engineering and occupational medicine

2.9 Safety features of the instrument

Heating bath	<ul style="list-style-type: none"> • electronic and mechanical overheating protection • electronic temperature control • temperature protection switch with programmable response temperature (Laborota 4002 and 4003 only)
Basic instrument	<ul style="list-style-type: none"> • clamps to secure evaporating and receiving flasks • adjustable immersion depth of the evaporating flask • evaporating flask may be lifted out of the heating bath manually • lift motor overcurrent protection (Laborota 4001, 4011 and 4003 only) • thermal overheating protection of drive motor
Glassware	<ul style="list-style-type: none"> • borosilicate glass • threaded connector fittings • condenser holder for G3, G4, G5 and G6 glassware
Optional	<ul style="list-style-type: none"> • Surlyn coating of the glass parts (Safecoat for G5 glassware) • guard hood • safety shield

2.10 Labels on instrument


Label	Description
	hot surface
	identification plate
<p>Disconnect the power plug before opening the instrument. Double-pole/neutral fusing. Replace fuse with same the and rating as specified</p> <p>  (2x) T 6,3A H 250VAC  </p>	fuse label (underside of instrument)

Table 2-1: Labels on instrument

3 Description of instrument

The Laborota 4000 / 4001 efficient, 4010 / 4011 digital, 4002 / 4003 control rotary evaporators are designed for distillation and evaporation of solutions.


Fig. 3-1: Laborota 4003 control rotary evaporator with G3 glassware

- | | |
|--------------------------|-------------------|
| 1 drive with vapour tube | 5 control panel |
| 2 evaporating flask | 6 receiving flask |
| 3 heating bath | 7 condenser |
| 4 instrument base | |

The solution to be evaporated is in the evaporating flask (2). The solution is evaporated by applying a suitable temperature of the heating bath (3) and vacuum pressure as appropriate to the thermodynamic properties of the solution. The settings can be adjusted on the control panel (5) of the instrument base (4) and the heating bath controls (3). The drive (1) rotates the evaporating flask to reduce the danger of boiling delays and accelerate evaporation by increasing the surface area of the solution. The vapour passes through the vapour tube into the condenser (7), where it condenses and flows into the receiving flask (6).

4 Setup and initial use

4.1 Scope of delivery

	Product	Quantity	Order no. 230 / 240 V 50 / 60 Hz	Order no. 120 V 50 / 60 Hz
	Basic instrument, Laborota 4000 efficient	1	519-00000-00	519-00000-01
or	Basic instrument, Laborota 4001 efficient	1	519-10000-00	519-10000-01
or	Basic instrument, Laborota 4010 digital	1	519-40000-00	519-40000-01
or	Basic instrument, Laborota 4011 digital	1	519-50000-00	519-50000-01
or	Basic instrument, Laborota 4002 control	1	519-20000-00	519-20000-01
or	Basic instrument, Laborota 4003 control	1	519-30000-00	519-30000-01
	Heating bath, 4000 digit	1	517-01002-00	517-01002-01
or	Heating bath, 4000 control	1	517-01001-00	517-01001-01
	Power cable	1	14-007-003-81	14-007-045-10
	Connecting cable - heating bath	1	14-007-045-08	14-300-009-28
	Laborota safety kit (4000 / 4001 / 4002 / 4003)	1	11-300-009-06	11-300-009-08
	Safety kit (Laborota 4010 / 4011)	1	11-300-009-10	11-300-009-12
	Product	Quantity	Order no. NS 29/ 32	Order no. NS 24 / 40
	Glassware G1	1	513-00100-00	513-00140-00
or	Glassware G3	1	513-00300-00	513-00340-00
or	Glassware G4	1	513-00400-00	513-00440-00
or	Glassware G5	1	513-00500-00	513-00540-00
or	Glassware G6	1	513-00600-00	513-00640-00
or	Glassware G1, coated	1	513-00110-00	513-00150-00
or	Glassware G3, coated	1	513-00310-00	513-00350-00
or	Glassware G4, coated	1	513-00410-00	513-00450-00
or	Glassware G5, coated	1	513-00510-00	513-00550-00
or	Glassware G6, coated	1	513-00610-00	513-00650-00
Included with the instrument				
	Screw coupling	1	23-09-03-01-03	
	Tension spring	1	22-03-02-01-05	
	Seal, PTFE 26	1	23-30-01-01-30	
	Collet 26	1	23-30-01-02-89	
	Vapour tube screw coupling	1	11-300-005-19	
	Heating bath control cable (LABOROTA 4002 / 4003 only)	1	14-300-003-36	
	Instruction manual	1	01-005-004-24	
	Warranty/declaration on absence of hazards	1	01-006-002-58	

Table 4-1: Scope of delivery

Setup and initial use

- Check on delivery**
- ➔ Unpack the rotary evaporator and ensure that it is complete and undamaged.
 - ➔ In case of damage, please inform Heidolph Instruments service (see chapter 8.4, page 125).

4.2 Transport

CAUTION

Vibration and mechanical shocks!

Damage to the housing and mechanism of the instrument.

- ➔ Avoid vibration and mechanical shocks.
-

Moving the unpacked instrument

- ✓ Ensure that the heating bath is empty.
 - ✓ Remove all hose and cable connections to other instruments, the cooling water/vacuum connection, and unplug the instrument from the AC power socket.
-
- ➔ Hold the underside of the instrument base.

4.3 Setting up the instrument base

- Where?** Chemicals may be released into the atmosphere via the pump when handling hazardous materials and distilling with solvents.
- ➔ Ensure that an exhauster is available for the exhaust air of the pump.
-

Do not set up the rotary evaporator in a corrosive atmosphere.

- ➔ Choose a level, firm base.

4.4 Initial operation

4.4.1 Attaching the heating bath

- ✓ Place the instrument base in a suitable location.


Fig. 4-1: Heating bath

1. Place the feet of the heating bath in the tracks (1). Ensure that the controls are facing forward.
2. If necessary, adjust the height of the feet.
 - Loosen the feet by hand or with an open-ended wrench and adjust them until the heating bath stands firmly without wobbling.

4.4.2 Filling the heating bath

The heating bath can be filled with a variety of liquid heating media.

When using to de-ionized or distilled water:

- ➔ Add 0.2% borax ($\text{Na}_2\text{B}_4\text{O}_7 \cdot 10 \text{H}_2\text{O}$) to the water.

The minimum/maximum volumes of the heating bath refer to the liquid level with the evaporating flask immersed.

1. Select a liquid suitable for your application such as tap water, oil with a low viscosity (40cP) and flash point (check safety data sheet) $> 260 \text{ }^\circ\text{C}$, water-soluble polyethylene glycol.
2. Lower the evaporating flask.
3. Fill the heating bath with medium. The level must be between the minimum and maximum marks.

4.4.3 Moving the heating bath

When using large evaporating flasks or spacers between the evaporating flask and vapour tube, it becomes necessary to increase the distance between the heating bath and the drive.

- Holding the heating bath by the grip, move it on the tracks and position it as required (see Fig. 4-1, page 83).

For applications for which the standard range is not enough, the range of movement can be increased to 210 mm using the extension plate (see chapter 10.2, page 133).

4.4.4 Connecting the instrument base

CAUTION

AC voltage and rated AC voltage printed on identification plate do not match!

Damage due to overvoltage.

- Ensure that your local AC voltage corresponds to that stated on the identification plate.
-

CAUTION

Unintentional activation of the drive and heating bath!

Injury and burning of the hands.

- Ensure that the power of the instrument base and heating bath is switched off.
-

By default the instrument intended for use with 230/240 V is equipped with a EURO plug (DIN 49441 CEE 7/VII 10/ 16 A 250 V), while those intended for 120 V have standard US plugs (NEMA Pub. No. WDI.1961 ASA C 73.1. 1961, p. 8, 15 A 125 V).

To use the instrument in a country with a different plug standard, use an approved adapter or have the plug replaced by a qualified electrician.

The plug sockets are located on the rear of the instrument.

- ✓ Instrument is switched off.
- ✓ Heating is installed on the instrument base.


Fig. 4-2: Connections

- | | |
|--|--|
| 1 heating bath connector cable | 4 temperature sensor/boiling temperature sensor (Laborota 4010 / 4011 / 4002 / 4003) |
| 2 power connection | 5 heating bath control cable (Laborota 4002 / 4003) |
| 3 remote control (Laborota 4001 / 4011 / 4003) | 6 VAC senso T / Rotavac vario control (Laborota 4002 / 4003) |

1. Insert the heating bath connector cable to the heating bath and the appropriate socket of the instrument base.

Laborota 4002 / 4003 control:

2. Connect the bayonet plug of the heating bath control cable (5-pin) to the heating bath and instrument base:
 - Grasp the plug by the knurled edge and turn it counter-clockwise to the stop.
 - Press the plug into the socket in such a manner that the locking pin fits into the groove.
 - Turn the knurled ring clockwise until the plug clicks into place.

The bayonet plug is now locked onto the connector.

3. Plug the power cable into the power input socket on the rear of the instrument.
4. Plug the power cable into an AC outlet.

4.4.5 Using the lift

The height of the drive can be adjusted in order to control the immersion depth of the evaporating flask.

The flask may not contact the bottom or side of the heating bath while immersed.

Manual lift Laborota 4000 efficient, 4010 digital, 4002 control (manual lift)


Fig. 4-3: Manual lift

1. Press the manual lift (1) downward and move left (lower) or right (higher) while pressed.
2. Release the manual lift (1) at the desired position.

The lift will be locked at the new position when released.

Motorised lift Laborota 4001 efficient, 4011 digital, 4003 control (motorised lift)

- ✓ Instrument base is installed.
- ✓ (for Laborota 4003 control) None of the **Set** modes are active (see chapter 6, page 109).
- ➔ Press the appropriate arrow button (up/down) (see Fig. 6-1, page 109 (9)) and release at the desired height.

The drive will remain at the desired height.

In the Laborota 4003 control, the motorised lift will not work as long as the instrument is in one of the **Set** modes.

4.5 Installing glassware

The following section describes the installation of the glassware. The rotary evaporator can be equipped with five different sets of glassware.

- ✓ Instrument base is set up and installed.
- ✓ Water connection is available (not for glassware G5).
- ✓ Lift is in highest position (see chapter 4.4.5, page 85).

WARNING

Danger of glass breakage!

Risk of severe injury due to broken glass.

- Check glass parts for damage (cracks, chips, etc.)
- Use only glass parts that are in flawless condition.
- Handle with care.

CAUTION

Drive rotating unintentionally!

Danger of injury to hands.

- Ensure that the rotation is turned off.

4.5.1 Installing the condenser holder (G3–G6)

Glassware sets with vertical condensers are supported by condenser holders to prevent them from tipping accidentally.


Fig. 4-4: *Installing the condenser holder*

CAUTION

Drive rotating unintentionally!

Danger of injury to hands.

- Ensure that the rotation is turned off.

1. Fasten stand rod (1) to drive (2) with Allen-head screw.

4.5.2 Installing the vapour tube

The vapour tube connects the evaporating flask to the condenser and permits the rotation of the evaporating flask by the drive.

WARNING

Danger of glass breakage!

Risk of severe injury due to broken glass.

- Check glass parts for damage (cracks, chips, etc.)
 - Use only glass parts that are in flawless condition.
 - Handle with care.
-

CAUTION

Drive rotating unintentionally!

Danger of injury to hands.

- Ensure that the rotation is turned off.
-

1. Remove vapour tube from its packaging.


Fig. 4-5: *Installing the vapour tube*

2. Remove the screw coupling (1) on the drive, take out the tension spring (2) and PTFE seal (3).


Fig. 4-6: Locking button

3. Press and hold the locking button (Fig. 4-6 (1)).
4. Unscrew the vapour tube screw coupling (Fig. 4-5 (6, 7)) from the drive.
5. Remove collet (Fig. 4-5 (5)).
6. Slide the vapour tube screw coupling (Fig. 4-5 (6, 7)) and collet (Fig. 4-5 (5)) onto the vapour tube with the ring facing forward (Fig. 4-5 (8)). Ensure that the ring clicks into place in the groove of the vapour tube (Fig. 4-5 (8)).
7. Insert the vapour tube (Fig. 4-5 (8)) and preassembled parts into the drive.
8. Press and hold the locking button (Fig. 4-6 (1)).
9. Hand-tighten the vapour tube (Fig. 4-5 (8)) in place at the drive button (Fig. 4-5 (4)) with the vapour tube screw coupling (Fig. 4-5 (6, 7)).


Fig. 4-7: Placing the PTFE seal on the vapour tube

Do not crease the PTFE seal or damage it with sharp fingernails.

10. Slide the PTFE seal (1) with the sealing lip forward (**Drive Side** label facing drive) onto the vapour tube (2).

4.5.3 Installing the condenser

- ✓ Screw coupling is removed from the drive head.
- ✓ Tension spring is removed from drive head.
- ✓ Vapour tube is installed in drive.
- ✓ PTFE seal is installed on vapour tube.

WARNING**Danger of glass breakage!**

Risk of severe injury due to broken glass.

- Check glass parts for damage (cracks, chips, etc.)
- Use only glass parts that are in flawless condition.
- Handle with care.

CAUTION**Drive rotating unintentionally!**

Danger of injury to hands.

- Ensure that the rotation is turned off.


Fig. 4-8: Installing glassware (here: glassware G1)

1. Version glassware sets G1 and G3:

- Slide screw coupling (2) over the condenser flange.
- Slide tension spring (3) over the condenser flange.
- Fasten condenser (1) to the drive (4) by tightening the coupling hand-tight.

2. Version **glassware G4**:

- Slide screw coupling (2) over the junction flange.
- Slide tension spring (3) over the junction flange.
- Fasten junction to the drive (4) by tightening the coupling (2) hand-tight.
- Insert condenser in the upper receptacle of the junction.

3. Version **glassware G5 (cooling water independent)**:

- Slide screw coupling (2) over the cooling trap jacket flange.
- Slide tension spring (3) over the cooling trap jacket flange.
- Fasten cooling trap jacket to the drive by tightening the coupling hand-tight.
- Insert the seal into the sealing groove of the cooling trap.
- Insert PTFE centring ring in the cooling trap and install the cooling trap insert.

The vacuum causes the seal to be pressed into place and the system is sealed hermetically by the inserted seal.

4. Version **glassware G6** (condenser for reflux distillation):


Fig. 4-9: Glassware G6

- Slide screw coupling (2) over the condenser flange.
- Slide tension spring (3) over the condenser flange.
- Fasten the screw coupling G6 (6) on the drive head (Fig. 4-8 (5)).
- Insert connector part G6 (5) with the seal G6 (4) in the screw coupling G6 (6).
- Fasten junction (1) to the screw coupling G6 (6) using the screw coupling (2).
- Insert condenser in the upper receptacle of the junction (1).

4.5.4 Mounting the vertical condenser (G3–G6) to the condenser holder

This step is not required for horizontal condensers (G1).

✓ Stand rod is installed.


Fig. 4-10: Condenser holder

Version **glassware G3, G4 and G6**:

1. Slide cross sleeve (3) onto the stand rod (1) and fix it with the clamping screw.
2. Place the condenser clamp (2) around the condenser and attach it to the stand rod with the clamping screws.
 - Adjust the stand rod to the angle of the glassware.
 - Fix the stand rod with the Allen screw (see also chapter 4.5.1, page 87).

1. Version **glassware G5** (cooling water independent):

- Insert the complete cooling trap holder into the transverse hole of the cross sleeve (2).
- Lock into place with the clamping screw.
- Adjust the stand rod (1) to the angle of the glassware.
- Place the cooling trap holder on the cooling trap jacket.
- Fix the cooling trap jacket with the clamping strap of the cooling trap holder.
- Fix the stand rod (1) with the Allen screw (see also chapter 4.5.1, page 87).

4.5.5 Installing the evaporating flask

WARNING

Danger of glass breakage!

Risk of severe injury due to broken glass.

- Check glass parts for damage (cracks, chips, etc.)
 - Use only glass parts that are in flawless condition.
 - Handle with care.
-

CAUTION

Drive rotating unintentionally!

Danger of injury to hands.

- Ensure that the rotation is turned off.
-

- Install the evaporating flask on the vapour tube with a conical-joint clamp.
-

4.5.6 Adjusting the angle of the evaporating flask

The angle of the evaporating flask can be adjusted by turning the drive.

CAUTION

Drive rotating unintentionally!

Danger of injury to hands.

- Ensure that the rotation is turned off.
-

WARNING

Danger of glass breakage!

Risk of severe injury due to broken glass.

- Check glass parts for damage (cracks, chips, etc.)
 - Use only glass parts that are in flawless condition.
 - Handle with care.
-

1. Grasp the condenser with your right hand.
-


Fig. 4-11: Angle and immersion depth of evaporating flask

2. With your left hand, turn the grey knob (1) at least one full turn anticlockwise toward **unlock**.
3. Adjust the angle of the evaporating flask by carefully moving the condenser.
4. With your left hand, turn the grey knob (1) toward **lock**.

4.5.7 Adjusting the immersion depth of the evaporating flask

CAUTION

Drive rotating unintentionally!

Danger of injury to hands.

→ Ensure that the rotation is turned off.

We recommend locking the immersion depth to ensure a consistent depth of the evaporating flask in the heating bath.

✓ Vertical stop (Fig. 4-24 (1)) in top position.

1. Version **without guard hood**:

- Turn black knob (Fig. 4-11 (2)) anticlockwise toward **unlock**.

The vertical stop (Fig. 4-24 (1)) will drop down.

- Move the lift down until the desired immersion depth has been reached.

2. Turn the black knob clockwise toward **lock**.

The immersion depth of the evaporating flask is now locked.

3. Version **with guard hood**:

- Move the lift down until the desired immersion depth has been reached.
- Hold the guard hood with your hand.
- Turn black knob (Fig. 4-11 (2)) anticlockwise toward **unlock**.

4. Guide the guard hood with your hand.
5. Turn the black knob clockwise toward **lock**.

The immersion depth of the evaporating flask is now locked.

4.5.8 Releasing the evaporating flask from the vapour tube

WARNING

Danger of glass breakage!

Risk of severe injury due to broken glass.

- Check glass parts for damage (cracks, chips, etc.)
 - Use only glass parts that are in flawless condition.
 - Handle with care.
-

CAUTION

Drive rotating unintentionally!

Danger of injury to hands.

- Ensure that the rotation is turned off.
-

- ✓ The evaporating flask must be cool.

1. Raise the evaporating flask from the bath (see chapter 4.4.5, page 85)
2. Hold the evaporating flask and release the clip.
3. Hold the locking button (see chapter 4.5.2, page 88)
4. Turn the nut on the vapour tube (Fig. 4-5 (7)) anticlockwise until the evaporating flask is released.
5. Let go of the locking button.

4.5.9 Installing the inlet pipe

The inlet pipe consists of:

- cock plug
- PTFE hose
- drain disk


Fig. 4-12: Installing the inlet pipe

1. Shorten the PTFE hose (1) to the desired length.
2. Grease the cock plug.
3. Position the drain disk (2) so that the condensed liquid can not flow back into the evaporating flask via the PTFE hose (1).
4. Insert the inlet pipe.

4.5.10 Installing the inlet hose


Fig. 4-13: Installing the inlet hose

1. Guide the inlet hose (4) through the cap (3) and seal (2) at the inlet tap (1).
2. Turn the cap (3) to the right.

The inlet hose (4) is not installed on the inlet tap (1).

4.5.11 Installing the receiving flask

1. Attach the receiving flask to the condenser/cooling trap jacket with spherical joint clip KS 35.
2. Secure with knurled screw.

The receiving flask is now installed on the condenser/cooling trap jacket/junction.

4.5.12 Connecting cooling water (not G5)

- ➔ Use a suitable hose with an inside diameter of 7–8 mm (e.g. hose set (see Accessories chapter 10.2, page 133))


Fig. 4-14: Connecting cooling water

1. Push the hose (cooling water inlet) (1) over the hose couplings (2) of the screw connection caps. Do this with the cooling water outlet (3) as well.
2. Secure the hoses with hose clips.

The cooling water inlet and outlet connections are now established.

4.5.13 Connecting the vacuum

Woulff bottle with vacuum valve

The Woulff bottle is connected to the vacuum valve (not Laborota 4002 / 4003 with Rotavac vario control).

WARNING

Imploding glass vessels!

Risk of severe injury due to broken glass.

- ➔ Check glass parts for damage (cracks, chips, etc.)
- ➔ Use only glass parts that are in flawless condition.


Fig. 4-15: Connecting Woulff bottle and vacuum valve

1. Remove screw connector cap (4) from Woulff bottle, take out hose clamp ring (5).
2. Insert included PTFE hose section (3) in cap nut (2) and tighten cap nut (2).
3. Slide screw connector cap (4) onto PTFE hose section (3); also slide hose clamp ring (5) on.
4. Connect to Woulff bottle.
The vacuum valve and Woulff bottle are now connected.
5. Attach the connected Woulff bottle and vacuum valve to the instrument base with the mounting clips (6) and two screws (7).

Connecting vacuum to glassware

- ✓ Suitable vacuum source available.
- ✓ Suitable vacuum hose (7–8 mm) available.


Fig. 4-16: Vacuum connection

1. Connect vacuum hose (3) to hose coupling (2) and screw connector cap (1).

Vacuum connections Hose routing of the individual vacuum systems


Fig. 4-17: Laborota 4000/4001 efficient, Laborota 4010/4011 digital with Rotavac valve control and Vac control automatic

- | | |
|--------------------------------|-----------------------|
| 1 condenser | 6 vacuum valve cable |
| 2 Y-piece | 7 power connection |
| 3 vacuum valve | 8 vacuum connection |
| 4 exhaust pump | 9 vacuum hose |
| 5 condensate cooler (optional) | 10 aeration/inert gas |


Fig. 4-18: Laborota 4000/4001 efficient, Laborota 4010/4011 digital with Rotavac valve control and Vac control automatic and Woulff bottle

- | | |
|--------------------------------|-----------------------|
| 1 condenser | 6 vacuum valve cable |
| 2 Woulff bottle | 7 power connection |
| 3 vacuum valve | 8 vacuum connection |
| 4 exhaust pump | 9 vacuum hose |
| 5 condensate cooler (optional) | 10 aeration/inert gas |


Fig. 4-19: Laborota 4002 / 4003 with Rotavac valve control and Vac Senso T

- | | |
|--------------------------------|----------------------|
| 1 condenser | 7 aeration/inert gas |
| 2 Y-piece | 8 vacuum connection |
| 3 vacuum valve | 9 power connection |
| 4 exhaust pump | 10 vacuum hose |
| 5 condensate cooler (optional) | 11 data cable |
| 6 vacuum valve cable | |


Fig. 4-20: Laborota 4002 / 4003 with Rotavac valve control and Vac Senso T and with Woulff bottle

- | | |
|--------------------------------|----------------------|
| 1 condenser | 7 aeration/inert gas |
| 2 Woulff bottle | 8 vacuum connection |
| 3 vacuum valve | 9 power connection |
| 4 exhaust pump | 10 vacuum hose |
| 5 condensate cooler (optional) | 11 data cable |
| 6 vacuum valve cable | |


Fig. 4-21: Laborota 4002 / 4003 with Rotavac vario control

- | | |
|--------------------------------------|----------------------|
| 1 condenser | 6 aeration/inert gas |
| 2 Y-piece | 7 vacuum connection |
| 3 exhaust pump | 8 power connection |
| 4 condensate cooler (optional) | 9 vacuum hose |
| 5 Rotavac vario control power supply | 10 data cable |


Fig. 4-22: Laborota 4002 / 4003 with Rotavac vario control with Woulff bottle

- | | |
|--------------------------------------|----------------------|
| 1 condenser | 6 aeration/inert gas |
| 2 Woulff bottle | 7 vacuum connection |
| 3 exhaust pump | 8 power connection |
| 4 condensate cooler (optional) | 9 vacuum hose |
| 5 Rotavac vario control power supply | 10 data cable |

4.6 Adding distillation material

Additional distillation material can be added during distillation and while under vacuum via the inlet hose and pipe.


Fig. 4-23: Inlet hose

1. Insert the inlet hose (1) in the additional material.
2. Turn the cock plug (2) parallel to the inlet hose (1).

The additional distillation material will be drawn into the evaporating flask.

4.7 Manual aeration

The system can be aerated via the cock plug (Fig. 4-23 (2)) of the inlet pipe:

➔ Turn the cock plug (Fig. 4-23 (2)) of the inlet pipe slowly upward.

The system will be aerated.

4.7.1 Installing/connecting accessories

Guard hood


Fig. 4-24: Installing the guard hood

1. Attach the guard hood holder (2) to the vertical stop (1) with the two included Allen-head screws (3).
2. Release the clamping screw (4) and align the guard hood (5) to the heating bath by moving it horizontally.
3. Retighten the clamping screw.

The guard hood is now installed. It can be opened using the grip.

Safety shield


Fig. 4-25: Installing the safety shield

1. Hang the safety shield (1) on the edge of the heating bath (2) and position it so that it does not interfere with the opening of the guard hood (if present).
2. Turn eccentric screw (3) by 180° to clamp the safety shield to the heating bath

Boiling temperature sensor (Laborota 4010 / 4011 and 4002 / 4003)


Fig. 4-26: Boiling temperature sensor

The boiling temperature sensor can not be used when using the temperature sensor T auto.

1. Remove the screw cap on the condenser.
2. Insert the boiling temperature sensor in the condenser so that the white PTFE side of the seal (2) faces the glass thread.
3. Secure the boiling temperature sensor with the screw connector cap.
4. Connect the bayonet plug (1) to the socket on the rear of the instrument (see chapter 4.4.4, page 84).

Temperature sensor T auto (Laborota 4002 / 4003)


Fig. 4-27: Temperature sensor Tauto

Only in conjunction with glassware sets G3, G4 and G6.

Adjust the depth to which the temperature sensor T auto extends into the condenser so that the lower edge of the sensor is positioned at 2/3 of the condenser height.

1. Remove the screw cap at the vacuum connection.
2. Insert the temperature sensor T(auto) into the vacuum connection so that the white PTFE side (3) of the seal faces the glass thread.
3. Slide the temperature sensor T(auto) to adjust its depth in the condenser.
4. Lock the temperature sensor T(auto) in place with the screw connector cap (1).
5. Connect the bayonet plug (2) to the socket on the rear of the instrument (see chapter 4.4.4, page 84).

Remote control
(Laborota
4001 efficient,
4011 digital,
4003 control)


Fig. 4-28: Remote control

- Connect the bayonet plug of the remote control to the socket on the rear of the instrument (see chapter 4.4.4, page 84).

Raising and lowering the lift

1. Press the arrow buttons (3) on the remote control.

Starting vacuum/rotation (Laborota 4003 only):

1. Press the **Vac Start Stop** (2) or **Rot Start Stop** (1) buttons on the remote control.

VAC senso T /
Rotavac vario
control power
supply (Laborota
4002 / 4003 control)


Fig. 4-29: VAC senso T / ROTAVAC vario control power supply

A	VAC senso T	6	vacuum connector
B	ROTAVAC vario control	7	aeration/inert gas
1	main power switch	8	power connection
2	LED	9	Laborota 4002 / 4003 connection
3	locking screws	10	ROTAVAC vario control vacuum pump
4	VAC senso T / Rotavac vario control connection	11	vacuum valve
5	bolt	12	vacuum switch box

1. Remove the positioning screws (3) from the locking tubes on the back of the instrument.
2. Insert the VAC senso T / ROTAVAC vario control power supply and bolts (5) into the locking tubes.
3. Secure the power supply with the positioning screws (3).
4. Establish connections (4, 6, 8, 9, 10, 11, 12).

5 Using the Laborota 4000 / 4001 efficient; 4010 / 4011 digital

This chapter covers the basic steps in using the Laborota 4000 / 4001 efficient and Laborota 4010 / 4011 digital rotary evaporators.

5.1 Setting the rotating speed

CAUTION

Drive rotating unintentionally!

Danger of injury to hands.

→ Ensure that the rotation is turned off.

- ✓ Rotation knob is set to left stop.
- ✓ Heating bath power switch is off.
- ✓ Instrument is switched on.

Raising the rotating speed permits the distillation rate to be increased.


Fig. 5-1: Laborota 4000 / 4001 efficient; 4010 / 4011 digital control panel

- | | |
|---|-----------------------------------|
| 1 display (left)
(Laborota 4010 / 4011
digital only) | 3 instrument base power
switch |
| 2 display (right)
(Laborota 4010 / 4011
digital only) | 4 rotation knob |

1. Switch on the power switch of the instrument base (3).

The instrument base power switch will be lit green (3).

2. Set the Rotation knob (4) to the desired speed.

The flask will start rotating. The Laborota 4010 / 4011 digital shows the selected speed on its display (1).

5.2 Setting the heating bath temperature

- ✓ Heating bath knob is set to left stop.
- ✓ Heating bath is filled with heating medium.
- ✓ Rotary evaporator is operational and the power switch is on.

WARNING

Surface temperatures in excess of 60°C during heating bath operation!

Danger of burns.

- Do not touch the inside and upper edge of the heating bath, the evaporating flask and the heating bath medium.
 - Wear protective gloves when changing the evaporating flask.
-

CAUTION

Overheated heating bath!

Damage and changes to the appearance of the heating bath.

- Never use the heating bath without liquid.
-

CAUTION

Damage to the guard hood!

- Set a maximum setpoint of 120 °C for operation with a guard hood.
-

The temperature difference between the heating bath and the vapour should be at least 20 K to ensure an adequate distillation rate.

As a rule of thumb, doubling the temperature difference will double the distillation rate.

Use a suitable oil as a heating bath medium for setpoint temperatures in excess of 100°C (observe safety data sheets, see chapter 4.4.2, page 83).

The message **OIL** will appear on the display of the heating bath (4) for one second when setting temperatures higher than 100°C.


Fig. 5-2: Heating bath (Laborota 4010/4011 digital)

1. Switch on the heating bath with its power switch (3).
The heating bath power switch will be lit green (3).
2. Adjust the setpoint temperature with the heating bath knob (4).
The heating bath display (1) will display the setpoint temperature flashing (Laborota 4010/4011 digital only). After 5 seconds without changes, the heating bath display (1) will switch to the actual temperature (not flashing).
The heating bath signal lamp (2) indicates that the heating bath is warming up.

5.3 Determining boiling temperature (Laborota 4010 / 4011, optional)

- ✓ Instrument is switched off.
- ✓ The boiling temperature sensor (optional) is connected.

→ Switch the instrument on.

The boiling temperature will be shown on the display of the rotary evaporator when the boiling temperature sensor is connected.

6 Using the Laborota 4002 / 4003 control

This chapter covers the basic steps in using the Laborota 4002 / 4003 control rotary evaporators.

The parameters for the distillation are entered and monitored on the control panel.

The selected parameters are saved. The parameters last used will appear on the display after a restart or an extended period in which the instrument is not used.


Fig. 6-1: instrument base / heating bath control panel

6.1 Setting the rotating speed

- ✓ Instrument is switched on.
- ✓ Rotation is switched off (**Rot Start Stop** button (7)).

At the left-hand stop the speed is not zero, but corresponds to the minimum speed value.

- ➔ Switch off rotation (**Rot Start Stop** button (7)).

1. Press **Select** (15) button until **Act Speed** LED is lit.
2. Press **Rot Start Stop** (7) button.

The drive rotates at minimum speed when the Rotation knob (11) is turned all the way to the left.

3. Turn the Rotation knob (11) to set the desired speed.

The set speed is shown on the display (14).

6.2 Setting the heating bath temperature

- ✓ Heating bath is filled with heating medium.
- ✓ Rotary evaporator is switched on at the instrument base power switch (Fig. 6-1 (10)) and operational.

WARNING

Surface temperatures in excess of 60°C during heating bath operation!

Danger of burns.

- ➔ Do not touch the inside and upper edge of the heating bath, the evaporating flask and the heating bath medium.
 - ➔ Wear protective gloves when changing the evaporating flask.
-

CAUTION

Overheated heating bath!

Damage and changes to the appearance of the heating bath.

- ➔ Never use the heating bath without liquid.
-

CAUTION

Damage to the guard hood!

- ➔ Set a maximum setpoint of 120 °C for operation with a guard hood.
-

Use a suitable oil as a heating bath medium for setpoint temperatures in excess of 100°C (observe safety data sheets, see chapter 4.4.2, page 83).

The message **OIL** will appear on the display of the heating bath (Fig. 6-1 (8)) for one second when setting temperatures higher than 100°C.

1. Switch on the heating bath with its power switch (Fig. 6-1 (18)).
The heating bath power switch will be lit green (Fig. 6-1 (18)).
2. Use the heating bath knob (Fig. 6-1 (16)) to set the temperature at least 20°C higher than the required temperature.
3. Press **Select** (Fig. 6-1 (6)) until **Set Bath** LED is lit.
4. Use the arrow buttons (Fig. 6-1 (9)) to set the setpoint temperature.
The setpoint temperature is visible on the display (Fig. 6-1 (8)). After 5 seconds, the current actual temperature appears on the display and the **Act Bath** LED is lit.

The lit yellow heating bath signal lamp (Fig. 6-1 (17)) over the heating bath power switch (Fig. 6-1 (18)) indicates that the heating bath is warming up.

6.3 Determining the boiling temperature (with optional boiling temperature sensor)

Connecting a boiling temperature sensor during operation:

- Switch the instrument off and back on to activate the program section.

The boiling temperature will be shown on the display of the control panel when the boiling temperature sensor is connected.

- Press **Select** (Fig. 6-1 (6)) until **Act Vap** LED is lit.

The boiling temperature is visible on the display (Fig. 6-1 (8)).

6.4 Controlling the vacuum

Together with the VAC senso T vacuum controller and the vacuum valve, the Laborota 4002 / 4003 control is a unit that permits the vacuum to be controlled via the opening and closing of a valve.

Together with the ROTAVAC vario control vacuum pump, the Laborota 4002 / 4003 control is a unit that permits the vacuum to be controlled via the speed of the vacuum pump.

The aeration valve is always open when the vacuum controller is switched off.

WARNING

Imploding glass vessels!

Risk of severe injury due to broken glass.

- Check glass parts for damage (cracks, chips, etc.)
- Use only glass parts that are in flawless condition.

Switching the vacuum control on and off

- Press **Vac Start Stop** (Fig. 6-1 (13)) button.

The LED next to the **Vac Start Stop** (Fig. 6-1 (13)) button is lit when active.

Press the **Vac Start Stop** (Fig. 6-1 (13)) button again to switch vacuum control off.

6.4.2 p const mode

Use the **Mode** button to choose between three vacuum control modes. A LED shows the currently active mode.

In **p const** mode, the vacuum is held constant at a given value. This setting is especially suitable for all standard processes.

For distillation, that means:

- No temperature sensor required.
- The pressure must be adjusted manually during mixture distillation.
- Knowledge of thermodynamic material data is required.

Selecting p const mode

Select the setpoint vacuum so that the boiling point of the substance to be distilled is between the bath temperature and the temperature of the cooling medium. Ideally, the difference between the bath temperature and boiling temperature, and that of the boiling temperature and temperature of the cooling medium, should each be 20 K.

Setting a value for Δp is superfluous when using the ROTAVAC vario control. The display (Fig. 6-1 (14)) shows **Set Δp : ----**. Vacuum control is realised via the speed of the pump.

The value of the differential gap (hysteresis) Δp indicates how much the pressure may rise (in mbar) before the vacuum valve is opened.

High hysteresis leads to high loss rates due to increased gas transport.

➔ Set Δp so that the valve does not activate too often.

For low boilers we recommend Δp values of 5–10 mbar, for high boilers 1–5 mbar.

1. Press **Mode** (Fig. 6-1 (12)) until **Mode p const** LED is lit.
2. Press **Select** (Fig. 6-1 (15)) until **Set Vac** LED is lit.
3. Use the arrow buttons (Fig. 6-1 (9)) to set the required vacuum.
The vacuum setpoint is visible on the display (Fig. 6-1 (14)).

VAC Senso T:

4. Press **Select** (Fig. 6-1 (15)) until **Set Δp** LED is lit.
5. Use the arrow buttons (Fig. 6-1 (9)) to set the desired differential gap (hysteresis) Δp .

6.4.3 T auto mode

In **T auto** mode, the vacuum is controlled via the temperature sensor T auto (accessory).

This has the following advantages:

- Only small amounts of the solvent are drawn off by the vacuum pump.
- The boiling point can be found automatically and corrected as required for mixtures.
- Mixture distillation can be performed without manual adjustment.
- No knowledge of thermodynamic material data is required.

Selecting T auto mode

- ✓ Temperature sensor T auto (accessory) is installed (see chapter 4.7.1, page 102).
- ✓ Cooling medium flowing for 15 minutes.

1. Press **Mode** (Fig. 6-1 (12)) until **Mode t auto** LED is lit.
2. Press **Select** (Fig. 6-1 (6)) until **Act Vap** LED is lit.
3. Read off the Act-Vap initial temperature.
4. Press **Select** (Fig. 6-1 (6)) until **Set Vap** LED is lit.

The set setpoint temperature has to lie 2–15 K above the Act-Vap initial temperature.

Values for **Set Vap** only slightly higher than the constant temperature (2–15 K) result in distillation with optimal recovery rates. Lower values lead to slower distillation, higher rates to faster.

5. Use the arrow buttons (Fig. 6-1 (9)) to set the required temperature.
The setpoint temperature is visible on the display (Fig. 6-1 (8)).
After five seconds, **Endp** appears in the display (Fig. 6-1 (14)).

If distillation is not to be aborted after a certain pressure is reached, enter the value 1 for **Endp**.

6. Use the arrow buttons (Fig. 6-1 (9)) to set the desired final pressure at which distillation is to be aborted.
The selected final pressure is visible on the display (Fig. 6-1 (8)).

6.4.4 p auto mode

The **p auto** mode combines the **T auto** and **p const** vacuum modes. **T auto** is used to reach the boiling point of a solution and the pressure determined for the continuation of the distillation is used with **p const**.

This has the following advantages:

- The boiling point is found automatically.
- No knowledge of thermodynamic material data is required.

Selecting p auto mode

- ✓ Temperature sensor T auto (accessory) is installed (see chapter 4.7.1, page 102).
-

Setting a value for Δp is superfluous when using the ROTAVAC vario control. The display shows **Set Δp** : ----. Vacuum control is realised via the speed of the pump.

The value of the differential gap (hysteresis) Δp indicates how much the pressure may rise (in mbar) before the vacuum valve is opened.

High hysteresis leads to high loss rates due to increased gas transport.

➔ Set Δp so that the valve does not activate too often.

For low boilers we recommend Δp values of 5–10 mbar, for high boilers 1–5 mbar.

Values for **Set Vap** only slightly higher than the constant temperature (2–15 K) result in distillation with optimal recovery rates. Lower values lead to slower distillation, higher rates to faster.

1. Press **Mode** (Fig. 6-1 (12)) until **Mode p auto** LED is lit.
 2. Press **Select** (Fig. 6-1 (6)) until **Set Vap** LED is lit.
 3. Use the arrow buttons (Fig. 6-1 (9)) to set the required temperature.
The setpoint temperature is visible on the display (Fig. 6-1 (8)).
 4. Press **Select** (Fig. 6-1 (15)) until **Set Δp** LED is lit.
 5. Use the arrow buttons (Fig. 6-1 (9)) to set the desired differential gap (hysteresis) Δp .
-

After the start of distillation in vacuum control mode **p auto**, the **Mode p auto** LED is lit continuously. Once the boiling point has been determined, the **Mode p const** LED also flashes.

6.4.5 Interrupting evacuation

To quickly interrupt the evacuation, e.g. in case of boiling delay or frothing:

→ Press **Hold** (Fig. 6-1 (1)) button.

The vacuum valve closes or the vacuum pump stops (ROTAVAC vario control). The LED next to **Vac Start Stop** (Fig. 6-1 (13)) button flashes.

Reactivating the vacuum

→ Press **Hold** (Fig. 6-1 (1)) button.

The LED next to **Vac Start Stop** (Fig. 6-1 (13)) button once again flashes.

6.4.6 Applying the current value as the setpoint (p const mode only)

To apply the current value as the setpoint during evacuation:

→ Press **Accept** (Fig. 6-1 (2)) button.

The current value is used as the setpoint for **p const** mode.

6.5 Time functions

6.5.1 Duration

The duration function records the time from switching on to switching off the rotation.

Setting the duration function

1. Press **Select** (Fig. 6-1 (15)) until **Set Time** LED is lit.
2. Use the arrow buttons (Fig. 6-1 (9)) to set the display (Fig. 6-1 (14)) to zero.

Recording of the duration starts when rotation is started by pressing the **Rot Start Stop** (Fig. 6-1 (7)) or **Auto Start Stop** (Fig. 6-1 (3)) buttons.

Recording of the duration ends when rotation is stopped by pressing the **Rot Start Stop** (Fig. 6-1 (7)) or **Auto Start Stop** (Fig. 6-1 (3)) buttons.

The recorded time can be displayed by selecting **Set Time** within 30 seconds after the rotation stops.

To display the duration during distillation:

1. Press **Select** (Fig. 6-1 (15)) until **Set Time** LED is lit.

The duration of distillation appears in the display (Fig. 6-1 (14)) for five seconds.

6.5.2 Timer

Setting the timer function

The timer function stops distillation after a specified time.

1. Press **Select** (Fig. 6-1 (15)) until **Set Time** LED is lit.
2. Use the arrow buttons to set the required time.

The timer starts counting down when rotation is started by pressing the **Rot Start Stop** (Fig. 6-1 (7)) or **Auto Start Stop** (Fig. 6-1 (3)) buttons. The distillation is aborted when the set time has elapsed.

- The system is aerated.
- The rotation is switched off.
- The lift is raised (Laborota 4003 control only).

30 seconds after distillation has been aborted, the preset time will be restored to the memory.

The instrument responds with five acoustic signals and a flashing display when a time is pre-programmed and rotation is started with the **Rot Start Stop** (Fig. 6-1 (7)) or **Auto Start Stop** (Fig. 6-1 (3)) buttons.

To switch between duration and timer functions during operation:

1. Press **Select** (Fig. 6-1 (15)) until **Set Time** LED is lit.
2. Use the arrow buttons (Fig. 6-1 (9)) to set the required time for the timer function, or 0 to select the duration function.

6.6 Managing distillation parameters

Up to 9 different data records can be stored in memory.

Each data record contains the following parameters:

- Set Time
- Set Vac
- Set Δp
- Set Bath
- Set Vap
- Endp
- Mode
- Pressure ramp if applicable

6.6.1 Storing data records

- ✓ Vacuum control is not active.
- ✓ Parameters are set.
- 1. Press **Store** (Fig. 6-1 (2)) button.
The number of the next free memory location appears on the display (Fig. 6-1 (8)).
- 2. To overwrite a memory location, use the arrow buttons (Fig. 6-1 (9)) to select it.
Flashing numbers indicate that the memory location is in use.
If all memory locations are in use, ---. will appear on the display (Fig. 6-1 (8)).
- 3. Use the arrow buttons to select a memory location number.
- 4. Press **Store** (Fig. 6-1 (2)) button.
An acoustic signal indicates that the data record has been saved.

6.6.2 Recalling data records

- ✓ Vacuum control is not active.
- 1. Press **Recall** (Fig. 6-1 (1)) button.
The data record with the lowest memory location number is recalled.
- 2. Use the arrow buttons to select a memory location number.
- 3. Press **Recall** (Fig. 6-1 (1)) button.
An acoustic signal indicates that the values have been transferred to the instrument's main memory.

6.7 Ramp function

Use the ramp function to program a pressure ramp in order to apply a range of pressure values in the **p const** vacuum control mode.

6.7.1 Programming a pressure ramp

- ✓ p const mode is selected.
- ✓ Timer settings have been checked.

The starting point is always at the time value 0 and a pressure of 950 mbar.

1. Press **Ramp Set** (Fig. 6-1 (4)) button.
The display (8) will indicate **t-1** for the first time value.
2. Use the arrow buttons (Fig. 6-1 (9)) to select the required value (shown on the display (Fig. 6-1 (14))).
3. Press **Ramp Set** (Fig. 6-1 (4)) button.
The time value is saved. The display (Fig. 6-1 (8)) will indicate **p-1** for the first pressure value.
4. Use the arrow buttons (Fig. 6-1 (9)) to set the required value.
The pressure value is visible on the display (Fig. 6-1 (14)).
5. Set additional pressure ramp values if required. 9 time and pressure values are available (up to **t-9** and **p-9**).
6. Press **Ramp Active** (Fig. 6-1 (5)) button.

Programming is complete when the final pressure value has been entered. The LED next to the **Ramp Active** (Fig. 6-1 (5)) button indicates that the next distillation started in **p const** mode will use the programmed ramp.

6.7.2 Recalling a pressure ramp

To recall a stored pressure ramp:

- ➔ Select the appropriate data record (see chapter 6.6.2, page 117).

6.7.3 Saving a pressure ramp

- ✓ Vacuum control is not active.
- ✓ A pressure ramp is programmed.
- ✓ **Ramp Active** (Fig. 6-1 (5)) has been pressed (LED lit).

1. Press **Store** (Fig. 6-1 (2)) button.
The number of the next free memory location appears on the display (Fig. 6-1 (8)).
2. To overwrite a memory location, use the arrow buttons (Fig. 6-1 (9)) to select it.
Flashing numbers indicate that the memory location is in use.
If all memory locations are in use, ---. will appear on the display (Fig. 6-1 (8)).
3. Use the arrow buttons (Fig. 6-1 (9)) to select a memory location number.
4. Press **Store** (Fig. 6-1 (2)) button.

An acoustic signal indicates that the data record, including the pressure ramp values, has been saved.

The programmed pressure ramps can be edited at any time by changing the relevant parameters (see chapter 6.6.1, page 117, and see chapter 6.7.1, page 117).

6.8 Auto Start Stop function

Use the Auto Start Stop function to automatically start and stop distillation.

- ✓ Instrument is operational.
- ✓ Evaporating flask immersion depth has been set (see chapter 4.5.7, page 94).
- ✓ Desired settings have been made.

1. Press **Auto Start Stop** (Fig. 6-1 (3)) button.

Lift moves downward (Laborota 4003 only); rotation, vacuum control and timing start. The LED next to the **Auto Start Stop** (Fig. 6-1 (3)) button is lit.

If one function is already active, the remaining functions will be activated.

All distillation parameters can be changed while the Auto Start Stop function is active.

Stopping the Auto Start Stop function

1. Press **Auto Start Stop** (Fig. 6-1 (3)) button.

Lift moves upward (Laborota 4003 only); rotation, vacuum control and timing stop. The LED next to the **Auto Start Stop** (Fig. 6-1 (3)) button goes out.

7 Faults and troubleshooting

Fault/Message	Cause	Remedy
Instrument can not be switched on.	Power cable not connected to an AC outlet.	→ Plug power cable into an AC outlet.
	Fuse defective	→ Replace fuse (see chapter 8.2, page 124)
No heating function	Main power switch turned off	→ Turn on main power switch
	Fuse defective	→ Replace fuse (see chapter 8.2, page 124)
	Main power switch defective	→ Contact service
	Heating bath connector cable is not plugged in (LABOROTA 4002 / 4003 control only)	→ Plug in heating bath connector cable
	Heating bath heater defective	→ Contact service
	Set Bath temperature higher than heating bath control	→ Use heating bath control to set temperature 20°C higher than Set Bath temperature
	Temperature limiter tripped	→ If medium in heating bath: let heating bath cool and reset temperature limiter → If no medium in heating bath: contact service
Drive does not turn	Main power switch defective	→ Contact service
	Speed control is set to left stop.	→ Turn speed control to right
	Drive defective	→ Contact service
	Fuse defective	→ Replace fuse (see chapter 8.2, page 124)
Message in heating bath display: Er1	Heating bath defective	→ Contact service
Message in heating bath display: Er2	Heating bath defective	→ Contact service
Message in heating bath display: Er3	Heating bath heated without liquid heating medium	→ Switch instrument off → Check whether temperature limiter has tripped. If so, reset (see chapter 8.1.2, page 123) → Add liquid heating medium

Fault/Message	Cause	Remedy
Motorised lift does not work	Main power switch turned off	→ Turn on main power switch
	Main power switch defective	→ Contact service
	Fuse defective	→ Replace fuse (see chapter 8.2, page 124)
	Lift is in end position	→ Use other arrow button
	Mechanical or motor defect	→ Contact service
	Vertical stop not correctly set	→ Adjust height (see chapter 4.5.7, page 94)
No evacuation	VAC Senso T / Rotavac vario control power not switched on	→ Turn power switch on
	Fuse defective	→ Replace fuse (see chapter 8.2, page 124)
Inadequate vacuum	System not sealing properly	→ Check seals and connections → Check ground connections, grease if required
	Vacuum pump defective	→ Consult manufacturer's documentation of vacuum pump
Instrument switches off unexpectedly (Laborota 4002 / 4003)	Timer is programmed	→ Check or disable timer (see chapter 6.5.2, page 116)
	Fuse defective	→ Replace fuse (see chapter 8.2, page 124)
Instrument switches off unexpectedly (Laborota 4002 / 4003), Display shows SAFE	Emergency shutdown of Laborota 4002 / 4003 control when a pressure of 1100 mbar is detected	→ Press any button to restore the Laborota 4002 / 4003 control to its initial state → Calibrate vacuum sensor (see chapter 8.1.1, page 122)

Table 7-1: Troubleshooting

8 Maintenance, cleaning, service

8.1 Maintenance

The vapour tube and PTFE seal must be serviced regularly.

1. Remove the vapour tube and PTFE seal (see chapter 4.5.2, page 88) and clean.
2. Check PTFE seal for wear or damage and replace if necessary.

8.1.1 Calibrating the pressure sensor

The pressure sensors of the ROTAVAC vario control and VAC senso T can be calibrated.

- ✓ Vacuum measurement equipment is available.
- ✓ Instrument is connected.

1. Insert the vacuum measurement equipment near the condenser (e.g. using a Y-piece).
2. Switch on the vacuum pump and VAC senso T / ROTAVAC vario control.
3. Switch off the rotary evaporator at the main power switch.
4. Press the **Rot Start Stop** (Fig. 6-1 (7)) and **Select** (Fig. 6-1 (15)) buttons at the same time and switch the rotary evaporator on at the main power switch.
5. Hold the **Rot Start Stop** (Fig. 6-1 (7)) and **Select** (Fig. 6-1 (15)) buttons for approx. 4 seconds.

The measured atmospheric pressure is shown in the left-hand display (Fig. 6-1 (14)). The right-hand display shows (Fig. 6-1 (8)) **01**.

6. Use the arrow buttons (Fig. 6-1 (9)) to set the value measured by the vacuum measurement instrument.
7. Press **Rot Start Stop** (Fig. 6-1 (7)) button.
02 appears in the display. Automatic evacuation starts.
8. When the measured value of the vacuum measurement instrument remains constant, enter the measured value with the arrow buttons.
9. Press **Rot Start Stop** (Fig. 6-1 (7)) button.
10. Switch the instrument off with the main power switch (Fig. 6-1 (10)).

Calibration is now complete.

8.1.2 Resetting the temperature limiter

The temperature limiter is designed to turn off the heating bath in case the temperature controller fails and the temperature of the heating bath exceeds the maximum by approx. 10%.

✓ Heating bath is cool.

WARNING

Hot surfaces!

Danger of burns.

→ Let heating bath cool.

1. Remove the liquid from the heating bath.


Fig. 8-1: Temperature limiter

2. Press the temperature limiter (1), using a pointed object if necessary.

8.1.3 Adjusting vacuum seal (glassware G6)


Fig. 8-2: Adjusting the vacuum seal, glassware G6

The seal in the valve (1) can be adjusted with a screwdriver.

CAUTION

Danger of glass breakage!

Danger of cuts.

→ Take care not to turn the adjuster screw too far.

1. Tighten the screw of the valve in the junction G6 (2) to the stop.
2. Turn the adjuster screw clockwise until a white ring can be seen on the glass shaft.

8.2 Replacing fuses

The fuses are located on the bottom of the instrument.


Fig. 8-3: Fuse on bottom of instrument

CAUTION

Unintentional activation of the drive and heating bath!

Risk of severe injury due to electrical shock. Danger of burns to hands.

- Ensure that the power of the instrument base and heating bath is switched off.

1. Unplug the instrument from its AC power socket.
2. Remove glassware.
3. Allow heating bath liquid to cool and remove the heating bath from the instrument base.
4. Carefully tip the instrument onto its side and turn over.
5. Use a flat screwdriver to remove the defective fuse.
6. Insert a new fuse and tighten with the flat screwdriver.

The fuse is now installed. Return the instrument to its original state.

8.3 Cleaning

CAUTION

Damage due to corrosive cleaning products!

Damage to the surface of the instrument.

- Use only mild soap solution. Do not use chlorine bleach or other chlorine-based cleaning products. Do not use abrasive cleanser, ammonia, steel wool or other cleaning products with metallic components.

- Wipe the surface of the instrument with a cloth dampened in a mild soap solution

8.4 Service

Your instrument is not working?

1. Please call Heidolph Instruments or your authorised Heidolph Instruments dealer.

**Addresses and
telephone numbers**

In Germany,
Austria and
Switzerland:

Tel: 00800-HEIDOLPH
+49 (0) 91 22 99 20 69

Fax: +49 (0) 91 22 99 20 65

E-mail: sales@heidolph.de

2. After consultation with a Heidolph service technician:
 - Copy the declaration on absence of hazards in this manual and complete it
 - Pack the instrument securely and send it together with the declaration on absence of hazards to:

Heidolph Instruments GmbH & Co. KG

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Walpersdorfer Str. 12

91126 Schwabach, Germany

9 Disassembly, storage, disposal

9.1 Disassembly

9.1.1 Removing glassware sets

CAUTION

Drive rotating unintentionally!

Danger of injury to hands.

→ Ensure that the rotation is turned off.

WARNING

Hot surfaces!

Danger of burns.

→ Let heating bath and glassware cool.

CAUTION

Danger of glass breakage!

Danger of cuts.

→ Handle with care.

→ Remove glassware sets in reverse order of their installation (see chapter 4.5, page 86).

9.1.2 Removing the Woulff bottle

CAUTION

Danger of glass breakage!

Danger of cuts.

→ Handle with care.

→ Remove the Woulff bottle in the reverse order of its installation.

9.1.3 Disconnecting cooling water/vacuum

1. Ensure that the cooling water/vacuum supplies are shut off.
2. Disconnect the hoses from the condenser.

9.1.4 Removing the heating bath

WARNING

Hot surfaces!

Danger of burns.

→ Let heating bath cool.

→ Remove the heating bath in the reverse order of its installation.

9.2 Storage

Where? → Store the instrument in a dry environment.

How? → Store the instrument in its original packaging.

→ Close the packaging with adhesive tape.

9.3 Disposal

→ Dispose of the instrument in accordance with the applicable regulations of your country.

10 Accessories, spare parts

10.1 Glass parts

Glassware G1


Fig. 10-1: Glassware G1

Order designation	Position	Order no.	Delivery quantity
Screw cap GL 18	1	23-09-03-01-24	1
Flask clamp NS 29/32	2	515-40000-00	1
Flask clamp NS 24/40	2	515-40000-01	(1)
Evaporating flask 1000 ml, NS 29/32	3	514-74000-00	1
Evaporating flask, 1000 ml, NS 24/40	3	514-74000-05	(1)
Vapour tube, NS 29/32	4	514-00000-01	1
Vapour tube, NS 24/40	4	514-00020-03	(1)
Receiving flask 1000 ml, S 35/20	5	514-84000-00	1
Receiving flask 1000 ml, S 35/20, coated	5	514-84000-04	(1)
Flask clamp S 35/20	6	515-42000-00	1
Inlet hose (PTFE-) 3.5/4.5x600	7	23-30-01-04-72	1
Screw connector cap GL10 red	8	23-30-01-04-69	1
O-ring seal 3.2x2.5	9	23-08-06-03-26	1
Inlet pipe	10	514-51000-00	1
Screw connector cap GL 14	11	23-09-03-01-27	3
Hose coupling for screw connector cap GL 14	11	11-300-005-22	3
Condenser G1	12	514-00100-00	1
Condenser G1B, coated	12	514-00110-00	(1)

Table 10-1: Glassware G1

Glassware G3


Fig. 10-2: Glassware G3

Order designation	Position	Order no.	Delivery quantity
Screw cap GL 18	1	23-09-03-01-24	1
Vacuum attachment, NS 29/32	2	514-00001-00	1
Condenser G3	3	514-00300-00	1
Condenser G3B, coated	3	514-00310-00	(1)
Flask clamp NS 29/32	4	515-40000-00	1
Flask clamp NS 24/40	4	515-40000-01	(1)
Evaporating flask, 1000 ml, NS 29/32	5	514-74000-00	1
Evaporating flask, 1000 ml, NS 24/40	5	514-74000-05	(1)
Vapour tube, NS 29/32	6	514-00000-01	1
Vapour tube, NS 24/40	6	514-00020-03	(1)
Flask clamp S 35/20	7	515-42000-00	1
Receiving flask 1000 ml, S 35/20	8	514-84000-00	1
Receiving flask 1000 ml, S 35/20, coated	8	514-84000-02	(1)
Inlet hose (PTFE-) 3.5/4.5x600	9	23-30-01-04-72	1
Screw connector cap GL10 red	10	23-30-01-04-69	1
O-ring seal 3.2x2.5	11	23-08-06-03-26	1
Inlet pipe	12	514-51000-00	1
Screw connector cap GL 14	13	23-09-03-01-27	3
Hose coupling for screw connector cap GL 14	13	11-300-005-22	3

Table 10-2: Glassware G3

Glassware G4


Fig. 10-3: Glassware G4

Order designation	Position	Order no.	Delivery quantity
Screw cap GL 18	1	23-09-03-01-24	2
Vacuum attachment, NS 29/32	2	514-00001-00	1
Screw connector cap GL 14	3	23-09-03-01-27	3
Hose coupling for screw connector cap GL 14	3	11-300-005-22	3
Condenser G4	4	514-55000-00	1
Condenser G4B, coated	4	514-55000-02	(1)
Junction G4	5	514-00401-00	1
Junction G4B, coated	5	514-00411-00	(1)
Flask clamp NS 29/32	6	515-40000-00	1
Flask clamp NS 24/40	6	515-40000-01	(1)
Evaporating flask, 1000 ml, NS 29/32	7	514-74000-00	1
Evaporating flask, 1000 ml, NS 24/40	7	514-74000-05	(1)
Vapour tube, NS 29/32	8	514-00000-01	1
Vapour tube, NS 24/40	8	514-00020-03	(1)
Flask clamp S 35/20	9	515-42000-00	1
Receiving flask 1000 ml, S 35/20	10	514-84000-00	1
Receiving flask 1000 ml, S 35/20, coated	10	514-84000-02	(1)
Inlet hose (PTFE-) 3.5/4.5x600	11	23-30-01-04-72	1
Screw connector cap GL10 red	12	23-30-01-04-69	1
O-ring seal 3.2x2.5	13	23-08-06-03-26	1
Inlet pipe	14	514-51000-00	1

Table 10-3: Glassware G4

Glassware G5


Fig. 10-4: Glassware G5

Order designation	Position	Order no.	Delivery quantity
Cooling trap insert	1	514-00501-00	1
Centering ring	2	23-30-01-04-88	1
Seal G5 (silicone)	3	23-30-01-01-88	(1)
Seal G5	3	23-30-01-01-39	1
Flask clamp NS 29/32	4	515-40000-00	1
Flask clamp NS 24/40	4	515-40000-01	(1)
Evaporating flask, 1000 ml, NS 29/32	5	514-74000-00	1
Evaporating flask, 1000 ml, NS 24/40	5	514-74000-05	(1)
Vapour tube, NS 29/32	6	514-00000-01	1
Vapour tube, NS 24/40	6	514-00020-03	(1)
Flask clamp S 35/20	7	515-42000-00	1
Receiving flask 1000 ml, S 35/20	8	514-84000-00	1
Receiving flask 1000 ml, S 35/20, coated	8	514-84000-04	(1)
Inlet hose (PTFE-) 3.5/4.5x600	9	23-30-01-04-72	1
Screw connector cap GL10 red	10	23-30-01-04-69	1
O-ring seal 3.2x2.5	11	23-08-06-03-26	1
Inlet pipe	12	514-51000-00	1
Screw connector cap GL 14	13	23-09-03-01-27	3
Hose coupling for screw connector cap GL 14	13	11-300-005-22	3
Screw cap GL 18	14	23-09-03-01-24	1
Cock plug, NS 19/38	15	15-003-003-24	1
Holder, complete	16	591-00061-00	1
Cooling trap jacket	17	514-00500-00	1
Cooling trap jacket, coated	17	514-00510-00	(1)

Table 10-4: Glassware G5

Glassware G6


Fig. 10-5: Glassware G6

Order designation	Position	Order no.	Delivery quantity
Screw cap GL 18	1	23-09-03-01-24	2
Vacuum attachment, NS 29/32	2	514-00001-00	1
Screw connector cap GL 14	3	23-09-03-01-27	3
Hose coupling for screw connector cap GL 14	3	11-300-005-22	3
Valve	4	514-48000-00	1
Junction G6	5	514-00601-00	1
Junction G6B, coated	5	514-00611-00	(1)
Screw coupling, drive flange	6	-	-
Screw coupling G6	7	515-62000-00	1
Seal G6	8	23-30-01-01-35	1
Connector part G6	9	22-30-01-05-02	1
Flask clamp NS 29/32	10	515-40000-00	1
Flask clamp NS 24/40	10	515-40000-01	(1)
Evaporating flask, 1000 ml, NS 29/32	11	514-74000-00	1
Evaporating flask, 1000 ml, NS 24/40	11	514-74000-05	(1)
Vapour tube, NS 29/32	12	514-00000-01	1
Vapour tube, NS 24/40	12	514-00020-03	(1)
Receiving flask 1000 ml, S 35/20	13	514-84000-00	1
Receiving flask 1000 ml, S 35/20, coated	13	514-84000-02	(1)
Flask clamp S 35/20	14	515-42000-00	1
Inlet hose (PTFE-) 3.5/4.5x600	15	23-30-01-04-72	1
Screw connector cap GL10 red	16	23-30-01-04-69	1
O-ring seal 3.2x2.5	17	23-08-06-03-26	1
Inlet pipe	18	514-51000-00	1
Condenser G6	19	514-23000-00	1
Condenser G6B, coated	19	514-23000-02	(1)
Condenser holder		591-00060-00	1

Table 10-5: Glassware G6

10.2 Accessories

Order designation	Order no. 230 / 240 V 50 / 60 Hz	Order no. 120 V 50 / 60 Hz	Remarks
Guard hood	591-00010-00		
Safety shield	591-00020-00		
Boiling temperature sensor	591-00030-00		
Temperature sensor T auto	591-00040-00		
Remote control	591-00050-00		For motorised lift
Extension plate	591-00090-00		
Vacuum pump, ROTAVAC valve control	591-00130-00	591-00130-01	For all instruments
Vacuum pump, ROTAVAC vario control	591-00140-00	591-00140-01	For Laborota 4002 / 4003
VAC senso T vacuum regulator	591-00350-00	591-00350-01	For Laborota 4002 / 4003
VAC control automatic vacuum controller	591-00340-00	591-00340-01	For all except Laborota 4002 / 4003
Vacuum switch box for 3 consumers	591-00400-00	591-00400-01	
Woulff bottle	591-00071-00		
Vacuum valve	591-24000-00		
Condensate cooler for ROTAVAC valve control	591-00081-00		
Condensate cooler for ROTAVAC vario control	591-00082-00		
Heating bath liquid (up to 180 °C)	515-31000-00		
Hose set	591-35000-00		

Table 10-6: Accessories

11 Appendix

11.1 Specifications

	Laborota 4000 efficient	Laborota 4001 efficient	Laborota 4010 digital	Laborota 4011 digital	Laborota 4002 control	Laborota 4003 control
Dimensions (WxDxH) Drive in lowest position, without glassware	396 x 440 x 360 mm footprint with lift lever	396 x 440 x 360 mm footprint	396 x 440 x 360 mm footprint with lift lever	396 x 440 x 360 mm footprint	396 x 440 x 360 mm footprint with lift lever	396 x 440 x 360 mm footprint
Connector cable	2 m					
Power consumption	1400 W/ 1500 VA	1430 W/ 1500 VA	1410 W/ 1500 VA	1440 W/ 1500 VA	1400 W/ 1500 VA	1430 W/ 1500 VA
Voltage	230 V / 50 / 60 Hz or 120 V / 50 / 60 Hz					
Drive	capacitor motor with electronic speed control					
Lift drive	manual	motor	manual	motor	manual	motor
Bath temperature display	digital					
Boiling temperature display (with optional boiling temperature sensor)	–	–	digital	digital	digital	digital
Speed range	20–270 rpm					
Speed display	scale	scale	digital	digital	digital	digital
Heating bath diameter	240 mm					
Weight without glassware	17 kg	18 kg	17 kg	18 kg	19 kg	20 kg
Heater rating	1300 W					
Lifting speed	–	approx. 18 mm/s	–	approx. 18 mm/s	–	approx. 18 mm/s
Lift distance	140 mm					
Cooling area	1200 cm ² , glassware G3					
Heating bath material	V4A (1.4404) AISI 316L					
Heating bath range	20 - 180°C					
Heating bath regulator precision	±1 K					
Heating bath temperature regulation	electronic				μ processor	
Protection class	IP 20					
Evaporation performance	approx. 1200 ml H ₂ O/h					
Permissible environmental conditions	0–40 °C at 80 % rel. humidity					
Remote control (optional)	–	yes	–	yes	–	yes
Fuses	120 V~/ 50/60 Hz T 12.5 A H 250 V~					
	230 V~/ 50/60 Hz T 6.3 A H 250 V~					

Table 11-1: Specifications

Laborota 4002 / 4003 control only:

- Aeration valve control with VAC senso T or ROTAVAC vario control accessories
- Timer
- Vacuum controller integrated with VAC senso T or ROTAVAC vario control
- Vacuum control mode p auto with optional VAC senso T or ROTAVAC vario control
- Vacuum control mode p const with optional VAC senso T or ROTAVAC vario control
- Vacuum control mode T auto with optional VAC senso T or ROTAVAC vario control and T auto temperature sensor
- Adjustable heating bath limiter

11.2 VAC senso T specifications

	Order no. 591-00350-00	Order no. 591-00350-01
Voltage	230 V / 50 / 60 Hz	120 V / 50 / 60 Hz
Power consumption	24 VA	24 VA
Fuse (2x)	T 0.125 A H 250 V~	T 0.25 A H 250 V~

Table 11-2: VAC senso T specifications

11.3 Solvent data

The nomogram shows the relationship between pressure and boiling temperature of a selection of solvents.

The difference between the temperature of the vapour and the cooling medium should be 20 K to ensure adequate condensation.

The temperature difference between the heating bath and the vapour should be at least 20 K to ensure an adequate distillation rate.

As a rule of thumb, doubling the temperature difference will double the distillation rate.


Fig. 11-1: Nomogram

Conversion Torr to mbar: [mmHg] \approx 3/4 [mbar]

Solvent data

Solvent	Empirical formula	MW [g/mol]	Boil. [°C]	ΔH_{vap} [J/g]	Vacuum for boil.at 40 °C	
					[mbar]	[mm(Hg)]
Acetone	C ₃ H ₆ O	58.08	56.5	550	556	387
Acetonitrile	C ₂ H ₃ N	41.05	81.8	833	230	173
Benzene	C ₆ H ₆	78.11	80.1	549	236	177
n-butanol (butyl alcohol)	C ₄ H ₁₀ O	74.12	117.5	619	25	19
Tert. butanol (tert. butyl alcohol)	C ₄ H ₁₀ O	74.12	82.9	588	130	98
2-butanone (methyl ethyl ketone)	C ₄ H ₈ O	72.11	79.6	473	243	182
Chlorobenzene	C ₆ H ₅ Cl	112.60	132.2	375	36	27
Cyclohexane	C ₆ H ₁₂	84.16	80.7	389	235	176
1,2-dichloroethane	C ₂ H ₄ Cl ₂	98.96	82.4	336	210	158
1,2-dichloroethylene (cis)	C ₂ H ₂ Cl ₂	96.94	59.0	320	479	134
1,2-dichloroethylene (trans)	C ₂ H ₂ Cl ₂	96.94	47.8	313	751	563
Dichloromethane (methylene dichloride)	CH ₂ Cl ₂	84.93	40.7	373	atm.	atm.

Solvent	Empirical formula	MW [g/mol]	Boil. [°C]	ΔH_{vap} [J/g]	Vacuum for boil.at 40 °C	
Diethyl ether	C ₄ H ₁₀ O	74.12	34.6	392	atm.	atm.
Diisopropyl ether	C ₆ H ₁₄ O	102.20	67.5	318	375	281
Dimethyl formamide	C ₃ H ₇ NO	73.09	153.0		11	8
1,4-dioxan	C ₄ H ₈ O ₂	88.11	101.1	406	107	80
Ethanol	C ₂ H ₆ O	46.07	78.4	879	175	131
Ethyl acetate	C ₄ H ₈ O ₂	88.11	77.1	394	240	180
Heptane	C ₇ H ₁₆	85.09	98.4	439	120	90
Hexane	C ₆ H ₁₄	86.18	68.7	370	335	251
Methanol	CH ₄ O	32.04	64.7	1225	337	253
3-methyl-1-butanol (isoamyl alcohol)	C ₅ H ₁₂ O	88.15	130.6	593	14	11
Pentachloroethane	C ₂ HCl ₅	202.30	160.5	203	13	10
Pentane	C ₅ H ₁₂	72.15	36.1	382	atm.	atm.
n-pentanol (amyl alcohol)	C ₅ H ₁₂ O	88.15	137.8	593	11	8
1-propanol (n-propyl alcohol)	C ₃ H ₈ O	60.10	97.8	787	67	50
2-propanol (isopropyl alcohol)	C ₃ H ₈ O	60.10	82.5	701	137	103
1,1,2,2-tetrachloroethane	C ₂ H ₂ Cl ₄	167.90	145.9	247	35	26
Tetrachloroethylene	C ₂ Cl ₄	165.80	120.8	233	53	40
Tetrachloromethane (carbon tetrachloride)	CCl ₄	153.80	76.7	225	271	203
Tetrahydrofuran	C ₄ H ₈ O	72.11	66.0	–	357	268
Toluene	C ₇ H ₈	92.14	110.6	425	77	58
1,1,1-trichloroethane	C ₂ H ₃ Cl ₃	133.40	74.1	251	300	225
Trichloroethylene	C ₂ HCl ₃	131.40	86.7	265	183	137
Trichloromethane (chloroform)	CHCl ₃	119.40	61.3	263	474	356
Water	H ₂ O	18.02	100.0	2259	72	54
Xylene (isomer mixture)	C ₈ H ₁₀	106.20	137–143	390	25	19

Tabelle 11-3: Solvent data

11.4 EU Declaration of Conformity

We, Heidolph Instruments GmbH & Co KG, hereby declare that the products Laborota 4000 / 4001 efficient, 4010 / 4011 digital and 4002 / 4003 control comply with the following standards and normative documents:

EMC directives:

- EN 61326: 1997 + A1:1998 + A2:2001+ A3 2003
- EN 61000-3-2: 2000
- EN 61000-3-3: 1995 + 1997 + A1:2001
- EN 61000-4-2:1995 +A1:1998+A2:2001
- EN 61000-4-3:2002 +A1:2002
- EN 61000-4-4:1995 +A1:2001 + A2:2001
- EN 61000-4-5:1995 +A1:2001
- EN 61000-4-6:1996 +A1:2001
- EN 61000-4-11:1994 + A1:2001

Low-voltage directive:

- EN 61010-1 + EN 61010-2-010

Modifications to the instrument without our prior approval will void this Declaration.

Date: 26.02.2007

Signature

Name of the signatory

Günter Vollgold

11.5 Warranty

Heidolph Instruments provides a three-year warranty on the products described herein (with the exception of glass and consumable parts) subject to registration using the enclosed warranty card or via the Internet (www.heidolph.com). The warranty period begins at the time of registration. Without registration, the serial number of the instrument shall apply. This warranty covers defects in materials and workmanship. Transit damage is excluded from this warranty.

To obtain such warranty service, contact Heidolph Instruments (phone: +49 – 9122 - 9920-69) or your local Heidolph Instruments Dealer. If defects in material or workmanship are found, your item will be repaired or replaced at no charge.

Misuse, abuse, neglect or improper installation are not covered by this warranty promise.

11.6 Declaration on absence of hazards

Should your instrument require repairs, please complete this declaration and mail or fax it to Heidolph Instruments GmbH & Co KG (see chapter 8.4, page 125).

1. Information about the instrument

Type code _____
 Serial number _____
 Reason for sending _____

2. Was the instrument cleaned, decontaminated or disinfected?

yes _____ no _____

3. Is the instrument in a state that would not represent a health hazard to service personnel?

yes _____ no _____

If not, to what substances has the instrument been exposed?

4. Information about sender

Surname, first name _____
 Company _____
 Department _____
 Street address _____
 City, postcode/ZIP _____
 Country _____
 Telephone _____
 E-mail _____

5. Legally binding declaration

The client hereby accepts liability for damages suffered by the client arising from incomplete or incorrect information.

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