
ROTAVAPOR

R-134 R-144

96530 A4 4000 0193 FAE

Hinweise

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Notes

- *No liability is accepted for damage and breakdowns resulting from a failure to observe the operating instructions.*
- *Illustrations and information in these operating instructions are subject to technical alterations which are necessary for improving an apparatus and its operation.*
- *The copyright of the texts, drawings and diagrams in the operating instructions is owned by BÜCHI LABORATORIUMS TECHNIK AG.*

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Contents	Page	Annex	Page
Setting up the apparatus			
1. The rotary evaporator	6	11. Tables	39
2. Technical data	7	11.1 Choosing the distillation temperatures	39
3. Installation	9	11.2 Solvent table	40
3.1 Notes and warnings	9	12. Troubleshooting	42
3.2 Standard equipment	11	12.1 Error messages	42
3.3 Installing the bath	13	12.2 Other malfunctions	42
3.4 Flanged screw couplings	13	13. Cleaning	45
3.5 Installing the seals	13	14. Assembling and dismantling	45
3.6 Bath immersion angle	14	14.1 Converting the sealing system	45
3.7 Setting up the stand	14	14.2 Replacing the electronic head	45
3.8 Combi-Clip	15	14.3 Replacing the drive unit	45
3.9 Connections at the back of the apparatus	17	15. Accessories and spare parts	47
3.10 Bath-Temp. probe	17	15.1 Quick-action jack	47
3.11 Hose connections	19	15.2 Baths	49
R-134			
4. Technical description	19	15.2.1 Water bath B-480	49
5. Supply schedule	19	15.2.2 Oil bath B-485	51
6. Starting up	19	15.3 Drive unit and sealing system	53
6.1 Controls	19	15.4 Glass assemblies	55
6.2 Functional description	21	15.4.1 Assembly A	55
6.2.1 Switching on the apparatus	21	15.4.2 Assembly V	57
6.2.2 Nominal pressure	21	15.4.3 Assembly C	59
6.2.3 Hysteresis	21	15.4.4 Assembly S	61
6.2.4 Distillation	21	15.4.5 Assembly E	63
6.2.5 Stop	23	15.4.6 Assembly CR	65
6.2.6 Hold state	23	15.5 Glass parts	66
6.2.7 Rotation	23	15.6 Electrical parts and probes	71
6.3 Working example	25	15.7 Miscellaneous	75
R-144			
7. Technical description	27	15.8 Wiring diagram	80
8. Supply schedule	27		
9. Installation	27		
9.1 Dual-Temp. probes	27		
10. Starting up	29		
10.1 Controls	29		
10.2 Functional description	31		
10.2.1 Switching on the apparatus	31		
10.2.2 Nominal pressure	31		
10.2.3 Hysteresis	31		
10.2.4 Manual distillation	31		
10.2.5 Stop	33		
10.2.6 Automatic distillation	33		
10.2.7 Hold state	35		
10.2.8 Rotation	35		
10.3 Working examples	35		
10.3.1 Manual distillation	37		
10.3.2 Automatic distillation	39		
10.3.3 Semiautomatic distillation	39		

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1. The rotary evaporator

The rotary evaporator makes it possible to carry out simple distillations rapidly and under mild conditions. In combination with a variety of accessories, it is capable of handling almost any task involving the evaporation of solvents.

- e.g.*
- *distillation of solvents*
 - *concentration of solvents*
 - *recrystallizations*
 - *Soxhlet extractions*
 - *powder and granule drying*

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2. Technical data

	R-134	R-144	B-480	B-485
<i>Voltages</i>	100, 120, 230, 240 V		120, 230 V	
<i>Tolerance</i>	-15% +10%		-15% +10%	
<i>Frequency</i>	50 and 60 Hz		50 and 60 Hz	
<i>Power consumption</i>	70 W	70 W	1300 W	
<i>Drive unit</i>	stepped motor			
<i>Rotation</i>	5 – 240 rpm			
<i>Vacuum measuring range</i>	1 – 1400 mbar			
<i>Vacuum range</i>	1 – 1400 mbar			
<i>Temperature range</i>			20 – 100°C	30 – 180°C
<i>Room-temperature</i>	10 – 35°C		10 – 35°C	
<i>Bath volume</i>			5 l	5 l
<i>Overtemperature switch</i>			140°C	260°C
<i>Power at apparatus socket</i>			max. 100 W	

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3. Installation

3.1 Notes and warnings

This unit is constructed in accordance with IEC 1010 and UL 1262 – Protective measures for electronic measuring instruments – and left the factory in a condition fully compliant with the safety regulations. To maintain this condition and ensure safe operation, the user must observe the notes and warnings contained in these operating instructions.

Before switching on, it is necessary to make sure that the operating voltage indicated on the identification plate corresponds to the mains voltage.

The mains plug must only be inserted into an earthed socket. The protective effect must not be reduced by an unearthed extension lead.

Warning

a) Any interruption of the earth wire inside or outside the unit, or loosening of the earth wire connection, can make the unit dangerous. Deliberate interruption is not permissible.

b) When opening covers or removing parts, except where this can be done by hand, live parts may be exposed. Connection points may also be live.

c) Before making adjustments, servicing, overhauling or replacing parts, the unit must be isolated from all supply points if it needs to be opened. Accordingly, if it is not possible to avoid making adjustments, servicing or making repairs with the unit open and live, this must only be carried out by a skilled person familiar with the dangers involved.

d) If it is to be assumed that the unit can no longer be operated safely, it should be taken out of service and secured against unintentional use.

It is to be assumed that the unit can no longer be operated safely

- if it has visible damage,*
- if it no longer works,*
- after prolonged storage under unfavourable conditions,*
- after severe stresses due to transport.*

Before opening the unit, always ensure first if:

- Watertap is closed**
- "On/Off" switch is in Off position**
- The unit is disconnected from the mains**

If it should be necessary to engage Büchi Service, hazardous substances (e.g. in terms of EC Directive L 360, 1976/79 or VBG 16) present on or in the unit have to be indicated and declared in the interest of working safety and environmental protection. If no declaration is made, Büchi must assume that the unit does not contain such substances.

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3.2 Standard equipment

The rotary evaporators R-134 and R-144 consist of 3 units.

- 1. Quick-action jack with drive and control unit*
- 2. Glass assembly*
- 3. Heating bath*

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

The assembly procedure is shown on the following pages

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
3.3 Installing the bath

The bath is hung in the guide rail on the Rotovapor. If a flask larger than 2 l is used, the bath simply has to be slid along this rail.

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3.4 Flanged screw couplings

All flanges (condensers, distributors, evaporating flasks, adapters) can be secured without the screw cap having to be removed. It only has to be undone until the flange can be pushed through.

In the RE variant, there is only the condenser screw coupling.

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3.5 Installing the seals

When the seals are inserted, the adjacent sequence (designated by numbers) must be followed. Conversion of the sealing system is explained in detail on page 45 under item 14.1.


The seals KD-22 and KD-26 must never be greased.

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3.6 Bath immersion angle


The angle of immersion of the evaporating flask in the bath can be changed as follows:

- Loosen the handle screw G
- Move the drive unit A into the desired position
- Retighten the handle screw G

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



3.7 Setting up the stand

- Insert the stand rod S into the orifice T
- Screw the stand rod tight with the handle screw L through the orifice R
- Secure the clamp B at the desired height with the boss K (not available with assembly A)

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3.8 Combi-Clip

In addition to hold a evaporation flask save and tight (A), the COMBI-CLIP does also:


- (B) *Open a frozen joint of the evaporation flask.*
- (C) *Remove a jammed vapour duct.*
- (D) *Hold flasks with flanges other than the original BÜCHI ones if the lever stopper is attached in the positioning holes 2.*

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
3.10 Bad-Temp. Sonde

3.9 Connections at the back of the apparatus

- 1 *Pressure sensor,
Vacuum hose to the system respectively to the
Woulff'bottle*

- CW *Cooling water, cooling water valve*

- 4 *Ventilation valve,
Possibility of nitrogen ventilation*

- T *Temperature probes,
Steam- or bath-probe*

- 3 *Mains supply,
Direct mains cable or from the bath*

- 2 *Fuses,
2 AT for 100/120 V
1 AT for 230/240 V*

- S *Auto.-probe,
Dual-Temp.-probe A
Dual-Temp.-probe B*

- P *Peripheral equipment,
Connection of further control-units*

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
Tips and notes:

Cooling water:

- For safety reasons, the hoses should be secured with commercially available hose clamps or cable binders
- To save cooling water, there is the possibility of using a recirculation condenser (e.g. B-700)
- Brittle hoses are to be replaced

Vacuum:

- The hoses used should all have the same internal diameter.
- The vacuum hoses should all be kept as short as possible
- It is recommended to insert a Woulfe's bottle between the vacuum source and the Rotavapor
- The hoses do not need to be secured
- Brittle hoses are to be replaced

R-134

4. Technical description

The rotary evaporator R-134 is the logical combination of a laboratory rotary evaporator and the tried-and-tested vacuum controller.

With a choice of 6 glass assemblies, it is capable of solving almost any problem involving the evaporation and condensation of solvents and requiring a regulated vacuum.

- e.g. – distillation of solvents
- concentration of solvents
 - powder and granule drying
 - reflux reactions
 - recrystallizations

5. Supply schedule

Quick-action jack with control unit R-134
Mains cable
Bath temperature probe
Glass assembly of your choice
Bath of your choice

6. Starting up

6.1 Controls

POWER	Apparatus on/off
DISPLAY	Switching of the display to pressure, temperature or speed of rotation
ROTATION	Adjustment of the speed of rotation with starting of the rotation
STOP	Interruption of any process
START	Distillation is started
p	Call up nominal pressure
Δp	Call up hysteresis
(AERATE)	Raising of the displayed value (with aeration of the system during a distillation)
	Lowering of the displayed value (with lowering of the pressure during a distillation)

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6.2.1 Switching on the apparatus

Key	Display	Function
POWER ON	A001 957	<ul style="list-style-type: none"> - Display of the software version - Current pressure display - Aeration valve open - Cooling water valve open for ca. 5 min, but can be closed with the STOP key

6.2.2 Nominal pressure

Key	Display	Function
p	147	<ul style="list-style-type: none"> - Display of the nominal pressure, which can be altered with the \cdot and (AERATE) keys (1 - 1400 mbar) - If no entry is made for 3 s, the displayed value is automatically stored and the display reverts to the previous state - The nominal value can be changed at any time

6.2.3 Hysteresis

Key	Display	Function
Δp	H 16	<ul style="list-style-type: none"> - Display of the hysteresis, which can be altered with the \cdot and (AERATE) keys (1 - 100 mbar und A)
	H A	<ul style="list-style-type: none"> - Automatic hysteresis, i.e. ± 60 mbar at 1000 mbar, reducing steadily to ± 2 mbar at 1 mbar pressure - If no entry is made for 3 s, the displayed value is automatically stored and the display reverts to the previous state

6.2.4 Distillation

Key	Display	Function
START	957	<ul style="list-style-type: none"> - System is evacuated to the preselected nominal value and is regulated at that value within the set hysteresis
	↓	
	147	<ul style="list-style-type: none"> - START LED lights up - Cooling water valve open - Rotation on (unless rotary knob is on 0) - Aeration valve closed - The START key can be actuated at any time except when the STOP LED is lit (no reaction)

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




6.2.5 Stop


Key	Display	Function
STOP	957	<ul style="list-style-type: none"> - System is aerated (aeration valve remains open) - STOP LED lights up - The cooling water valve remains open for ca. 5 min. - The rotation switches off - The STOP key can be used to return to the neutral state at any time - Pressing the STOP key twice switches the cooling water off immediately

6.2.6 Hold state

If foaming occurs during the distillation, it is possible to intervene in the process → **Hold state.**

It is possible to change from the distillation mode to the HOLD STATE by briefly pressing the  or  (AERATE) key.

- The current pressure with the stored hysteresis is regulated after the  or  (AERATE) key has been released
- The START LED flashes
- Pressing the START key restores the corresponding mode

6.2.7 Rotation

In the distillation mode, the rotation proceeds at the set speed as soon as START is pressed.
For rotation only, without starting a distillation, it suffices to adjust the Rotation knob slightly.

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6.3 Working example

The values in the following application example apply to the distillation of water. For the distillation of other solvents, the appropriate pressure values can be taken from the table on page 40.

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- cooling water connected
- all electrical connections correct
- all hose connections correct
- vacuum source connected
- heating medium in the bath

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POWER ON A001 – Switch on apparatus
 967 – Display of the system pressure

p	157	- Display of the stored nominal pressure
/	72	- Set desired nominal pressure
Δp	H 7	- Display of the stored hysteresis
/	H 5	- Set desired hysteresis

START 967 – System is evacuated to the desired nominal pressure and then keeps it within the desired hysteresis
 ↓
 72

Lower flask into the bath with the quick-action jack.

If foaming or violent splashing occurs, it is possible to intervene directly with the / keys.

(AERATE) Aerate system while key depressed
 Evacuate system while key depressed

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STOP 72 – Distillation is interrupted and system aerated
 ↓
 967

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
R-144

7. Technical description

The rotary evaporator R-144 is the combination of a laboratory rotary evaporator and the tried-and-tested vacuum controller with automatic distillation monitoring. With a choice of 6 glass assemblies, it is capable of handling tasks such as distillations with regulated and monitored vacuum.

e.g. – distillation of solvents

- concentration of solvents
- powder and granule drying
- reflux reactions
- recrystallizations

With the built-in controller, distillations can be carried out from beginning to end without the need for supervision. The built-in automation makes it possible to work quickly under optimal conditions with minimal solvent losses (optimized recovery rate).

Of course, the controller can also be used for simple pressure monitoring.

The R-144 is supported by an extensive range of accessories, making it a reliable laboratory aid.

8. Supply schedule

Quick-action jack with control unit R-144

Mains cable

Bath temperature probe

Dual-Temp probe

Glass assembly of your choice

Bath of your choice

9. Installation

9.1 Dual-Temp probes

- 1 Insert PTFE stop cock of the Dual-Temp probe into condenser.
- 2 Insert the PTFE hose with the sensors between the condenser wall and the cooling coils. The PTFE hose for the probe should be on top or side but not on the bottom, where the distillate runs back.
- 3 Tighten the sensor hose with the screw cap and plug the sensor cable in the R-144 or B-168 (socket S).

Note: The aeration of the ROTAVAPOR with diagonal condenser is performed with the aeration valve. If continuous feed of sample solution is also required, then, a glass assembly with the vertical condenser, type "V", must be used.

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
Vertical condensers «V» und «S»

- 1 Insert vacuum port of Dual-Temp probe into condenser.
- 2 Insert the PTFE hose with the sensors between the condenser wall and the cooling coils.
- 3 Tighten the sensor hose with the screw cap and plug the sensor cable in the R-144 or B-168 (socket S).
- 4 Vacuum connection.
- 5 Verify that the PTFE hose with the sensors is inserted straight, without bows or loops.
If necessary, rotate the vacuum port with the PTFE hose back and forward to position the hose correctly.
- 6 If the distance between the condenser wall and the coil does not allow the probe to be introduced easily, the probe can alternatively also be introduced in the middle, along the cold finger.

10 Starting up

10.1 Controls

POWER	Apparatus on/off
DISPLAY	Switching of the display to pressure, temperature or speed of rotation
ROTATION	Adjustment of the speed of rotation with starting of the rotation
STOP	Interruption of any process
MANUAL	Manual distillation is started
AUTO	Automatic distillation is started if Dual-Temp probe is connected
p	Call up nominal pressure
Δp	Call up hysteresis
(AERATE)	Raising of the displayed value (with aeration of the system during a distillation)
	Lowering of the displayed value (with lowering of the pressure during a distillation)

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10.2 Functional description

The following numerical examples represent an arbitrary choice.

10.2.1 Switching on the apparatus

Key	Display	Function
POWER ON	A001	- Display of the software version
	957	- Current pressure display - Aeration valve open - Cooling water valve open for ca. 5 min, but can be closed with the STOP key

10.2.2 Nominal pressure

Key	Display	Function
p	147	- Display of the nominal pressure, which can be altered with the \uparrow and \downarrow (AERATE) keys (1 – 1400 mbar) - If no entry is made for 3 s, the displayed value is automatically stored and the display reverts to the previous state - The nominal value can be changed at any time

10.2.3 Hysteresis

Key	Display	Function
Δp	H 16	- Display of the hysteresis, which can be altered with the \uparrow and \downarrow (AERATE) keys (1 – 100 mbar und A).
	H A	- Automatic hysteresis, i.e. ± 60 mbar at 100 mbar, reducing steadily to ± 2 mbar at 1 mbar pressure - If no entry is made for 3 s, the displayed value is automatically stored and the display reverts to the previous state

10.2.4 Manual distillation

Key	Display	Function
START	957	- System is evacuated to the preselected nominal value and is regulated at that value with in the set hysteresis
	147	- MANUAL LED lights up - Cooling water valve open - Rotation on (unless rotary knob is on 0) - Aeration valve closed - The MANUAL key can be actuated at any time except when the STOP LED is lit (no reaction)

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10.2.5 Stop

Key	Display	Function
STOP	957	<ul style="list-style-type: none">- System is aerated (aeration valve remains open)- STOP LED lights up- The cooling water valve remains open for ca. 5 min- The flask is lifted out of the bath- The rotation switches off- The STOP key can be used to return to the neutral state at any time- Pressing the STOP key twice switches the cooling water off immediately

10.2.6 Automatic distillation

Conditions:

- Dual-Temp probe connected
- Rotation not on 0 (-> E-8)

In the AUTO mode, the manually selected nominal values are not active. However, they can be further adjusted at any time without affecting the AUTO mode.

Key	Display	Function
AUTO	957	<ul style="list-style-type: none">- The apparatus searches for the optimal distillation conditions and regulates the established pressure under these conditions within the automatic hysteresis- AUTO LED lights up- Cooling water valve open- Rotation on- Aeration valve closed- The AUTO key can be actuated at any time except when the STOP LED is lit (no reaction)

Description of the phases of automatic distillation:

- Phase 1:
 - Stepwise evacuation until distillation begins
- Phase 2:
 - Optimization of the distillation pressure
- Phase 3:
 - Monitoring of the distillation
- Phase 4:
 - End of distillation (rotation off, flask lifted out of the bath, aeration valve remains open, cooling water valve closes after ca. 5 min, beep at minute intervals)
 - Cooling water and beep can be switched off by pressing the STOP key

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


10.2.7 Hold state

If foaming occurs during manual or automatic distillation, it is possible to intervene in the process → **Hold state**.

It is possible to change from the manual and automatic mode to the HOLD STATE by briefly pressing the  or  (AERATE) key.

- The current pressure with the stored hysteresis is regulated after the  or  (AERATE) key has been released
- The MANUAL or AUTO LED flashes
- Pressing the AUTO or MANUAL key restores the corresponding mode

10.2.8 Rotation

In the AUTO and MANUAL mode, the rotation proceeds at the set speed as soon as AUTO or MANUAL is pressed.

For rotation only, without starting a corresponding distillation run, it suffices to adjust the Rotation knob slightly.

10.3 Working examples

The values in the following application examples are approximate values for the distillation of water. For the distillation of other solvents, the appropriate pressure values can be taken from the table on page 40.

The known 20/40/60°C rule (see page 39) also applies here.

The following points should be checked before every distillation:

- cooling water connected
- all electrical connections correct
- all hose connections correct
- vacuum source connected
- heating medium in bath
- Dual-Temp probe connected (only for automatic and semiautomatic distillation)

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10.3.2 Automatic distillation

Set bath to the desired temperature (60°C).

POWER ON A001 – Switch on apparatus
967 – Display of the system pressure

Lower flask into the bath with quick-action jack.

AUTO 967 – System automatically searches for the optimal distillation pressure, regulates it within the automatic hysteresis and finally ends the distillation

If foaming or violent splashing occurs, it is possible to intervene directly with the \uparrow / \downarrow keys.

(AERATE) Aerate system while key depressed
Evacuate system while key depressed

The distillation can be continued by pressing the START-key.

STOP 72 – Distillation is interrupted and
↓
967 system aerated


10.3.3 Semiautomatic distillation

The distillation is started manually as described under 10.3.1. When the distillation is running, the AUTO key is actuated. The distillation is then completed independently.

11. Tables

11.1 Choosing the distillation temperatures

As a rule the cooling water temperature cannot be freely chosen. The distillation conditions therefore depend inter alia on the temperature of the available cooling water. This results in the following rule of thumb:

Advantages of 60°C bath temperature:

- The evaporating flask can be changed without scalding the hands
- The rate of evaporation of water from the heating bath is not very great
- The distillation capacity of the Rotavapor is efficiently utilized

If distillation is carried out according to the rule of thumb 20 – 40 – 60°C, the pressure is to be chosen so that the vapour temperature (and hence the boiling point) is 40°C. The required vacuum can be taken from the solvent table on the next page or determined fully automatically by the apparatus with the Dual-Temp probe.

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12. Troubleshooting

12.1 Error messages

<i>Message</i>	<i>Cause</i>	<i>Remedy</i>
<i>E—1</i>	<i>Dual-Temp probe faulty</i>	<i>Replace</i>
<i>E—2</i>	<i>Dual-Temp probe faulty or not connected</i>	<i>Have replaced at service point</i>
<i>E—3</i>	<i>Pressure sensor faulty</i>	<i>Notify service point</i>
<i>E—4</i>		<i>Notify service point</i>
<i>E—5</i>	<i>Leak rate too high</i>	<i>Remedy leak</i>
<i>E—6</i>		<i>Notify service point</i>
<i>E—7</i>		<i>Notify service point</i>
<i>E—8</i>	<i>No rotation during automatic operation</i>	<i>Rotation on</i>
<i>E—9</i>	<i>Optional equipment incorrectly connected at the back of the apparatus</i>	<i>Check connections</i>
<i>E—10</i>	<i>Hardware fault</i>	<i>Notify service point</i>

12.2 Other malfunctions

<i>Fault</i>	<i>Cause</i>	<i>Page</i>
<i>– Apparatus does not work</i>	<i>– Apparatus not connected to the mains supply</i>	<i>17</i>
	<i>– Fuse faulty</i>	<i>17</i>
<i>– Quick-action jack does not lift the flask out of the bath on its own</i>	<i>– Motor faulty</i>	
<i>– Bath does not heat up</i>	<i>– Bath not connected to the mains supply</i>	<i>48-51</i>
	<i>– Overtemperature cut-out has responded</i>	<i>48-51</i>
	<i>– Fuse faulty</i>	<i>48-51</i>
<i>– Flask does not rotate</i>	<i>– Rotation on 0</i>	<i>19/23</i>
<i>– C in the display</i>	<i>– Calibration mode</i>	<i>29/35</i>
	<i>– Press STOP to quit this mode</i>	

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

13. Cleaning

It is possible to clean the unit with a moist cloth. For the cleaning of the bath it is also possible to use a soft detergent.

Do not clean the unit with solvents.

14. Assembling and dismantling

14.1 Converting the sealing system

Conversion from the RE to the EL sealing system:

- Push out covering ring (6) with plastic spatula
- Screw on EL adapter (3):
 - Push in stop button (8)
 - Screw EL adapter into orifice (2)
 - Tighten EL adapter (3) with bolt (7)
- Screw on flanged screw coupling (4) with spring (5)

Conversion from the EL to the RE sealing system:

- Unscrew flanged screw coupling (4) with spring (5):
 - Push in stop button (8)
 - Unscrew EL adapter with bolt (7)
 - Insert covering ring (6) into orifice (2)

14.2 Replacing the electronic head

Removing head (1):

- Switch off apparatus and pull out mains plug
- Unscrew head fixing screw (6)
- Carefully lift head (1) at the back and pull it out forwards
- Pull out the two plugs (4) + (5)

Fitting head (1):

- Push in the two plugs (4) + (5)
- Push in head (1) forwards and then lower at the back
- Screw up head (1) with head fixing screw (6)
- Reconnect apparatus to the mains

14.3 Replacing the drive unit

Removing drive unit (2):

- Switch off apparatus and pull out mains plug
- Loosen handle screw (3)
- Carefully lift out drive unit (2)
- Pull out drive unit plug (4) on the head (1)

Fitting drive unit (2):

- Push in drive unit plug (4) on the head (1)
- Screw drive unit (2) tight in the desired position with handle screw (3)

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32167 Bellow frame, bottom

Operation:

- The flask can be moved into and out of the bath by applying the same force in any position of the quick-action jack
- To minimize the force applied when operating the quick-action jack, the energy accumulator can be matched to the weight of the glass assembly
- Adjust the screw in the orifice (S) with a screwdriver (turn to the right to increase force and to the left to decrease force)
- The quick-action jack has an additional support at the back, at the bottom of the tower. This support can be pulled out by loosening the Allen screw (I)
- Any rocking of the quick-action jack on an uneven surface can be compensated by adjusting the levelling screw (N)

Warning!

Because the spring is under tension, the energy accumulator may only be opened by authorized service technicians.

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


Controls

POWER Main switch for switching on the bath

HEATING Pilot lamp indicates when the bath is heating

TEMPERATURE °C Rotary knob for setting the desired bath temperature (20 - 100°C)

A built-in overtemperature cut-out prevents the heater from burning out. This safety device resets itself automatically when the heater has cooled down.

If deionized or distilled water is used in the water bath, 0,2% of borax ($\text{Na}_2\text{B}_4\text{O}_7 \cdot x 10\text{H}_2\text{O}$) must be added.

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00591 Supporting ring for distributor
 00633 Retaining spring
 00636 Seal KD-22
 10179 Seal KD-26
 11623 Back-up ring, stainless steel
 11629 Viton O-ring
 11635 Protecting ring with 2 O-rings
 23375 Combi clip, STJ 29.2
 25424 Pin
 27267 Combi clip, STJ 24
 30635 Motor
 30809 Hall sensor
 32004 Rubber buffer for motor securing device
 32005 Sliding ring
 32006 Covering ring
 32007 EL adapter
 32014 Flanged screw coupling
 32017 Insert spring
 32018 Rubber buffer for motor mounting
 32026 Drive belt
 32034 Half-dish
 32035 Cover dish
 32193 EPDM O-ring
32370 Bearing cpl. (00568+00569
 +00633+32005+32009+32153
 +32155+32157)
32369 Retention pin cpl.

32299 Conversion kit RE → EL cpl.

00591 Supporting ring for distributor
 10179 Seal KD-26
 11623 Back-up ring, stainless steel
 11162 Seal extractor
 11629 Viton O-ring
 11635 Protecting ring with 2 O-rings
 32007 EL adapter
 32014 Flanged screw coupling
 32017 Insert spring
 27773 Bolt

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
15.4.1 Glass assembly A

32300 Glass assembly A cpl.

(as shown in drawing with 1l evaporating flask and 1l receiving flask)

- 00636 Seal KD-22
- 00637 Glass stopcock STJ 18,8/38
- 00643 PTFE tube, 460 mm
- 03275 SGJ clip
- 23375 Combi clip STJ 29,2
- 23875 Screw cap GL-10
- 23900 O-ring 2.7 x 3 mm
- 28096 PTFE tube, 600 mm
- 32001 Vapour duct, STJ 29,2/32, 280 mm
- 32014 Flanged screw coupling
- 32022 Diagonal condenser A (inkl. 3 x 18913+23875+23900)
- 32302 Diagonal condenser A cpl.** (0063+00637+00643+28096+32022)
- 32340 Vapour duct with Combi clip, STJ 29,2/32, 280 mm**
- 32885 Hose nipple GL-14, bent, cpl.

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15.4.4 Glass assembly S

32319 Glass assembly S cpl.

(as shown in drawing with 1l evaporating flask and 1l receiving flask without vacuum adapter 1006)

- 00591 Supporting ring for distributor
- 00637 Glass stopcock STJ 18,8/38
- 00646 PTFE tube, 300 mm
- 00659 Vertical condenser S
(incl. 2 x 18913+25435)
- 00661 Vacuum connector
- 00935 PTFE disc (draining disc)
- 00936 Shut-off valve
- 01006 Vacuum adapter with cock
- 03275 SGJ clip
- 10179 Seal KD-26
- 11162 Seal extractor KD-26
- 11623 Back-up ring, stainless steel
- 11635 Protecting ring with 2 O-rings
- 23375 Combi clip STJ 29,2
- 23733 Adapter STJ 29,2/32
- 23875 Screw cap GL-10
- 23900 O-ring 2.7 x 3 mm
- 25419 Adapter with Combi clip, STJ 29,2/32**
- 25435 Screw cap GL-14
- 28096 PTFE tube, 600 mm
- 32014 Flanged screw coupling
- 32021 Distributor (incl. 23875+23900)
- 32331 Distributor cpl.**
(00591+00637+00646+00935+00936+28096+32021)
- 32885 Hose nipple GL-14, bent, cpl.

Plastic + Glas (PLG)

32320 Glass assembly S cpl. PLG

(as shown in drawing with 1l evaporating flask, uncoated, and 1l receiving flask, coated, without vacuum adapter 1006)

- 23531 Vertical condenser S PLG
(incl. 2 x 18913+25435)
- 32310 Distributor PLG
(incl. 23875+23900)
- 32332 Distributor cpl. PLG**
(00591+00637+00646+00935+00936+28096+32310)

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
15.4.5 Glass assembly E

32323 Glass assembly E cpl.

(according to drawing with 1l evaporating flask and 1l receiving flask)

- 00591 Supporting ring for distributor
- 00637 Glass stopcock STJ 18,8/38
- 00646 PTFE tube, 300 mm
- 00935 PTFE disc (draining disc)
- 01002 Expansion vessel
- 01003 U-tube (incl. 25435)
- 01004 Vacuum adapter (incl. 18913)
- 03275 SGJ clip
- 10179 Seal KD-26
- 10183 Condenser (incl. 2x18913)
- 11151 Condenser clamp
- 11162 Seal extractor KD-26
- 11623 Back-up ring, stainless steel
- 11635 Protecting ring with 2 O-rings
- 23375 Combi clip STJ 29,2
- 23733 Adapter STJ 29,2/32
- 23875 Screw cap GL-10
- 23900 O-ring 2.7x3 mm
- 25419 Adapter with Combi clip, STJ 29,2/32**
- 25435 Screw cap GL-14
- 27343 Wing screw
- 27344 Cross piece cpl.
- 27347 Cross piece
- 28096 PTFE tube, 600 mm
- 32008 Rod, bent
- 32014 Flanged screw coupling
- 32019 Handle screw
- 32049 Distributor (incl. 23875+23900)
- 32885 Hose nipple GL-14, bent, cpl.

Plastic + Glas (PLG)

32324 Glass assembly E cpl. PLG

(as shown in drawing with 1l evaporating flask, uncoated, and 1l receiving flask, coated)

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

Evaporating flasks (pear shaped)
For EL systems only with adapter
flange STJ

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

Flanged evaporating flask
for EL systems

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
Flanged drying flask
for EL systems

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
**Special flasks for drying
powdery substances**

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
**"Soxhlet" extraction unit
for EL models**

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
<i>Solid/liquid, complete</i>	<i>Code</i>
<i>Capacity: 200 ml</i>	<i>11744</i>
<i>Capacity: 500 ml</i>	<i>11745</i>
<i>Extraction cartridges</i>	
<i>Capacity: 200 ml</i>	<i>08560</i>
<i>Capacity: 500 ml</i>	<i>00989</i>

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


15.6 Electrical parts and probes

23349 Water jet pump B-164 24 V

11289	Reduction piece 1/2" – 3/8"
11290	Reduction piece 1/2" – 3/4"
11514	Sieve \varnothing 18.5 for water tap
11515	Gasket 18.5 x 13 x 1 mm
16475	Water jet pump
18913	Hose adapter GL-14 cpl.
20859	Water valve B-162, 24 V DC
20882	Vacuum valve cable B-162
23347	Gasket "Kalrez" 10.5 x 1.78
20814	Connection GL-14
20815	Back hold valve

23349 Water jet pump B-167 24 V with cooling water valve

01308	Needle of needle valve
03560	O-ring 6 x 1.5 mm
11289	Reduction piece 1/2" – 3/8"
11290	Reduction piece 1/2" – 3/4"
11514	Sieve \varnothing 18.5 for water tap
11515	Gasket 18.5 x 13 x 1 mm
16475	Water jet pump
18913	Hose adapter GL-14 kpl.
23347	Gasket "Kalrez" 10.5 x 1.78
27104	Core with spring large
27107	Core with spring small
27108	Coil
30073	Vacuum valve cable B-167
20814	Connection GL-14
20815	Back hold valve

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
15.7 Miscellaneous**22699 Cooling water valve 24 V**

- 01308 Needle valve for reg. nozzle
- 03560 O-ring 6 x 1.5 mm
- 11289 Reduction piece 1/2" – 3/8"
- 11290 Reduction piece 1/2" – 3/4"
- 11514 Sieve ø 18.5 for water tap
- 11515 Gasket 18.5 x 13 x 1 mm

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
25519 Woulff bottle 3-neck cpl.

- 18913 Hose adapter, GL-14, cpl.
- 18914 Screw cap, GL-14
- 18916 Hosester, GL-14, Ø 8 mm
- 25518 Woulff bottle (glass)


25519 Flacon de Woulff cpl.

- 18913 Olive, GL-14, cpl.
- 18914 Ecrou chapeau, GL-14
- 18916 Olive coudee, GL-14, Ø 8 mm
- 25518 Flacon de Woulff (verre)


11152 Bride Ø 100 mm*for scaffolding mounting of the Woulff bottle***11152 Bride Ø 100 mm***pour le montage stativ du flacon de Woulff*


- 11606 Wasserregulierdüse 1/2"**
 (01308+03560+11514+11515)
- 01308 Regulierdüse
 - 03560 O-Ring, 6,0 x 1,5
 - 11514 Einsatzsieb
 - 11515 Dichtung zu Regulierdüse


- 18797 Vakuummeter-Satz**
 (00394+10056+18826)
- 00394 Y-Stück O, 10 mm
 - 10056 Stativstange, 12 x 300
 - 18826 Vakuummeter

- 17622 Vakuumschlauch 16/6**
- 04133 Kühlschlauch 9/6**

- 11289 Reduziergewinde 1/3" - 3/8"**
- 11290 Reduziergewinde 1/2" - 3/4"**

- 01330 Glisseal Laborfett 60 g**

11606 Water reg. valve 1/2"
(01308+03560+11514+11515)
01308 Needle valve
03560 O-ring, 6,0 x 1,5
11514 Sieve for needle valve
11515 Gasket for needle valve

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



Only an authorized technician from your BÜCHI agency may open the cover plate of the electronic module and carry out servicing work.

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