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Primus 25 advanced[®] Innovative PCR Technology / Innovative PCR-Technologie

CONTENTS

1	SAFETY INSTRUCTIONS	2
2	LICENSE NOTICE	2
3	SYSTEM OVERVIEW	3
4	INSTALLATION	4
5	MAIN MENU	4
5.1	RUN	4
5.2	SELECT + RUN PROGRAM	4
5.3	RUN INCUBATION	5
5.4	MODE/EDIT/UTILS	5
5.4.1	PROGRAM	5
5.4.2	VIEW	5
5.4.3	PRINT	5
5.4.4	SETUP	5
6	PROGRAMMING	6
6.1	Program steps	6
6.2	Sample program	7
7	TECHNICAL SUPPORT AND ORDERING INFORMATION	8
8	TECHNICAL SPECIFICATIONS	9
9	MAINTENANCE AND REPAIR OF THE THERMAL CYCLER	10
9.1	Cleaning, maintenance and repair of the thermal cycler	10
9.2	Temperature check-up	10
9.3	Calibration	10
10	WARRANTY	10
11	DISPOSAL INSTRUCTION	11

1 SAFETY INSTRUCTIONS

Before the first use of the thermal cycler please read the instruction manual entirely.

Special notice should be taken of the following

Caution of dangerous voltage

Please ensure the voltage indicated on the device exactly matches your local electrical supply

Caution of dangerous explosive material

Explosive or reactive material mustn't be heated in the thermal cycler.

Caution - liquids

Ensure that no liquids can enter the device.

Caution - hot surfaces

The thermal plate and the heated lid quickly reach temperatures above 50°C. There is danger of burning! Keep the lid closed until the temperatures reaches 30°C or less.

Only use materials (tubes, tube caps, plates and sealing films) which are heat resistant to 120°C


Caution – environment

The ventilation of the device must not be covered.

2 LICENSE NOTICE

This instrument is licensed for research and development and for uses other than human in vitro diagnostics under one or more of the following patents of Applied Biosystems Corporation: U.S. Patent Nos. 5,656,493, 5,038,852, 5,333,675, 5,475,610 (claims 1-159 and 164-166), 6,703,236 (claims 7-10) and 7,238,517 or corresponding claims in their non-U.S. counterparts. No right is conveyed expressly, by implication or by estoppel under any other patent of Applied Biosystems, including but not limited to U.S. Patent No. 6,814,934 and its non-U.S. counterparts, which describe and claim thermal cyclers capable or real-time detection.

3 SYSTEM OVERVIEW

- | | | |
|----|-------------------------|---|
| 1 | Heated lid | |
| 2 | Universal block | With automatic height adjustment |
| 3 | Display | for 25 x 0.2 ml tubes or 13 x 0.5 ml tubes with flat caps |
| 4 | Cursor key pad | Displays menus and programs |
| 5 | Alphanumeric key pad | Moves the cursor on the display |
| 6 | Alpha key | For entering parameters |
| 7 | Shift lock key | Switches between numbers and letters |
| 8 | Delete key | Switches between lower and upper case |
| 9 | Insert key | Deletes parameters or functions |
| 10 | Mode key | Inserts parameters or functions |
| 11 | Stop / Esc key (red) | Selects edit and display options |
| 12 | Run / Enter key (green) | Exits a menu option to return to a higher menu level |
| | | Selects menus and program options |

4 INSTALLATION

Before using the Primus 25 advanced® for the first time, check the voltage range and be certain it corresponds to your power supply.

Place your Primus 25 advanced® in an adequate position in the lab. The unit should always be surrounded by at least 25 cm of free space to ensure proper cooling. Place several units beside each other, not behind each other.

Please avoid running the Primus cyler without having tubes placed inside the block.

5 MAIN MENU

The software for the Primus 25 advanced® is designed for simple and intuitive use. Generally, the arrow keys are used to move the cursor on the display screen. Text is entered using the alphanumeric key pad. The position of a letter on the key tells how many times the key must be pressed (in rapid sequence) to obtain that letter. Use the Run/Enter key to access the various function menus, and the Stop/Esc key to exit the function menus and return to the main menu.

When switching on the instrument, the screen first displays an initialization routine, then the main menu appears from which all other functions of the Primus 25 advanced® are accessible.

5.1 RUN

To start the last run program select RUN. The program starts immediately.

Attention: The modifications in the program after the last run will be ignored. In this case, select the program with SELECT + RUN PROGRAM and press enter to start.

Programs can be interrupted at any time during a run with the Stop/Esc key. Follow the instructions on the display.

5.2 SELECT + RUN PROGRAM

To run a previously stored program select SELECT + RUN PROGRAM. A list of all stored programs (up to 90) is displayed. Select the desired program with the arrow keys, or enter the program number directly. Press enter to start. Before each run, a user name may be entered via the alphanumeric key pad. At this point you can decide whether a GLP report shall be printed via the parallel printer interface (GLP ONLINE PRINTS; GLP = Good Laboratory Practice). Switch between YES und NO with the (↑↓) cursor keys. A predicted program run time will be displayed before the start.

Programs can be interrupted at any time during a run with the Stop/Esc key. Follow the instructions on the display.

5.3 RUN INCUBATION

The Primus 25 advanced® may be used as incubator with the INCUBATION function. Choose the desired temperature (INCUB) and switch the lid heating on or off (LIDHT). The current temperature, the status of the lid heating and the total time are displayed. Use the Stop/Esc key to exit the INCUBATION function.

5.4 MODE/EDIT/UTILS

With this menu option you can edit and review all programming and operating functions.

5.4.1 PROGRAM

With the submenu PROGRAM programs may be created and changed..

5.4.2 VIEW

With the VIEW option, an overview of the stored programs may be viewed without risking an accidental alteration of the data. The GLP report of the last program run can be viewed on the display. Information about setup, system, and block are available here.

5.4.3 PRINT

With the PRINT menu programs, program lists, GLP reports and equipment information may be printed via the parallel printer interface.

5.4.4 SETUP

The basic setup of the Primus 25 advanced® is entered via SETUP. To enter the time and date use CLOCK. If several Primus thermal cyclers are connected via the lab bus (RS485), enter the machine's address here.

6 PROGRAMMING

Thanks to the convenient user interface of the Primus 25 advanced®, it is easy to create PCR programs of various degrees of complexity.

In the main menu, move the cursor to MODE/EDIT/ UTILS and press Enter twice, or go to the PROGRAM menu with the Mode key.

The following options are available:

NEW	To create a completely new program.
EDIT	To change an existing program.
COPY	To copy a previously stored program (e.g. to modify a program).
DELETE	To delete one or all programs.

6.1 Program steps

The following program steps are available for programming the Primus 25 advanced®:

TEMP	<ul style="list-style-type: none"> a) Temperature b) Temperature change within a program loop. The sign states if the temperature is increased or decreased. c) Period of time for which the preset temperature is to be held. Using the „minus“ key "FOREVER" can be chosen d) Time change within a program loop. The sign indicates whether the time is increased or decreased.
RAMP	The ramp always starts at the temperature defined last and runs to the set temperature.
PAUSE	The process is paused for the stated time at the preset temperature. The pause may be interrupted at any time by pressing a key. The thermal cycler emits an audible signal when the pause is reached. After pressing the key Enter, the program continues with the next step.
LOOP[Start of a program loop which is closed by a subsequent Loop] step. Input the number of cycles. No "loop within a loop" (nested loops) are possible!
LOOP]	End of a program loop.
LIDHT	The heated lid may be switched on or off. The stated temperatures are kept constant within the program after the Lidheat step. Temperatures from 70 °C to 120 °C are valid.
STORE	Cooling to 8°C, the time is programmable.
LIDOP	The program is stopped and the user is asked to open the lid.
LIDCL	The program is stopped and the user is asked to close the lid.
END	The END step cannot be programmed or deleted, but is always found at the end of a program. When reaching this step, the instrument is actively run to ambient temperature and the operator is called to the instrument .

6.2 Sample program

Select **MODE**, then **PROGRAM** and create a new program with **NEW**. Enter a program name and confirm with **Enter** to access the programming functions.

01 LIDHT ON 120 °C Heated lid on, temperature 120 °C

Press **Enter** to insert a new program step. Move from each parameter to the next with the cursor keys (←→) or with **Enter**.

02 TEMP 94.0 °C +0.0 °C Denaturation at 94 °C
0:02:00 +0:00 For 2 minutes

Enter another program step with **Enter**. Choose one of the steps listed above by using the (↑↓) cursor keys, e.g., **LOOP**].

03 LOOP[30x Open a program loop - 30 cycles

Continue by entering a temperature profile.

04 TEMP 94.0 °C +0.0 °C Denaturation at 94 °C
0:00:30 +0:00 For 30 seconds
(Press 0 to enter the time in seconds)

05 TEMP 50.0 °C -0.5 °C Annealing at 50 °C with stepwise
0:00:30 +0:00 temperature decrement by 0.5 °C (Touch-
down)

06 TEMP 72.0 °C +0.0 °C Elongation with stepwise time increment
0:00:30 +0:05 by 5 seconds

End the cycling step with **LOOP]** and with the next step switch the heated lid off.

07 LOOP] End the program loop
08 LIDHT OFF Lid heat off

Ramps may be inserted into any program. The ramp starts from the last programmed temperature (here 72 °C) and decreases temperature to the set temperature (here 8 °C) at a constant rate.

09 RAMP 4.0 °C Target temperature 8 °C with
1.0 °C/s a speed of 1 °C per second

10 STORE 8.0 °C +0.0 °C Final temperature 8 °C
Forever +0.00 (Minus key!)

11 End End of program

With the cursor keys (↑↓) the **HEADER of the program is accessible**. Here, a comment (**CMNT**), author (**AUTH**), and user (**USER**) may be input. Also, a denaturation step in case of power failure (**PWRFAIL-DENAT**) may be programmed.

Exit the program editor with **Stop/Esc**. Edited or newly created PCR programs are tested in a simulation run before being saved and before running. You can only exit the editing mode if the program makes sense. If not, the editor will indicate possible error sources. A calculated program run time can be displayed (**CALCTIM**), and you can save your newly created PCR program with **Enter**.

7 TECHNICAL SUPPORT AND ORDERING INFORMATION

For technical questions and more detailed information on PEQLAB's products please visit www.peqlab.com to find the respective contact person.

This manual is supplied with the following products:

Product	Cat. No.
Primus 25 advanced® Personal cycler with universal block for 25 x 0.2 ml tubes of 13 x 0.5 ml tubes with flat caps	95-4002
110 V version for USA	95-4002-US

As reaction vessels you can use 0.2 ml tubes or 0,5 ml tubes with flat caps. For optimum results, only use thin-walled tubes with a slim cone.

Item	Description	Cat. No.
PCR Tubes 0.2 ml (domed cap)	natural, 1 bag of 1000	82-0337-A
	natural, 5 bags of 1000	82-0337-B
	assorted colours, 5 bags of 200	82-0491-A
	assorted colours, 25 bags of 200	82-0491-B
PCR Tubes 0.2 ml (flat cap)	natural, 1 bag of 1000	82-0620-A
	natural, 5 bags of 1000	82-0620-B
	assorted colours, 5 bags of 200	82-0622-A
	assorted colours, 25 bags of 200	82-0622-B
PCR Tubes 0.5 ml (flat cap)	natural, 1 bag of 1000	82-0350-A
	natural, 5 bags of 1000	82-0350-B
	assorted colours, 5 bags of 200	82-0533-A
	assorted colours, 25 bags of 200	82-0533-B

8 TECHNICAL SPECIFICATIONS

Temperature

Temperature range of block	4.0°C to +105.0°C
Control accuracy of block	± 0.1 °C
Block Uniformity (at 72°C)	± 0.7°C
Cooling rate of block	2°C/sec
Heating rate of block	2°C/sec
Temperature range of heated lid	70°C to 120°C
Minimal setting steps heated lid	1 °C

Block

Capacity	25 x 0.2 ml tubes or 13 x 0.5 ml tubes with flat caps
Heated lid	Automatic height adjustment for various tube heights

Software

Memory	90 Programs with a maximum of 99 Steps/Program
Programming steps	Temp, Ramp, Pause, Loop[, Loop], LidHt, LidOp, LidCl,
Store	
Time (h:min:sec)	0:00:01 to 9:59:59, infinite ('forever')
Time increment/decrement (min:sec)	0:01 to 9:59
Temperature increment / decrement	0.1 °C to 9.9 °C
GLP report	
Instant incubation	
Autocalibration	
Comment, author, user for each program	
Remote control for use with PC	
Power failure recognition and auto restart	
Internal real time clock	
Service control system	

Interfaces

RS232 Interface (9600 Baud)
Parallel printer interface

General

Power consumption	120 VA
Power supply	230 V/115 V AC, 50/60 Hz
Dimensions (W x H x D)	225 x 250 x 280 mm
Weight including block	approx. 6.3 kg

9 MAINTENANCE AND REPAIR OF THE THERMAL CYCLER

9.1 Cleaning, maintenance and repair of the thermal cycler

For safety reasons the device must be switched off and the mains plug pulled before cleaning is carried out!

The surface of the thermal cycler can be wiped off with a damp cloth or ethanol. Do not use aggressive or scouring cleaners or organic solvents for cleaning. The device should be protected from aggressive chemicals. Make sure that no liquid penetrates into the interior of the device.

If sample liquids penetrates into the drill holes of the thermal cycler clean them immediately with mild soap followed by distilled water. You can use mild disinfectants for decontamination.

Maintenance (temperature check-up see 9.2) should be carried out regularly to detect possible technical failures.

The electric fuses are placed at the back of the device between the power switch and the mains plug.

After pulling off the mains plug and pushing upwards a safety catch the fuses can be changed. Use only fuses with correct values (indications about the fuse type are placed at the back of the device).

The device may only be opened by a qualified specialist. Unauthorized work on the device voids the warranty. All kinds of repairs may only be carried out by authorized persons using original replacement parts. The replacement of single components of the device should not be done by the user but exclusively by authorized specialists. Therefore the device must be sent to the manufacturer. The surface of the device should be decontaminated with a lint free cloth soaked with 70 % ethanol before sending. The drill holes of the thermoblock shall be cleaned with a 1.5 % hypochloride solution with the help of a cotton bud. Please fill in and sign the decontamination certificate and send it to the manufacturer together with the device.

9.2 Temperature check-up

Maintenance of the temperature at regular intervals (e.g. biannual) is recommended. The thermometer for use should be equipped with wired sensing devices (e.g. Pt1000) and should have an accuracy of measurements of 1/10 class B ($\Delta T = \pm 0.03^\circ\text{C}$) in the range of 0 to 100°C .

While processing the temperature control the room temperature should be 25°C at maximum. Three different measurements at the right and left border and in the middle of the thermal plate should be carried out in tubes filled with glycerine.


9.3 Calibration

If there is a reasonable suspicion that the device has deviated from the calibration range, because the variations of the temperature check-up are considerably out of the acceptable tolerance, a new calibration must be carried out by the manufacturer. Please contact your local distributor.

10 WARRANTY

The period of warranty for the product at hand is 24 month from date of purchase. The warranty does not apply for defects caused by incorrect, improper or inappropriate use. The customer has to check up on obvious defects promptly upon receipt of the device and send a written report. Otherwise the claim for warranty expires. The manufacturer reserves the right to carry out changes to the thermal cycler or to the technical documentation of the thermal cycler on the basis of product development or improvement without prior notice.

11 DISPOSAL INSTRUCTION

This product is subject to the WEEE disposal regulations and cannot be disposed of with regular waste.

Please contact PEQLAB Biotechnologie GmbH or your local distributor for disposal of Primus 25 advanced® thermal cyclers. Waste electronic equipment is not useless waste. Environmentally friendly disposal can help to retrieve valuable resources. Thereby you contribute to the protection of the environment and human health.

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
This instrument is licensed for research and development and for uses other than human in vitro diagnostics under one or more of the following patents of Applied Biosystems Corporation: U.S. Patent Nos. 5,656,493, 5,038,852, 5,333,675, 5,475,610 (claims 1-159 and 164-166), 6,703,236 (claims 7-10) and 7,238,517 or corresponding claims in their non-U.S. counterparts. No right is conveyed expressly, by implication or by estoppel under any other patent of Applied Biosystems Corporation, including but not limited to U.S. Patent No. 6,814,934 and its non-U.S. counterparts, which describe and claim thermal cyclers capable of real-time detection.

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D PEQLAB Biotechnologie GmbH, 91052 Erlangen, Freecall (D): 0800 100 20 16, info@peqlab.de, www.peqlab.de
AT PEQLAB Biotechnologie GmbH, 6404 Polling, Tel: +43 (0) 5238 84 169, info@peqlab.at, www.peqlab.at
UK PEQLAB Ltd., Southampton SO31 7ZN, Freephone (UK): 0808 202 1302, info@peqlab.co.uk, www.peqlab.co.uk
USA PEQLAB LLC, Wilmington, DE 19810, Toll-Free (US): 877 737 5220, info@peqlab.us, www.peqlab.us