50012338

D

F

GB

E

**Heraeus**  
INSTRUMENTS

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**Baureihe 6000**

**Brutschränke** B 6030, B 6060, B 6120, B 6200, B 6420, B 6760

Operating Instructions

**Series 6000**

**Incubators** B 6030, B 6060, B 6120, B 6200, B 6420, B 6760

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
**Série 6000**

**Incubateurs** B 6030, B 6060, B 6120, B 6200, B 6420, B 6760

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**Serie 6000**

**Incubadores Microbiológicos** B 6030, B 6060, B 6120, B 6200, B 6420, B 6760

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




Stand: 2. Februar 1993

## Table of Contents

<b>1. Introduction</b>	<b>3</b>
<b>2. Installation and Connections</b>	<b>4 - 7</b>
<b>3. Operating Elements</b>	<b>8 - 20</b>
Air-flap adjustment	
Temperature registration	
Temperature control	
Thermal protection	
Timers	
Master switch	
<b>4. Commissioning</b>	<b>21 - 22</b>
<b>5. Maintenance and Repairs</b>	<b>23 - 24</b>
<b>6. Appendix</b>	
Table 3 Classification of Heating Ovens acc. to DIN-Standard 12 880, Part 1	<b>25</b>
Technical data	<b>26 - 28</b>
Dimensioned sketch	<b>29</b>
Controller adjustment	<b>30</b>
Eurotherm controller with interface RS 232	<b>31</b>
Eurotherm controller with interface RS 485	<b>32</b>
Pt 100	<b>33</b>

## 1. Introduction

Your new incubator is characterized by a high degree of economic efficiency and outstanding quality.

Please read the following instructions carefully so that you'll be in a position to operate your incubator without any difficulty.

### Special Notes

**The incubators meet the requirements of DIN-Standard 12 880 and VDE-Regulation 0700, Part 1.**

**Radio-interference suppression: acc. to VDE-Regulation 0875, Part 1.**

**Type of protection: IP 20 acc. to DIN-Standard 40050.**

**The thermal safety class as per section "Thermal protection" (Page 14) is of particular importance for the installation and operation of the incubator.**

**In addition to the above, please observe all relevant local safety regulations.**

**The incubators must not be used for drying processes or other heat treatments releasing combustible vapors which might form explosive mixtures with the air. Furthermore, the incubators are not designed for heat treatments of dangerous substances.**

If you wish to use the incubator for purposes other than those for which it has been designed, please contact the nearest Heraeus Service representative.

## 2. Installation and Connections

**Incubators of thermal safety class 0 (see nameplate) must be installed so as to permit permanent attendance, cf. section "Thermal protection" (Page 14).**

If you frequently work with cultures or test materials which are particularly important, difficult to replace or even irreplaceable, appropriate precautions should be taken (**protection by temperature limit controller**).

### Transport

Care should be taken not to lift the incubator by the door or the door handle.

### Connection to the Power Supply

The incubators come with a power supply cable and a shock-proof plug (16 A).

For line fusing, a 16 A delay-action fuse is required.

Prior to connecting the unit to the line, make sure that the power supply ratings are those stated on the nameplate. The connection must be made in accordance with the relevant local safety regulations.

Rated voltage, rated frequency, rated power and current are listed under "Technical Data" and on the nameplate.

### Place of installation

The incubator must be placed on a level surface.

### Adjustable feet

Its position is locked in if the adjustable feet are turned. This requires two wrenches (Size 24).


### Articulated casters

Models B 6420 and B 6760 are fitted with 4 articulated casters. The lever which releases the caster is situated above the locking lever. After the incubator has been positioned, restrain the locking levers of the casters.

### Stability (models 6420, 6760)

In order to ensure optimum stability, turn the front casters forward as soon as the incubator is in position, and depress the locking levers (see fig. 1).

Close the door(s) when the incubator is being moved!

**Ambient temperature**

The ambient temperature should be at least 5°C lower than the intended working temperature.

Avoid installation in drafty locations. The ambient temperature should not permanently exceed 35°C.

The rear wall of the incubator is fitted with 2 tubular sockets. When the air flap is set for fresh air operation, cold air passes through the lower socket and hot air exits through the upper one.

Air inlet pipe

An air inlet pipe can be connected to supply the incubator with fresh air.

Air outlet pipe


The vapors produced in the incubator during the heating phase can be expelled through an air outlet pipe.

**Vapors and gases which are detrimental to health and the environment must be carried off by all means.**

For this purpose, a temperature-resistant and corrosion-proof air outlet pipe should be used, which can be connected to the air outlet socket.

We recommend that you install a draft interruptor if several incubators in a row are connected to a central exhaust manifold (see illustration).

Bench installation

Incubators (B 6060 and B 6120) can also be built into laboratory benches. For safety reasons, the laboratory benches must be of non-combustible materials in accordance with DIN-Standard 4102. Installation is carried out as follows:

- \* Connect the incubator to the line.
- \* Fit it underneath the bench, taking care not to damage the power supply cable.

A spacer provides for clearance between the rear wall of the laboratory bench and the incubator. Clearance all the way around the unit must be at least 25 mm.


## Stacking

### Accessories required:

Models:	6060	6120	6200
Stock no.:	50027660	50027661	50027661

Incubator models 6060, 6120 and 6200 can be stacked as follows:

A max. of 2 on top of each other and 2 stacks side-by-side; clearance between stacks 25 mm.


The procedure is carried out in the following manner:

Remove the plastic rivets.


Put down the facing.

Secure the stacking elements and angular facing with head screws.

Place the second incubator on top of the first one. The stacking elements fit into the hollow feet of the upper incubator to secure it in place.

3. Operating Elements from page 8 - 20

3.1) Air-flap adjustment

3.2) Temperature registration

3.3) Temperature control

3.4) Upper/lower limit safeguard

3.5) Timers

3.6) Master switch

### 3. Operating Elements

#### 3.1 Air-flap adjustment


open

closed

#### 3.2 Temperature registration


Temperature indicator

dial thermometer

20 to 50°C

Temperature recorder

round-chart  
temperature recorder,  
measuring range:

0 to 100°C

The round chart recorder works without ink or color ribbon. The curve is recorded as the sharp needle point presses periodically against the pressure sensitive record chart. But for the replacement of the record chart, this system is absolutely maintenance-free.

The record charts are graduated from 0 - 100°C. Future orders require that the article number be stated.

#### For 50 Hz

1-day record chart, art. no. : 50029351

7-day record chart, art. no. : 50029352

30-day record chart, art.no. : 50029353

#### For 60 Hz

20h record chart, art. no. : 50029581

140h record chart, art. no. : 50029582

25-day record chart, art. no.: 50029583

Minimum quantity: 1 set (100 record charts).

#### Replacement of the record chart:

A slight clockwise rotation detaches the plexiglass cover. After removal of the record chart, the following operating elements become visible:

#### Stroke adjustment:

Setscrew for adjustment of the stroke; clockwise rotation increases the stroke.

### Feed selection

This slide is used to set the desired feed.

- 1 day / 7 days / 30 days (at a line frequency of 50 Hz).
- 20 h / 5.8 days / 25 days (at a line frequency of 60 Hz).


The pointer tip is visible in the narrow slot beneath the stroke-adjustment screw, and appears every 6 seconds to register the measuring point. Depending on when it is disconnected, the pointer could stop in the recording position.

**In order to avoid damage to the pointer when it is in the recording position, make sure that a record chart has been inserted before the plexiglass cover is closed.**

After the desired feed has been set, the new record chart is inserted and pushed onto the four guide pins in the center. Clockwise rotation of the driving shaft sets the time at the circumference of the record chart to the white mark at the edge of the instrument. A slight clockwise turn will lock on the plexiglass cover. (Note the inscription "ELMES 12" in the center.)

**3.3 Temperature Control**  
Mechanical temp. controller  
"Heratherm"  
Models 6030 - 6200 only

Temperature range:  
20 - 50°C

**Setting:**

Use a screwdriver to release the locking device of the temperature controller.

**Note the marking !**

Select the desired working temperature with the control knob.

Tighten the locking device again.


**Note the marking !**

The set working temperature, once reached, is kept constant with the controller opening and closing the heating circuit as required. The white pilot lamp indicates that the heating is working.

If the chamber temperature and dial readout do not match, the setting must be corrected. (see Appendix).

Electronic temperature controller "Digicon"

electronic controller w. digital temp. display for temp. range: 20 - 70°C

Setting:

The electronic controller features digital display of the set and actual temperatures.

Press the control button and select the desired working temperature. The set value is indicated in the digital display.

When the control button is released, the actual value appears.

The LED (top left) indicates that the heating system is working.

The "Digicon S" temperature controller controls the temperature of the incubator in accordance with the setpoint. The actual temperature is indicated digitally. When the (blue) setpoint button is pressed, the display will indicate the desired working temperature.

The interface module is located above the "Digicon S" temperature controller (in models B 6030, 6420 and 6760, it is next to it).


Connect connecting cable (diode plug) to socket "W<sub>ext</sub> (0..20mA)".

0 mA = 0°C      20 mA = 70°C  
Admissible load R<sub>i</sub> = 500 Ω

This connection converts from the internal to the external setpoint program. There must be a bridge between pins 3 and 4 in the diode plug.

The actual temperature can be measured at the two "X<sub>ext.</sub>" sockets as voltage of 0-10V.

0 Volt = 0°C      10 Volt = 70°C  
Admissible internal resistance R<sub>i</sub> = 1 kilo-Ω

- 1 + 0...20 mA
- 2 - 0...20 mA
- 3 setpoint
- 4 ext. setpoint
- 5
- 6 + 0...10 V
- 7 - 0...10 V

50012338

Digicon S with interface

Operation with external program

Measurement of the actual temperature

View of module with pin assignment



Connect only power sources / measuring amplifiers which are fed from a safety circuit (cf. here VDE-Regulation 0100/410).

(Example: via safety transformer and VDE-Regulation 0551.)

Microprocessor-based temperature controllers

"Thermicon T"

"Thermicon P"  
(programmable)

For description and start-up, see enclosed operating instructions.

"Eurotherm  
type 818"

For description and start-up, see enclosed operating instructions.

### 3.4 Thermal protection

The requisite precautions concerning thermal safety are stipulated in DIN-Standard 12 880, Part 1, Table 3 (see Appendix).

The thermal safety class which applies to these incubators is listed on the nameplate.

#### No overtemperature protection

##### Thermal safety class 0:

No direct protection for the incubator, its environment and load.

##### **Only admissible for non-dangerous loads.**

- If the temperature control system fails, the subsequent rise of the temperature inside the incubator will lead to damage.
- Only admissible for attended operation.
- The operation of these incubators must be checked at sufficiently short intervals.
- Incubators of this thermal safety class may not be built into laboratory furniture or stacked.

Supplementary equipment can be added to the incubators of safety class 0, if a temperature limit controller or an upper/lower limit controller is installed.

A Heraeus Serviceman can carry out such an addition at any time.

Temperature limit controller for excess temperature


Thermal safety class 3.1

Protection of the incubator, its environment and the load against inadmissible excess temperature.

**Admissible only for non-dangerous loads.**

Function:

The temperature limit controller is operationally and electrically independent of the temperature control system. In case of failure, it takes over the control function.

Temperature limit controller

When the control knob is set to the limit stop, the temperature limit controller functions as incubator protection. If the temperature limit controller is set to a value somewhat higher than that selected on the controller, it functions as load protection.

Once the temperature limit controller begins to function and the red warning light comes on, proceed as follows:

1. Disconnect the incubator from the line.
2. Investigate and eliminate the cause of the failure.
3. Turn on the incubator again as described.

Setting:

To ascertain the temperature at which the temperature limit controller operates, and thus to be able to adjust it properly, the incubator must be put into operation. Set the control knob to the limit stop (for protection of the incubator). After the selected operating temperature has been reached, reset the temperature limit controller to the switch point. If the controller comes into operation, the red warning lamp lights up.

You'll find the optimum setting of the temperature limit controller by turning back the control knob by approximately one graduation mark. The red warning lamp will go out.

Operational test:

The temperature limit controller should, at reasonable intervals, be tested for proper operation. It is recommended that this test also be carried out routinely by the authorized operating personnel, e.g. prior to the start of a longer work process.

Connection for central monitoring system (potential-free contact)

Provisions have been made to connect customer-installed warning systems (e.g. telephone system, building control systems).

Function:


When the thermal safeguard (temperature limit cut-out, temperature limit controller) comes into action, or a timer switches the incubator on or off, the appropriate message will be transmitted to the customer-installed warning system.

The potential-free contacts (1 make and 1 break contact) are designed for the following circuits:

- \* Power supply circuits:  
max. 250 V A.C., 3 A
- \* SELV/SELV-E circuits  
(cf. VDE-Regulation 0100/410):

25/50 V A.C., 4 A  
60/120 V D.C., 0.5 A

see also connection diagram

Upper/lower temperature limit controller functioning in case of excess and insufficient temperature


Thermal safety class 3.3

Protection of the incubator, its environment and the load against inadmissible excess and insufficient temperature.

**Only admissible for non-dangerous loads.**

Function:

This controller complements the function of the temperature limit controller as per thermal safety class 3.1. With the flip switch, it can be turned on/off as desired to offer additional protection against insufficient load temperature if the controller fails.

Lower limit safeguard

Upper limit safeguard

Temperature limit controller with upper and lower limit safeguard.

Setting of the upper limit safeguard:

Set the desired temperature for take-over in the event of overheating as described under "temperature limit controller" (excess temperature).

Setting of the lower limit safeguard:

Turn on the white switch after the set temperature has been reached.

Pilot lamp in the switch lights up.

A slow turn of the knob will reveal the trigger position on the dial of the lower limit controller (the orange warning lamp lights up).

After that, turn back the knob a little so that the orange lamp will go out.

Operational test:

The upper/lower limit controller should, at reasonable intervals, be tested for proper operation. It is recommended that this test also be carried out routinely by the authorized operating personnel, e.g. prior to the start of a longer work process.


### 3.5 Timers: Synchronous timer

#### Function:

Synchronous time switch  
for switching the oven on and off

6 h duration at 50 Hz line frequency  
5 h duration at 60 Hz line frequency

24 h duration at 50 Hz line frequency  
20 h duration at 60 Hz line frequency

#### 24h-Clock

#### Setting:

Function of the mode switch

Position I: To turn on the incubator after the set time has expired, or for operation without a synchronous timer.

Position II: To turn off the incubator after the set time has expired.

#### 6h-Clock

For shut-off only.


Set the timer to the desired time by turning the knob to the right. Make sure to exceed the minimum adjustment angle of approx. 30°.

Shorter switching times are obtained if the timer is first turned to the right and then turned back to the desired time.

#### Caution:

The premounted dial (24 h) applies to a line frequency of 50 Hz. A line frequency of 60 Hz requires that the enclosed 20 h dial be substituted.

The button is released from the dial by impact from below (snug fit).


Mounting of the dial:

- Set timer to 0 position
- Mount dial in 0 position
- Fasten the control button.

**Daily program timer**

**Function:**  
To switch the incubator on and off.

**Setting:**  
At the start of operation, set the actual time on the program dial by turning it to the right.


Select the desired switching times by locking down the riders.

The locked red dial segments indicate the on-time.

Minimum switching frequency 15 minutes.

**Weekly program timer**

**Function:**  
To switch the incubator on and off.

**Setting:**  
Set the actual time and weekday on the program dial by turning it to the right.

1 = Mon., 2 = Tues., 3 = Wed., 4 = Thurs.,  
5 = Fri., 6 = Sat., 7 = Sun.

Select the desired switching time by locking down the riders.

Minimum switching frequency 2 hours.

**Mode converter:**

- ☉ = automatic operation with timer function.
- ☒ = manual operation: no timer function.

**Digital weekly program timer**


See separate operating instructions.

### 3.6 Master switch:

#### Function:

Turn the knob to switch the incubator on/off at all poles.


Turn on the master switch  
- The green pilot lamp lights up.

#### 4. Commissioning

##### Shelves and supports

The shelves and their supports are secured inside the incubator for transport. The transport guards must be removed before the incubator is put into operation.

The shelf supports can be inserted at random in the perforated carrier racks. The shelves slide over the supports in such a way that the fork-shaped stabilizer catches the lower side of the support.

**The working space proper is defined as that space which is 1/10 less than the overall chamber of the respective incubator (cf. also DIN-Standard 12 880, Part 2).**

Only this portion of the chamber may be loaded.  
**Do not load the bottom of the incubator!**

Turn on the master switch.

- The green pilot lamp lights up.

Set the temperature controller to the desired working temperature.

After the desired temperature has been reached, set the temperature safeguard.

**Note on  
loading the unit**

**Loading the incubator correctly greatly influences  
the temperature conditions inside the chamber.**

In practice, incubators are often loaded not only beyond the lateral boundaries of the working space, but also with stacks of Petri dishes which are too high.

If culture vessels are stacked, their surfaces come into contact with each other (e.g. Petri dishes). This means the stated temperature variation will be exceeded if the stacks consist of three or more single items.

However, it is not only important that the preset temperatures be kept inside the chamber; in some instances, temperature must recover within an appropriate period of time after the test materials have been introduced. This point is less important when dealing with long-term incubation projects (over 72 hours) so that more dishes can be crammed into the available space.

With relatively short incubation projects, on the other hand, (up to 72 hours), overloading must be avoided, because the temperature recovers too slowly even if there is only one layer of Petri dishes. Here in particular, attention must be paid to the fact that good temperature recovery is guaranteed only if the test material is distributed evenly throughout the entire working space.

Don't arrange the load too closely on the shelves (load only 70% of their area), to not impair air circulation and to ensure uniform heating.

## 5. Maintenance and Repairs


The incubator should, at reasonable intervals, be tested for proper operation. This applies especially to the temperature limit controller.

**The safety of the incubator and its conformity with the relevant standards will be guaranteed only if maintenance work, repairs and/or modifications are carried out by authorized personnel.**

Before any maintenance or repair work is carried out, the incubator must be disconnected from the line.

### Readjusting the door

If the door no longer closes properly, adjust the bearing pin and the catches:

#### 1. Adjusting the catch

- Loosen the catch nut M6 with a wrench (size 10) and turn the catch clockwise for at least one complete rotation. After this, return the catch to the same position. Tighten catch nut M6 again.

#### 2. Adjusting the bearing pin

- Loosen the Phillips screw of the bearing pin. Use a small spindle to prevent the bearing pin from rotating also. After you've loosened the screw, turn the bearing pin with the spindle, and determine the desired setting for height and depth. Retighten the Phillips screw, and check if the door is tight.

### Removal of the glass door

Open the door, slightly pressing upward. Note the withdrawal position at approx. 15°. The door can be removed and refitted only in this position.

### Replacing the gasket

Remove the defective gasket and install a replacement. After that, it is imperative that the unit be tested for tightness. Make sure to replace the clamps.

### Removal of the internal fittings

Use a Phillips screwdriver to remove the fittings inside the working chamber.

- Unscrew the bottom plate
- Unscrew the lateral cooling baffles
- Remove the individual internal fittings

Hints on care

Inner casing:

Commercially available detergents in small quantities should be used, but **no acids, chloric solvents or saline solutions.**

Control elements (modules):

It is best to wipe these only with a damp cloth.

Spare parts

If you have complaints or wish to order spare parts, please state the data on the nameplate.

The safety and correct operation of the incubator are guaranteed only with genuine Heraeus spare parts.

A list of spare parts can be requested either from your local Heraeus representative or directly from the factory. All you need do is to state the model type and serial number.

Wiring diagrams:

BR 6000	220 V/1 A.C.	500 11 620
BR 6000	110 V/1 A.C.	500 30 370

**6. Appendix - Table 3, from DIN-Standard 12 880**

Safety devices against malfunction of the temperature control circuit.

Classification of heating ovens according to DIN-Standard 12 880 Part 1, table 3.

The class describes the type and scope of protection, as well as the safety measures required to deal with a simple defect in the temperature control circuit.

Class	Purpose	Scope	Safety device acc. to Sub-section 5.5	Safety measures
0	-	-	-	Only attended operation*) with harmless materials is permitted; or overheating is excluded by constructional measures
1	Protection of oven	In the event of a defect, no danger emanates from the oven.	Temperature limiter or thermal fuse	Special safety measures depending on the intended use.
2	Protection of oven, environment and material being treated.	In the event of a defect, no danger emanates from the oven or from the material being treated.	Adjustable temperature limit cut-out	
3		In the event of a defect, the material being treated is protected against overheating and/or undercooling (e.g. incubator)	Adjustable temperature limit controller, coming into action in case of:	
3.1				
3.2	insufficient temperature			
3.3		excess temperature and insufficient temperature		

\*) If operation is supervised, the unit must be checked at reasonably short intervals to make sure it works correctly.

	B 6030	B 6060	B 6120	B 6200	B 6420	B 6760
<b>Materials:</b>						
Inner casing*)	1.4301	1.4301	1.4301	1.4301	1.4301	1.4301
Shelves*)	1.4301	1.4301	1.4301	1.4301	1.4301	1.4301
Door gasket	Silicone rubber	Silicone rubber	Silicone rubber	Silicone rubber	Silicone rubber	Silicone rubber
Heating elements*)	1.4435	1.4435	1.4435	1.4435	1.4435	1.4435
Outer casing**)						
*) Stainless and acid-resistant steel						
**) Galvanized sheet steel with temperature-resistant outer coating						
<b>Colours:</b>						
Outer casing, grey-white	RAL 9002	RAL 9002	RAL 9002	RAL 9002	RAL 9002	RAL 9002
Modules, agate-grey	RAL 7038	RAL 7038	RAL 7038	RAL 7038	RAL 7038	RAL 7038
Operating elements, light-blue	RAL 5012	RAL 5012	RAL 5012	RAL 5012	RAL 5012	RAL 5012
<b>Dimensions</b> see dimensioned sketches P.26						
<b>External dimensions:(mm)</b>						
Casing						
T <sub>A</sub> incl. door handle	610	610	610	790	790	790
B <sub>A</sub>	552	744	895	895	744	1200
H <sub>A</sub>	552	552	696	816	1813	1813
<b>Internal dimensions:(mm)</b>						
t <sub>i</sub>	370	370	370	550	550	550
b <sub>i</sub>	352	403	554	554	544	1000
h <sub>i</sub>	231	380	524	644	1366	1366
<b>Working-chamber volume (liters)</b>	30	57	107	196	409	751

50012338

**Technical Data**

	B 6030	B 6060	B 6120	B 6200	B 6420	B 6760
<b>Weights (kg)</b>						
Empty weight	41	52	69	96	157	239
max. charge	50	50	50	75	75	150
max. shelf load:						
surface	20	20	20	20	30	30
spot	10	10	10	10	15	15
<b>Electrical data:</b>						
Rated voltage V	220	220	220	220	220	220
Rated freq. Hz	50/60	50/60	50/60	50/60	50/60	50/60
Rated current A	0.95	1.6	1.6	2.4	3.4	5.1
Rated power kW	0.21	0.36	0.36	0.54	0.75	1.13
<b>Times: (minutes)</b>						
Rise time (air-flap closed) from 25°C to:						
70°C	70	35	45	70	80	84
50°C	35	30	35	60	72	83
37°C	20	25	35	40	65	67
Recovery time (door opened for 30 sec.) to 99% of the initial value						
at 70°C	5	3	2	5	22	15
at 50°C	3	3	2	5	17	6
at 37°C	2	3	2	5	13	6

50012338





**Technical Data**

	B 6030	B 6060	B 6120	B 6200	B 6420	B 6760
<b>Air changes per h</b> at rated temp. air flap open	5	5	5	5	5	5
<b>* Holding power (Wh/h)</b> recirculated air at 70°C	50	68	97	120	232	335
at 50°C	30	41	57	71	126	181
at 37°C	19	25	35	42	71	104
<b>Temperatures**:</b>						
	mech./electron controller	mech./electron controller	mech./electron controller	mech./electron controller	mech./electron controller	mech./electron controller
Rated temp. (°C)	50 / 70	50 / 70	50 / 70	50 / 70	50 / 70	50 / 70
Variation in space (°C), air flap <u>closed</u>						
at 70°C	± 1.7	± 2.2	± 2.5	± 2.5	± 3.4	± 4.7
at 50°C	± 0.7	± 1.4	± 1.6	± 1.6	± 2.4	± 2.9
at 37°C	± 0.5	± 0.8	± 1.1	± 0.9	± 1.4	± 1.8
Fluctuation at rated temperature (°C)	≤ 0.5	≤ 0.5	≤ 0.5	≤ 0.5	≤ 0.5	≤ 0.5
Accuracy of temp. display Telethermometer in % of full-scale reading	± 3	± 3	± 3	± 3	± 3	± 3
Digital display in % of full-scale reading	± 1	± 1	± 1	± 1	± 1	± 1

\* corresponds to the heat load of the installation area

\*\* measured in accordance with the DIN-Standard 12 880, Part 2.  
50012338

# Dimensioned Sketches

Baugrößen 6030, 6420, 6760

Baugrößen 6060, 6120, 6200


## Appendix

### Adjustment of the mechanical controller "Heratherm"

If the temperature in the working space does not correspond to that on the dial of the mechanical controller, the dial can be readjusted.


To do this, the following steps are necessary:

1. Remove the blue button with the dial
2. Push the blue button off the dial (see illustration below)
3. Reposition the dial on the module
4. Compare the temperature value with the dial value and readjust
5. Remount the blue button

## Eurotherm Controller with Interface RS 232

Assignment of the 25-pin D-Sub socket:

Pin	1	Shield
	2	Rx ( reception )
	3	Tx ( transmission )
	4	( connected to 5 )
	5	( connected to 4 )
	6	( connected to 8 and 20 )
	7	COM ( common, ground)
	8	( connected to 6 and 20 )
	20	( connected to 6 and 8 )


----- All unlisted pins are not wired -----

The interface pin assignment is designed to permit direct connection to an IBM or IBM-compatible computer if the computer is equipped with a 25-pin plug. If the computer interface has a 9-pin plug, a commercially available adapter (9-pin socket and 25-pin plug) can be used.

Upon delivery of the unit, the controller address is set to 00, and the Baud rate to 9600. For other addresses or other Baud rates, the controller must be reconfigured (see operating instructions for the EURO THERM controller).

## Eurotherm Controller with Interface RS 485

Assignment of the 25-pin D-Sub socket:

Pin	1	Shield
	2	( connected to 12 )
	3	Tx + ( transmission + )
	4	( connected to 5 )
	5	( connected to 4 )
	6	( connected to 8 and 20 )
	7	COM ( common, ground)
	8	( connected to 6 and 20 )
	12	Rx + ( reception + )
	13	Rx - ( reception - )
	16	Tx - ( transmission - )
	20	( connected to 6 and 8)


----- All unlisted pins are not wired -----

Upon delivery of the unit, the controller address is set to 00, and the Baud rate to 9600. For other addresses or other Baud rates, the controller must be reconfigured (see operating instructions for the EURO THERM controller).





**Caution:** In the case of some RS 422/485 interface cards for personal computers, it may be necessary to change the **polarity** of the signals Rx and Tx. If the controller does not respond to the computer interface, connect Rx+ of the PC interface to Tx- of the controller, and Tx+ of the PC interface to Rx- of the controller.

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