

Telstar BiOptima IT configuration

The Telstar solution for Hospital Pharmacies in applications that involve the use of pharmacy compounding software programs.

The majority of hospital pharmacies rely on software solutions to aid them in dose preparation and improve the safety of the medications and the pharmacy workflow.

These software solutions can include on-line (weight-based) gravimetric verification, barcode reading, automatic and electronic documentation of the preparations.

IT configuration is a recent addition to the Bioptima family of Microbiological Safety cabinets specifically engineered to address the current needs of pharmacies.

This configuration is valid for all Bioptima models: One, Cyto, Exa and Biopima.

Designed and manufactured according to EN 12469 and the DIN 12980 norms these units are equipped with all the necessary hardware for effortless integration in the daily routine of a modern day pharmacy or to upgrade existing pharmacies to a state of the art compounding facility.



Top Features



**Swiveling support
for keyboard**

**Support for
printer and PC**



Integrated 21.5" monitor



Telstar[®] BiOptima IT configuration

Available sizes

The IT configuration is available for sizes 4,5 & 6 in all BiOptima models: One, Cyto, Exa and BiOptima

Main features

BiOptima IT series are standard equipped with:

- 2 European sockets.
- UV light
- Segmented Stainless Steel work surface
- Support frame H 900mm
- LCD Monitor: 21.5" screen integrated in the back wall of the cabinet, with HDMI connection to the outside
- 2x USB Feed-through
- Full size pre-filter under the work area
- Electrically operated front window with unique back-up battery and safety mechanism

Options

- Support for printer and PC
- Swiveling support for keyboard

Accessories

- Components for ducting
- Taps
- IV bar
- Weighing stone (low vibration for balance)



Check for the nearest Telstar office at
<https://www.telstar.com/international/>

Headquarters

Av. Font i Sagué, 55 T +34 937 361 600
08227 Terrassa (Spain) F +34 937 861 380

www.telstar.com
sales.laboratory@telstar.com



ISO 9001: Certified Company

FL-BIOPTIMAIT-EN-0122
Telstar reserves the right to improvements and specifications changes without notice.