




## MAGNETIC SEPARATION

BIOVENDIS

XXL



### Type 1

Magnetic separator for manual separation using Kits in high volumes. Consists of a clear polyacril block with 2 holders for 50 ml tubes: one duct with one magnet for separation from a maximum 25 ml volume and one with two magnets for separation of a maximum 50 ml liquid volume.

Optimized for use with 50 ml tubes and Kits or other kits using magnetic particles. Transparent material allows for constant visual control of the separation.

### Type 2

Magnetic separator for manual separation using Kits in high volumes. Consists of a clear polyacril block with 2 holders for 50 ml tubes: each product has one magnet for separation from a maximum 25 ml liquid volume.

Optimized for use with 50 ml tubes and Kits or other kits using magnetic particles. Transparent material allows for constant visual control of the separation.

**IMPORTANT!** We recommend cleaning the Magnetic Stand M96 with deionized water only; alcohol-based cleaning agents may cause cracking of the stand material.

**For Research Use Only. Not for human or animal therapeutic or diagnostic use.**

	XXL Type 1	XXL Type 2
Cat. No.	209616	209617
Format	2x 50ml tubes	2x 50ml tubes
Material	clear polyacril	clear polyacril
Colour	transparent	transparent
Items per Unit	1	1
Criteria	up to 50ml volume	up to 50ml volume