

DATASHEET

For in vitro research only.
Not for diagnostic or clinical use.

Optical Clear Sealing Tape

100 sheets

Cat. No. 317500

BIOVENDIS

Biovendis Products GmbH
Joseph-Meyer-Str. 13-15
68167 Mannheim

Tel. + 49 - 621 - 43 00 52 52
Fax + 49 - 621 - 43 00 52 50
E-Mail info@biovendis.de

PRODUCT

Optical Clear Sealing Tape

ORDERING INFORMATION

Cat. No. 317500 - Optical Clear Sealing Tape

PLATE MATERIAL

Polypropylene (PP), Polyethylene (PE)

SEALING CONDITIONS

Bond strength depends on application pressure, so sufficient pressure must be applied.

SEAL INTEGRITY RANGE

-70 °C to 110 °C;

During PCR these seals can withstand higher temperatures (> 110 °C) provided adequate support/pressure is applied by the thermal cycler heated lid.

DESCRIPTION

The Optical Clear Sealing Tape is an adhesive plate seal, specially developed for qPCR and other fluorescence applications.

The strong adhesive prevents evaporation during thermal cycling.

No sticky residue on the plate when the seal is peeled off after PCR.

Certified free from DNase, RNase and human genomic DNA contamination.

STORAGE

Seals should be stored in a cool, dry place. Avoid direct exposure to sunlight. It is recommended to use seals within one year from date of purchase.

SHEET SIZE

141,3 x 79,4 mm (with tabs)

120,7 x 79,4 mm (without tabs)